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transcriptional regulation of myeloid differentiation. Experimental hematology, 36(11), pp.1449-1460.

- Groner, Y., Ito, Y., Liu, P., Neil, J.C., Speck, N.A., Andre Van Wijnen and Springerlink., 2017. RUNX Proteins in Development and Cancer. Singapore: Springer Singapore.
- Lin, S., Mulloy, J.C. and Goyama, S., 2017. RUNX1-ETO Leukemia. RUNX proteins in development and cancer, pp.151-173.
- Mackay, Ian M., 2004. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3: 190-212.

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AML1-ETO RQ

Autumn 2025, Version 1.2

For Real-Time PCR Quantitative Detection of
AML1-ETO Transcripts
For Research Use Only

 24 (Cat# AML1-ETORQ24)

 48 (Cat# AML1-ETORQ48)

 96 (Cat# AML1-ETORQ96)

 NG-WI-ASL-31-102

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. RNA Isolation	7
14. cDNA Synthesis	7
15. PCR Protocol	7
16. Devices and software	8
17. Programming RotorGene	8
18. Programming of StepOne	9
19. Programming Other Machines	10

20. Data Analysis: RotorGene.....	10
21. Data Analysis: StepOne	12
22. AML1-ETO% calculation	15
23. Analytical Sensitivity	15
24. Disposal Method	15
25. Technical Support	16
26. Contact Information.....	16
27. References	16
28. Symbols	17

1. Introduction

AML1-ETO RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying AML1-ETO Transcripts. All required reagents are included in the PCR Mix provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to calculate BCR-ABL% and prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

Important Note: *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

2. Intended Use

AML1-ETO RQ kit is intended for the detecting of AML1-ETO transcripts in peripheral blood and AML1-ETO percentage calculation in patients undergoing therapy. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, MIC and StepOne machines.

3. Background Information

AML1-ETO is the most common chromosomal aberration in de novo acute myeloid leukemia (AML) patients. This abnormality is resulted from t(8;21)(q22;q22) translocation. The translocation produces a fusion protein which inhibits myeloid transcriptional factors and as a result, cellular differentiation is blocked. This alteration occurs in approximately 7% of adult and 12% of pediatric with AML.

4. Test Principle

The Translocation is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time

PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual or a quick guide and the following reagents:

Label	Content	Quantity
AML1-ETO Mix*	Master mix for AML1-ETO	480 µl
ABL Mix*	Master mix for ABL	480 µl
AE1	AML1-ETO Standard 1: 100,000 copy/µl	150 µl
AE2	AML1-ETO Standard 2: 10,000 copy/µl	150 µl
AE3	AML1-ETO Standard 3: 1,000 copy/µl	150 µl
AE4	AML1-ETO Standard 4: 100 copy/µl	150 µl
AE5	AML1-ETO Standard 5: 10 copy/µl	150 µl
A1	ABL Standard 1: 100,000 copy/µl	150 µl
A2	ABL Standard 2: 10,000 copy/µl	150 µl
A3	ABL Standard 3: 1,000 copy/µl	150 µl
A4	ABL Standard 4: 100 copy/µl	150 µl
Water	PCR Grade Water	200

6. * 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**

- Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2 ml PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 48 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used. For optimum results, a sample should include 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for two days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, Germany)
- TRIzol isolation reagent (Cat. no. 15596026, Thermochemical, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

14. cDNA Synthesis

4ug of total RNA is required and should be reverse transcribed to cDNA using random hexamers. Different kits are available in the market for this purpose.

Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20ul of cDNA add 30ul of nuclease free water.

Upon successful RNA extraction and cDNA synthesis, ABL titer will reach above 20,000 copy/ul.

15. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

Each sample should be examined for both of AML1-ETO fusion gene and for ABL control gene. So, two sets of reactions are required. In AML1-ETO set, consider 1 tube for each sample as well as 6 tubes for the 5 standards (AE1 to AE5) and Negative control or NTC. In ABL set, consider 1 tube for each sample and 5 tubes for the 4 standards (A1 to A4) and Negative control or NTC. Place required number of tubes on cold block.

Pipette 20µl of AML1-ETO Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA or standard or control to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.

16. Devices and software

AML1-ETO RQ kit is designed to work with Rotor-Gene, StepOne and MIC.


17. Programming RotorGene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the p210 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); AML1-ETO 0.1 is for strip tubes and AML1-ETO 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the next page image.

Select a tube number containing bcr Mix for the Green channel and tube with ABL Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

18. Programming of StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. One negative control, five standards for AML1-ETO, four standards for ABL and a few samples are defined. You may change plate setup using right-click options

(copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished click on Start Run and save the experiment. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Both of AML1-ETO Mix and ABL Mix contain ROX with final concentration of 300nM in reaction.

20. Data Analysis: RotorGene

Before analyzing results, make sure that in the sample menu, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC”, respectively.

Analyze the data according to RotorGene manual. Perform quantitative analysis for both **AML1-ETO (the Green channel)** and **ABL (the Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double-click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold at 0.05. Figures 1 and 2 represent typical graphs for the RotorGene machine.


Fig 1. Typical AML1-ETO graph in Green channel for Rotor-Gene


Fig 2. Typical ABL graph in Green channel for Rotor-Gene

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the Green/AML1-ETO and Yellow/ABL channels with sigmoid graphs and CT of 20-40 for the Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the Green/AML1-ETO channel while it is positive in the Yellow/ABL channel with a sigmoid graph and CT of 20-27 and ABL titer above 20,000 copies/ul..
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green/AML1-ETO and the Yellow/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green/AML1-ETO channel while it is positive in the Yellow/ABL channel with sigmoid graph and CTs above 27 and ABL titer less than 20,000 copies/ul. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 1ug RNA has been used. CTs higher than 27 and ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **false negative** reports.*

21. Data Analysis: StepOne

Analyze the data according to StepOne manual. Briefly, click on "Analyze" and set the threshold for the **AML1-ETO/FAM** on 0.05 and for the **ABL/VIC** on 0.1.

Figures 3 and 4 represent typical graphs for the StepOne machine.


Fig3. Typical AML1-ETO graph in FAM channel for StepOne


Fig 4. Typical ABL graph in VIC channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the FAM/AML1-ETO and the VIC/ABL channels with sigmoid graphs and CT of 20-40 for the Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the FAM/AML1-ETO channel while it is positive in the VIC/ABL channel with a sigmoid graph and CT of 20-27 ABL titer above 20,000 copies/ul.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/AML1-ETO and VIC/ABL channels. Improper extraction or test setup could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM/AML1-ETO channel while it is positive in the VIC/ABL channel with sigmoid graph and CT of above 27 and ABL titer less than 20,000 copies/ul. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 1ug RNA has been used. CTs higher than 27 and ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **false negative** reports.*

22. AML1-ETO% calculation

To assess the response to therapy, AML1-ETO% value for each patient can be calculated. This kit uses NCN method for this purpose (*Beillard E. 2003, Leukemia 17:2474*). In this method AML1-ETO% value is the AML1-ETO transcripts (titer) normalized by the ABL transcripts (titer) and then multiplied by 100%.

Usually, titer of ABL gene transcripts is higher than for AML1-ETO fusion gene. Therefore, results of above calculation would be a number less than 100%. This value may fall even below 0.001% in optimum cases. However, at diagnosis and in case of relapse, AML1-ETO transcripts may surpass ABL gene transcripts. As a result, ratio of AML1-ETO/ABL would be above 100%. Such results have been reported previously too. For example, please note table 12 of the paper published by J. Gabert *et al* (2003, *Leukemia 17:2318*) for patients. He reported a range of 110% to 2000% with average of 420% for peripheral blood and range of 29% to 530% with average of 210% for bone marrow samples. Therefore, a ratio above 100% is not unexpected and simply means higher transcripts of the fusion gene compared to the control gene.

23. Analytical Sensitivity

Analytical sensitivity of the kit has been determined by examining dilution series of cloned cDNA and was estimated to be about 2 copeis/ul.

24. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

25. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

26. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393


+98 990 -1813124

Website: www.novingene.com

27. References

- Balkhi, M.Y., Christopeit, M., Chen, Y., Geletu, M. and Behre, G., 2008. AML1/ETO-induced survivin expression inhibits transcriptional regulation of myeloid differentiation. *Experimental hematology*, 36 (11), pp.1449-1460.
- Groner, Y., Ito, Y., Liu, P., Neil, J.C., Speck, N.A., Andre Van Wijnen and Springerlink, 2017. RUNX Proteins in Development and Cancer. Singapore: Springer Singapore.
- Lin, S., Mulloy, J.C. and Goyama, S., 2017. RUNX1-ETO Leukemia. RUNX proteins in development and cancer, pp.151-173.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3, 2004: 190-212.

28. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com