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Aspergillus RQ4 Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of Aspergillus DNA
For use with Rotor-Gene 4 plex or higher
For Research use only

 24 (Cat# Aspergillus4RQ24)

 48 (Cat# Aspergillus4RQ48)

 96 (Cat# Aspergillus4RQ96)

 NG-WI-ASL-57-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	5
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. Internal Control (IC)	7
14. DNA Isolation	7
15. PCR Protocol	8
16. Devices and software	8
17. Programming Rotor-Gene	8
18. Programming Other Machines	9
19. Data Analysis: Rotor-Gene	10

20. Sensitivity.....	13
21. Disposal Method	14
22. Technical Support.....	14
23. Contact Information.....	14
24. References	14
25. Symbols	15

1. Introduction

Aspergillus RQ4 kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting Aspergillus DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

Aspergillus RQ4 kit is intended for the detection of four Aspergillus species including *A. fumigatus*, *A. flavus*, *A. terreus* and *A. niger* by detecting the ITS region of the fungus (Without differentiation between *A. terreus* and *A. niger*). Detection is accomplished by polymerase chain reaction and is applicable with Rotor-Gene 4-plex and MIC machines.

3. Background Information

Aspergillus is an opportunistic type of fungal mold that exists everywhere. More than a few hundred Aspergillus species have been identified. Among them about 40 could be pathogenic for humans.

While they are not pathogenic for a healthy person, they may cause mild allergic reactions to significant and life-threatening infections in patients with compromised immune system. Most of these infections have been related to *A. fumigatus*, *A. flavus*, *A. terreus*

and *A. niger*. Early diagnosis of this disease is crucial for proper patient management.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
ASPG 4C Mix	PCR mix*	360 µl
<i>A.fum</i> Pos	Positive Control for <i>A. fumigatus</i>	100 µl
<i>A.flu</i> Pos	Positive Control for <i>A. flavus</i>	100 µl
<i>A.ter</i> Pos	Positive Control for <i>A. terreus</i>	100 µl
<i>A.nig</i> Pos	Positive Control for <i>A. niger</i>	100 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for in vitro diagnostics.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Proper samples to test for Aspergillus include whole blood, which should be collected in sterile tubes. We recommend EDTA or citrate as an anticoagulant. Samples such as plasma, cerebrospinal fluid (CSF) and bronchoalveolar lavage (BAL) can also be used. Samples can be stored at 2-8°C for a few hours or at -20°C or lower for up to a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To evaluate the possibility of DNA extraction failure and PCR inhibition and prevent false negative results, the Aspergillus RQ kit contains an internal control (IC). This IC can be used during extraction process or simply added to ASPG 4C Mix. To monitor both DNA extraction and PCR reaction, internal control should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of internal control is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of internal control should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to ASPG 4C Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC reagent should be added to each reaction. For example, for 10 PCR reactions, 10ul of the internal control should be added to 150ul of ASPG 4C Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the Internal control should generate a CT of 27-32 in the Yellow/VIC Channel.

14. DNA Isolation

DNA isolation can be performed using different fungal DNA extraction kits. For blood sample we recommend using:

- QIAamp® UCP PurePathogen Blood Kit (Cat. no. 50112, Qiagen GmbH, Hilden, Germany)

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus four for Positive control and one for the negative control.

If IC is introduced during the extraction process, pipette 15µl of [ASPG 4C Mix](#) to each PCR tube.

If IC is added to the ASPG 4C Mix, add 15ul of the [prepared mix](#) (as described in section 12) to each PCR tube.

Then add 10ul of isolated DNA, [Pos Ctrl](#), or water to each tube.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software

Aspergillus RQ4 kit is designed to work with Rotor-Gene 4-plex and MIC.


17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Aspergillus RQ4 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Aspergillus RQ4 0.1 is for strip tubes and Aspergillus RQ4 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain ASPG 4C Mix).

Aspergillus RQ4 (V1.1)

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for Green, Yellow, Orange and Red dyes.

19. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Analyze the data according to the manufacturer's recommendations.

Signal in the **Green** channels is due to **A. niger** or **A. terreus**. A signal in the **Yellow channel** is due to **IC**. Signal in the **Orange** channels is due to **A. flavus** and Signal in the **Red** channels is due to **A. fumigatus**.

Briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green. In the pop-up for Automatic Threshold, increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK or simply put the threshold on 0.05 for the Green channel.

Repeat the above for Cycling A. Yellow, Orange and Red but, cancel the Automatic Threshold and manually put the threshold on 0.1. Figures 1, 2, 3 and 4 represent typical graphs for the Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene


Fig 3. Typical Control graph in Orange channel for Rotor-Gene


Fig 4. Typical Control graph in Red channel for Rotor-Gene

Consider the following points when analyzing:

- A sample is **Positive** for *A. niger* or *A. terreus* if it is positive in the Green channel with a CT of less than 40.
- A sample is **Positive** for *A. flavus* if it is positive in the Orange channel with a CT of less than 40.
- A sample is **Positive** for *A. fumigatus* if it is positive in the Red channel with a CT of less than 40.
- A sample is **Negative** if it is negative in the Green, Orange and Red channels while it is positive in the Yellow channel with a CT of 27-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green, Orange, Red and Yellow channels.

The interpretation of results is summarized in the following Table.

Green	Yellow	Orange	Red	Result
+	+/-	-	-	Pos for <i>A. niger</i> or <i>A. terreus</i>
-	+/-	+	-	Pos for <i>A. flavus</i>
-	+/-	-	+	Pos for <i>A. fumigatus</i>
-	+	-	-	Negative
-	-	-	-	Inconclusive

20. Sensitivity

The analytical detection limit of the kit was assessed with cloned target and showed a limit of detection equal 2 copies/μl for *A. terreus* and *A. niger*, 7 copies/μl for *A. flavus* and 9 copies/μl *A. fumigatus*.

21. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

22. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

23. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124



Email: info@novingene.com

Website: www.novingene.com

24. References

- Fosses Vuong, M., Hollingshead, C.M. and Waymack, J.R., 2023. Aspergillosis.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212
- Samson, R.A. and Varga, J., 2009. What is a species in Aspergillus? Medical Mycology, 47(s1), pp.S13–S20.
- Van, H., Geert Cauwenbergh and Donald., 2013. Aspergillus and Aspergillosis. Springer Science & Business Media.

25. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com