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BRAF RQ Kit Manual

Autumn 2025, Version 1.2

For Real-Time PCR Detection of BRAF Mutations
For use with Rotor-Gene and StepOne
For Research Use Only

 24 (Cat# BRAFRQ24)

 48 (Cat# BRAFRQ48)

 NG-WI-ASL-40-102

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.

Table of Contents

1. Introduction 1

2. Intended Use 1

3. Background Information 1

4. Test Principle 1

5. Kit Contents 2

6. Packaging models..... 2

7. Storage and Stability 2

8. Product Use Limitations 2

9. Additionally Required Items 3

10. General Precautions 3

11. Specimen 4

12. DNA Extraction 4

13. Protocol: DNA Sample Assessment..... 4

14. Devices and software..... 5

15. Programming Real-time PCR..... 5

16. Analysis: Sample Assessment..... 7

17. Protocol: Detection of BRAF mutations..... 9

18. BRAF Mutation Detection Analysis 10

19. Analytical Sensitivity	17
20. Disposal Method	18
21. Technical Support.....	18
22. Contact Information.....	18
23. References	18
24. Symbols	19

1. Introduction

BRAF RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting BRAF V600 mutations in human genomic DNA. All required reagents are included in the PCR Mix provided in the kit. Control Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC).

This kit is intended for Research Use Only!

2. Intended Use

BRAF RQ kit is intended for detecting BRAF V600 mutations in human genomic DNA. Most of the BRAF mutations are located on codon 600 and constitute V600E, V600Ec, V600D, V600K, V600R. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene or StepOne machines.

3. Background Information

BRAF oncogene is among the most frequently mutated kinases in human cancer. Mutations in codon V600 have been reported in different types of cancers including 40-50% of melanomas, 10-70% of thyroid carcinomas, 10% of colorectal cancers and 3-5% of non-small cell lung cancers (NSCLC).

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-

amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and following reagents:

Label	Content	Quantity	Volume
BRAF Ctrl Mix	PCR Mix for quality control	2	480 µl
V600E Mix	PCR Mix to check V600E mutation	1	480 µl
V600Ec Mix	PCR Mix to check V600Ec mutation	1	480 µl
V600D Mix	PCR Mix to check V600D mutation	1	480 µl
V600K Mix	PCR Mix to check V600K mutation	1	480 µl
V600R Mix	PCR Mix to check V600R mutation	1	480 µl
BRAF Pos	1% Positive control	1	250 µl
BRAF Neg	Negative control	1	250 µl
Water	PCR Grade Water	1	200 µl

6. Packaging models

The kit is available in 24 and 48 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.

- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipment/items
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen

Sections of fixed tissue (paraffin block) or unfixed tissue can be examined. Note that tumor tissue is heterogeneous with uneven distribution of mutant cells; therefore, different parts of tissue may produce different results.

12. DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that the minimum volume of DNA required for the test is 40 microliters. We recommend extracting 100 microliters or a higher volume of DNA, enough to repeat the test if necessary.

Note: Extracted sample should contain 10-50 ng/ul DNA.

13. Protocol: DNA Sample Assessment

Before examining a sample for BRAF mutations, the quality of DNA should be assessed. If the results are within the desired range, then a second test for detecting BRAF mutations will be performed. To qualify DNA extraction, follow these steps:

Note! Check the volume of extracted DNA sample and make sure it is more than 40 microliters and preferably about 100 microliters. Volumes less than 40 microliters are not enough to proceed.

First, thaw the **BRAF Ctrl Mix** on ice completely, mix by inversion followed by a quick spin and store on crushed ice after. Place the required number of microtubes on a cold block including one for each sample, plus two for negative control and water.

Pipette 20µl of BRAF Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and Water to each tube.

Cap the tubes and inspect visually to ensure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

14. Devices and software


BRAF RQ kit is designed to work with Rotor-Gene and StepOne.

15. Programming Real-time PCR Machine

Rotor-Gene or StepOne could be setup using templates provided in the flash card, or accessible by kit QR code. For the Rotor-Gene, select “BRAF 0.1” is for strip tubes and “BRAF 0.2” is for 0.2ml tubes, Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the following image.

Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain BRAF Ctrl Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

You may also set up the machine as below.

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM/Green and VIC/Yellow dyes/channels. All mixes contain ROX with a final concentration of 300nM in the reaction.

Make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

16. Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. As the initial step, results of negative control will be evaluated.

Perform quantitative analysis for both **DNA sample (the FAM/Green channel)** and **internal control (the VIC/Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on "Cycling A. Green". Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Test is valid only if:

- 1) Water sample is **Negative** in the FAM/Green channel, and it is **Positive** in the VIC/Yellow channel with CT of 27-33.
- 2) Negative control is **Positive** in the VIC/Yellow channel with CT of 27-33.
- 3) Negative control is **Positive** in the FAM/Green channel with CT of 23-28.

Expected results for a valid test are summarized in the Table1.

Sample	VIC/Yellow	FAM/Green
NTC	Pos (CT: 27-33)	Neg
Negative Ctrl	Pos (CT: 27-33)	Pos (CT: 23-28)

Table1. Expected results for valid Sample Assessment.

If the above are met, then the test results are valid and may proceed to DNA sample assessment.

B) DNA sample Quality Analysis:

- **DNA sample is qualified only if** sample is positive in FAM/Green channel with CT of 22-30 and positive in VIC/Yellow channel with CT of 27-33. This sample can be further examined for BRAF mutations.

- **Sample should be diluted with water**, if a sample is positive in FAM/Green channel with CT of less than 22 and positive in VIC/Yellow channel with CT of 27-33. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the BRAF mutation test.

- **Result is invalid**, if a sample in the FAM/Green channel has CT above 30 and, in the VIC/Yellow channel CT of 27-33. Also, the result is invalid if a sample in the VIC/Yellow channel has CT above 33 and in the FAM/Green channel has CT of 22-30. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 2.

	FAM/Green	VIC/Yellow	Result
1	+ CT: 22-30	+ CT: 27-33	Valid
2	+ CT<22	+ CT: 27-33	Sample dilution
3	+ CT>30	+ CT: 27-33	Invalid

4	+	+	Invalid
	CT: 22-30	CT>33	

Table 2. DNA Quality Analysis

17. Protocol: Detection of BRAF mutations

To detect BRAF mutations, each sample should be examined with six mixes, one for DNA quantitation and one for each of five mutations. Therefore, to examine only one sample, 24 microtubes are required; six for the sample, 6 for the Positive control, six for the Negative control and six for NTC. Respectively for each extra sample, 6 microtubes will be added. Therefore, for two samples, 30 microtubes and for three samples 36 microtubes are required. Figure1 shows tube setup for three samples.

To start, Place required number of tubes on cold block organized in series of six each.

Pipette 20µl of BRAF Ctrl Mix to each tube in the first series.

Pipette 20µl of V600E Mix to each tube in the second series.

Pipette 20µl of V600Ec Mix to each tube in the third series.

Pipette 20µl of V600D Mix to each tube in the fourth series.

Pipette 20µl of V600K Mix to each tube in the fifth series.

Pipette 20µl of V600R Mix to each tube in the sixth series.

Continue by adding 5µl of extracted DNA, Positive control, Negative Control, and water to each tube. Consider the first and second tube in each series for positive control and Negative control and the third for water/NTC. The next tubes would be for samples.

Tubes setup for mixes and samples is shown in Fig 1.

Cap the tubes and visually inspect to ensure all are capped securely. Place tubes in the machine and program it according to the section 15.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.


Fig1. Tubes setup for Mixes and samples.

18. BRAF Mutation Detection Analysis

Before analyzing results for BRAF mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be examined. If the test is valid, then may proceed to BRAF mutation analysis.

Briefly, click on Analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase,

the sample is considered Negative, and CT, if present, is not reliable.

A) Test is valid only if:

- 1) NTC is Negative in FAM/Green channel with all Mixes, and
- 2) Positive Ctrl, Negative Ctrl and NTC are Positive in the VIC/Yellow channel with each Mix with a CT of 27-33, and
- 3) Positive and Negative controls are Positive in the FAM/Green channel with Control Mix with a CT of 22-27, and
- 4) Positive control is Positive in the FAM/Green channel **with each of specific Mixes** with a CT of 30-35, and
- 5) Sample is Positive in the VIC/Yellow channel with All mixes with CT of 27-33.
- 6) Sample is Positive in the FAM/Green channel with Control Mix with a CT of 22-30 and

Above steps for Test validation are summarized in Figure 2.

If all above are met, then the Test is valid, and results can be analyzed further for BRAF mutations.


Fig 2. Test validation flowchart

B) Data analysis for BRAF mutations

- 1) Select the patient's samples, which are positive in the FAM/Green channel with BRAF specific mixes with a CT of 20-40, and document the CTs in Table 3.
- 2) Calculate ΔCT for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the ΔCT in the Table 3 too.
- 4) If the ΔCT is within the valid range as mentioned in the Table 3, the sample has the BRAF mutation.
- 5) If a sample is positive for two or more BRAF mutations, it is considered positive for the mutation with the lowest ΔCT and negative for the other mutations.

Above steps for detection of BRAF mutations are summarized in Figure 3.

BRAF Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
Ctrl Mix	22-30		-		
600E Mix	20-40		≤ 7.3		
600Ec Mix	20-40		≤ 10.6		
600D Mix	20-40		≤ 9.6		
V600K Mix	20-40		≤ 9		
V600R Mix	20-40		≤ 8.3		

Table 3. Valid CT and ΔCT ranges for BRAF Mixes.

As an example, If CT of a sample is 24.7 with Control Mix, 35.2 with V600E Mix, 30.5 with V600Ec Mix, 40.4 with V600D Mix, 42.1 with V600K Mix and 43.5 for V600R mix, then ΔCT is 10.5 for V600E (35.2 - 24.7), 5.8 for V600Ec (30.5-24.7), 15.7 for V600D

and 17.4 for V600K, 18.8 for V600R (Table. 4). According to table 3, delta CT for V600Ec is in valid range, and patient is Positive for V600Ec.

BRAF Mix	Valid CT	Sample CT	Valid Δ CT	Sample Δ CT	Result
Ctrl Mix	22-30	24.7	-		Valid
V600E Mix	20-40	35.2	≤ 7.3	10.5	Neg
V600Ec Mix	20-40	30.5	≤ 10.6	5.8	Positive
V600D Mix	20-40	40.4	≤ 9.6	15.7	Neg
V600K Mix	20-40	42.1	≤ 9	17.4	Neg
V600R Mix	20-40	43.5	≤ 8.3	18.8	Neg

Table 4. BRAF mutation analysis for a specific sample.

Note that a normal sample may cross-react with some of the mutation specific mixes. However, the calculated Δ CT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in the FAM/Green channel for BRAF mutation mixes, will only be considered positive only if Δ CT is within the valid range.

Note: if sample is positive with V600E Mix, it could be positive for either V600E and V600K Mutation.


Fig 3. Sample analysis flowchart for BRAF mutation detection

Results can also be interpreted as below:

- Sample is **Negative** for BRAF mutations, if it is negative in the FAM/Green channel with all mixes or has CT above 40.
- Sample is **Negative** for BRAF mutations if it is positive in the FAM/Green channel with CT of 20-40 for one or more BRAF mixes, but ΔCT is higher than cutoff range (Table 3).

- Sample is **Positive** for BRAF mutations, if sample has a valid CT and Δ CT is lower than cutoff range (Table 3) in the FAM/Green channel for one of the five BRAF mixes. Briefly,
 - It is Positive for V600E or V600K if with V600E Mix has a $CT \leq 40$ and $\Delta CT \leq 7.3$.
 - It is Positive for V600Ec if with V600Ec Mix has a $CT \leq 40$ and $\Delta CT \leq 10.6$.
 - It is Positive for V600D if with V600D Mix has a $CT \leq 40$ and $\Delta CT \leq 9.6$.
 - It is Positive for V600K if with V600K Mix has a $CT \leq 40$ and $\Delta CT \leq 9$.
 - It is Positive for V600R if with V600R Mix has a $CT \leq 40$ and $\Delta CT \leq 8.3$.
- If a sample has a valid CT and Δ CT (Table 3) for more than one Mix, sample is Positive for the BRAF mutation, which has the lowest Δ CT.

Interpretation of results for BRAF mutation test are summarized in Table 5.

Note that, if a sample is Negative with this kit in terms of BRAF mutations, the following conditions should be considered:

- Sample is Negative only for BRAF mutations which was mentioned, and it could be Positive for other BRAF mutations.
- Sample is Positive for BRAF mutations in quantities below the sensitivity of the kit.

BRAF Mix	Sample CT	Sample Δ CT	Conclusion
600E Mix	>40	-	Neg for 600E mutation
	20-40	>7.3	Neg for 600E mutation
	20-40	≤ 7.3	Pos for 600E or 600K mutation
600Ec Mix	>40	-	Neg for 600Ec mutation
	20-40	>10.6	Neg for 600Ec mutation
	20-40	≤ 10.6	Pos for 600Ec mutation
600D Mix	>40	-	Neg for 600D mutation
	20-40	>9.6	Neg for 600D mutation
	20-40	≤ 9.6	Pos for 600D mutation
V600K Mix	>40	-	Neg for 600K mutation
	20-40	>9	Neg for 600K mutation
	20-40	≤ 9	Pos for 600K mutation
600R Mix	>40	-	Neg for 600R mutation
	20-40	>8.3	Neg for 600R mutation
	20-40	≤ 8.3	Pos for 600R mutation

Table 5. Interpretation of results for BRAF mutation.

19. Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated BRAF DNA that can be identified in the background of wild-type DNA and is mentioned in Table 6. The sensitivity depends on the DNA quantity in a sample. Maximum sensitivity is achieved when a sample with Control Mix has a CT of 22-27 in the FAM/Green channel; with CT of 28-30 in VIC/Yellow, the sensitivity will decrease for some mutations.

Reaction	Ctrl Mix CT: 22-27	Ctrl Mix CT: 28-30
V600E	1%	2%
V600Ec	0.5%	1%
V600D	0.5%	1%
V600K	1%	2%
V600R	2%	4%

Table 6. Sensitivity of the BRAF mutation detection assay.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124


Email: info@novingene.com

Website: www.novingene.com

23. References

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Sepulveda, A.R. and Lynch, J.P., 2013. Molecular Pathology of Neoplastic Gastrointestinal Diseases. Springer Science & Business Media.
- Sullivan, R.J. ed., 2014. BRAF targets in melanoma: biological mechanisms, resistance, and drug discovery (Vol. 82). Springer.
- Xing, M.B.R.A.F., 2005. BRAF mutation in thyroid cancer. Endocrine-related cancer, 12(2), pp.245-262.

24. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com