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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Corbel, M.J., 2006. World Health Organization, Food and Agriculture Organization of The United Nations and

International Office of Epizootics. Brucellosis in humans and animals. Geneva: World Health Organization.

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Moreno, E., Cloeckart, A. and Moriyón, I., 2002. Brucella evolution and taxonomy. Veterinary microbiology, 90(1-4), pp.209-227.
- Whatmore, A.M., 2009. Current understanding of the genetic diversity of Brucella, an expanding genus of zoonotic pathogens. Infection, Genetics and Evolution, 9(6), pp.1168-1184.
- Yu, W.L. and Nielsen, K., 2010. Review of Detection of Brucella sp. by Polymerase Chain Reaction. Croatian Medical Journal, 51(4), pp.306–313.

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Brucella RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of Brucella DNA
For Research Use Only

 24 (Cat# BrucellaRQ24)

 48 (Cat# BrucellaRQ48)

 96 (Cat# BrucellaRQ96)

 NG-WI-ASL-59-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction 3

2. Intended Use 3

3. Background Information 3

4. Test Principle 4

5. Kit Contents..... 4

6. Packaging models 4

7. Storage and Stability 4

8. Product Use Limitations 5

9. Additionally Required Materials 5

10. General Precautions..... 5

11. Specimen, Storage and Transport..... 6

12. Interfering Substances 6

13. Internal Control (IC)..... 6

14. DNA Isolation 7

15. PCR Protocol 7

16. Devices and software 8

17. Programming Rotor-Gene 8

18. Programming StepOne..... 9

19. Programming Other Machines..... 10

20. Data Analysis: Rotor-Gene	10
21. Data Analysis: StepOne	12
22. Sensitivity	14
23. Disposal Method.....	14
24. Technical Support	15
25. Contact Information.....	15
26. References	15
27. Symbols	16

1. Introduction

Brucella RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting Brucella DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

Brucella RQ kit is intended for detecting Brucella species including *Brucella melitensis*, *Brucella suis*, *Brucella abortus*, and *Brucella Canis*. The target genes are *BCSP31* and *IS711*. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Brucella is a gram-negative bacteria that causes Brucellosis, one of the most common zoonotic diseases. Brucellosis is most prevalent in areas where people live in close contact with infected animals. Annually, about 500,000 new cases of human Brucellosis are reported worldwide. Consumption of contaminated milk or milk products, direct contact with an infected animal, or inhalation of aerosols are the main routes of transmission. Twelve species of Brucella have been identified of which *Brucella melitensis*, *Brucella suis*, *Brucella abortus*, and *Brucella canis* have moderate to significant human pathogenicity. Initial symptoms include fever,

sweats, anorexia, headache, fatigue, and pain in muscles, joints and back.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
Brucella Mix	PCR mix*	360 µl
Pos Ctrl	Positive Control	100 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**

- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw the kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Proper samples to test for Brucella include whole blood, which should be collected in sterile tubes. We recommend EDTA or citrate as an anticoagulant. Samples such as plasma, serum, and bodily fluids can also be used. Samples can be stored at 2-8°C for 48 hours or at -20°C or lower for up to a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the Brucella RQ kit

contains an internal control (IC). This IC can be used during the extraction process or added directly to the Brucella Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to PCR Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 27-32 in the Yellow Channel for RotorGene and 28-34 in the VIC Channel for StepOne.

14. DNA Isolation

DNA isolation can be performed using different DNA extraction kits. For blood sample we recommend using:

- QIAamp DNA Blood Mini Kit (Cat# 51104, Qiagen, Germany)
- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold

block. Consider one tube for each sample, plus one for positive control and one for the negative control.

If the IC is introduced during the extraction process, pipette 15µl of the Brucella Mix to each PCR tube.

If the IC is added to the Brucella Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.

Then add 10ul of isolated DNA, Positive Control, or water to each tube.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software


Brucella RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Brucella template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Brucella 0.1 is for strip tubes and Brucella 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Brucella Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click on Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One positive control, one negative control, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The Brucella Mix contains ROX with the final concentration of 300nM in reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Analyze the data according to the Rotor-Gene manual. A Signal in the **Green** channel indicate **Brucella** and a signal in the **Yellow channel** indicate **IC**.

Briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green. In the pop-up for Automatic Threshold, increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK or simply set the threshold on 0.1 for the Green and Yellow channels.

Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 27-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of the Green and Yellow channels.

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for Brucella/FAM and IC/VIC channels at 0.1. Figures 3 and 4 represent typical graphs for the StepOne machine.


Fig 3. Typical Controls graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM channel with a CT of less than 40.
- A sample is **Negative** if it is negative in the FAM channel while it is positive in the VIC channel with a CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM and VIC channels.

The interpretation of the results is summarized in the following Table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

22. Sensitivity

The analytical detection limit of the kit was assessed with a cloned target and showed a limit of detection equal to 1 copy/ μ l.

23. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite

overnight and then discarded.

24. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

25. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com


Website: www.novingene.com

26. References

- Corbel, M.J. 2006, World Health Organization, Food and Agriculture Organization Of The United Nations and International Office Of Epizootics., 2006. Brucellosis in humans and animals. Geneva: World Health Organization.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212
- Moreno, E., Cloeckert, A. and Moriyón, I., 2002. Brucella evolution and taxonomy. Veterinary microbiology, 90(1-4), pp.209-227.
- Whatmore, A.M., 2009. Current understanding of the genetic diversity of Brucella, an expanding genus of zoonotic pathogens. Infection, Genetics and Evolution, 9(6), pp.1168-1184.

- Yu, W.L. and Nielsen, K., 2010. Review of Detection of Brucella sp. by Polymerase Chain Reaction. Croatian Medical Journal, 51(4), pp.306–313.

27. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com