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# EGFR RQ Kit Manual

Autumn 2025, Version 1.3

For Real-Time PCR Detection of EGFR Mutations  
For use with Rotor-Gene or StepOne  
For Research Use Only

 24 (Cat# EGFRRQ24)

 48 (Cat# EGFRRQ48)


 NG-WI-ASL-42-103

**RUO**

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

## Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	5
5. Kit Contents .....	5
6. Packaging models .....	5
7. Storage and Stability.....	6
8. Product Use Limitations .....	6
9. Additionally Required Items .....	6
10. General Precautions .....	7
11. Specimen.....	7
12. DNA Extraction .....	7
13. Protocol: DNA Sample Assessment.....	8
14. Programming Real-time PCR.....	8
15. Devices and software.....	8
16. Analysis: Sample Assessment.....	10
17. Protocol: Detection of EGFR mutations .....	12
18. EGFR Mutation Detection Analysis.....	14

19. Analytical Sensitivity .....	20
20. Disposal Method .....	21
21. Technical Support.....	21
22. Contact Information.....	21
23. References .....	22
24. Symbols .....	22

## 1. Introduction

EGFR RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detection of 34 somatic mutations in the EGFR gene. All required reagents are included in the PCR Mix provided in the kit. Control Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC).

This kit is intended for Research Use Only!

## 2. Intended Use

EGFR RQ kit is intended for the detection of 34 somatic mutations in EGFR exons 18, 19, 20 and 21. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene or StepOne machines.

## 3. Background Information

Epidermal Growth Factor Receptor (EGFR) is a tyrosine kinase receptor and is considered as an oncogene. EGFR is involved in the regulation of cellular proliferation, differentiation and survival. Mutations in EGFR exons 18, 19, 20 and 21 are associated with the development of different human cancers, especially Non-Small Cell Lung Carcinoma (NSCLC), and glioblastoma. Since, the choice of anti-EGFR therapies are highly dependent on the EGFR mutations, it is essential that patients are tested for these mutations.

This kit provides ready to use reagents for detection of 34 mutations. Detected mutations include three-point mutations in exon 18 (G719A, G719S and G719C without differentiating them), 24 deletions in exon 19 (without differentiation between them), three insertion (without differentiation between them) and two-

point mutation in exon 20 (S768I and T790M) and two point mutations in exon 21 (L858R, L861Q). For more information see table1.

Exon	Mutation	Cosmic* ID	Base change
18	G719A	6239	2156G>C
	G719S	6252	2155G>A
	G719C	6253	2155G>T
19	Deletions	6210	2240_2251 del 12
		6218	2239_2247 del 9
		6220	2238_2255 del 18
		6223	2235_2249 del 15
		6225	2236_2250 del 15
		6254	2239_2253 del 15
		6255	2239_2256 del 18
		12367	2237_2254 del 18
		12369	2240_2254 del 15
		12370	2240_2257 del 18
		12382	2239_2248TTAAGAGAAG>C
		12383	2239_2251>C
		12384	2237_2255>T
		12385	2235-2255delinsAAT
		12387	2239_2258>CA
		12403	2239-2256 delinsCCA
		12419	2238_2252>GCA
		12422	2238_2248>GC
		12678	2237_2251 del 15
		12728	2236_2253 del 18
		13550	2235-2248delinsAATTC
		13551	2235_2252>AAT
		13552	2235-2251delinsAATTC
		23571	Del 2238-2252
20	S768I	6241	2303G>T
	Insertion	12376	2307_2308insGCCAGCGTG
		12378	2310_2311insGGT
		12377	2319_2320insCAC

	T790M	6240	2369C>T
21	L858R	6224	2573T>G
	L861Q	6213	2582T>A

Table1: EGFR mutations detected by NovinGene kit.

\* COSMIC:Catalogue of somatic mutations in cancer: <http://cancer.sanger.ac.uk/>.

## 4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Lable	Content	Quantity	Volume
EGFR Ctrl Mix	PCR Mix for quality control	2	480 µl
G719X Mix	PCR Mix for G719A/G719SG719C	1	480 µl
19Del Mix	PCR Mix for exon 19 Deletions	1	480 µl
20Ins Mix	PCR Mix for exon 20 insertions	1	480 µl
S768I Mix	PCR Mix for S768I	1	480 µl
T790M Mix	PCR Mix for T790M	1	480 µl
L858R Mix	PCR Mix for L858R	1	480 µ
L861Q Mix	PCR Mix for L861Q	1	480 µ
EGFR Pos	Positive control	1	250 µl
EGFR Neg	Negative control	1	250 µl
Water	PCR Grade Water	1	200 µl

## 6. Packaging models

The kit is available in 24 and 48 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit and required equipment/items
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## 11. Specimen

Sections of formalin-fixed tissue (paraffin block) or unfixed tissue can be examined. Note that tumor tissue is heterogeneous with uneven distribution of mutant cells; therefore, different parts of tissue may produce different results.

## 12. DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPET DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 40 microliters of 10-50 ng/ul. We recommend extracting 100

microliters or higher volume of DNA, enough to repeat the test if necessary.

### 13. Protocol: DNA Sample Assessment

Before examining a sample for EGFR mutations, quality of DNA should be assessed. If the results are within the desired range, then the second test for detecting EGFR mutations will be performed. To qualify DNA extraction, follow steps below.

**Note! Check the volume of extracted DNA sample and make sure it is more than 40 microliters and preferably about 100 microliters. Volumes less than 40 microliters are not enough to proceed.**

First, thaw the **EGFR Ctrl Mix** on ice completely, mix by inversion followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for negative control and water.

**Pipette 20µl of EGFR Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and Water to each tube.**

Cap the tubes and inspect visually to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

### 14. Devices and software


EGFR RQ kit is designed to work with Rotor-Gene and StepOne.

### 15. Programming Real-time PCR

Rotor-Gene or StepOne could be setup using templates provided in the flash card, or accessible by kit QR code. For the Rotor-Gene, select “EGFR 0.1” is for strip tubes and “EGFR 0.2” is for 0.2ml tubes, ensuring that the locking ring is attached to the rotor.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain EGFR Ctrl Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine. You may also set up the machine as below.

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM/Green and VIC/Yellow dyes/channels.

All mixes contain ROX with final concentration of 300nM in the reaction.

Make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

## 16. Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. As the initial step, results of negative control will be evaluated.

Perform quantitative analysis for both **DNA sample (the FAM/Green channel)** and **internal control (the VIC/Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on "Cycling A. Green". Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold at 0.05.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

If the above are met, then the test results are valid and may proceed to DNA sample assessment.

### A) Test is valid only if:

- 1) Water sample is **Negative** in the FAM/Green channel, and it is **Positive** in the VIC/Yellow channel with CT of 27-33.
- 2) Negative control is **Positive** in the VIC/Yellow channel with CT of 27-33.

- 3) Negative control is **Positive** in the FAM/Green channel with CT of 23-28.

Expected results for a valid test are summarized in the Table2.

Sample	VIC/Yellow	FAM/Green
NTC	Pos (CT: 27-33)	Neg
Negative Ctrl	Pos (CT: 27-33)	Pos (CT: 23-28)

Table2. Expected results for valid Sample Assessment.

### B) DNA sample Quality Analysis:

- **DNA sample is qualified only if** sample is positive in FAM/Green channel with CT of 22-30 and positive in the VIC/Yellow channel with CT of 27-33. This sample can be further examined for EGFR mutations.

- **Sample should be diluted with water**, if a sample is positive in the FAM/Green channel with CT of less than 22 and positive in the VIC/Yellow channel with CT of 27-33. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the EGFR mutation test.

- **Result is invalid**, if a sample in the FAM/Green channel has CT above 30 and, in the VIC/Yellow channel CT of 27-33. Also, the result is invalid if a sample in the VIC/Yellow channel has CT above 33 and in the FAM/Green channel has CT of 22-30. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in the Table 3.

	FAM/Green	VIC/Yellow	Result
1	+ CT: 22-30	+ CT: 27-33	Valid
2	+ CT<22	+ CT: 27-33	Sample dilution
3	+ CT>30	+ CT: 27-33	Invalid
4	+ CT: 22-30	+ CT>33	Invalid

Table3. DNA Quality Analysis

### 17. Protocol: Detection of EGFR mutations

To detect EGFR mutations, each sample should be examined with eight mixes, one mix for DNA quantitation and seven mixes for detection of 34 mutations. Therefore, to examine only one sample, 32 microtubes are required: 8 for the Positive control, 8 for the Negative control and 8 for NTC and 8 for the sample. Respectively for each extra sample, 8 microtubes will be added. Therefore, for two samples, 40 microtubes and for three samples 48 microtubes are required. Figure1 shows tube setup for three samples.

To start, Place required number of tubes on a cold block organized in series of six each.

**Pipette 20µl of EGFR Ctrl Mix to each tube in the first series.**

**Pipette 20µl of G719X Mix to each tube in the second series.**

**Pipette 20µl of 19Del Mix to each tube in the third series.**

**Pipette 20µl of 20Ins Mix to each tube in the fourth series.**

**Pipette 20µl of S768I Mix to each tube in the fifth series.**

**Pipette 20µl of T790M Mix to each tube in the sixth series.**

**Pipette 20µl of L858R Mix to each tube in the sixth series.**

**Pipette 20µl of L861Q Mix to each tube in the sixth series.**  
**Continue by adding 5µl of Positive control, Negative Control, water and extracted DNA to each tube. Consider the first and second tube in each series for positive control and Negative control and the third for water/NTC. The next tubes would be for samples.**

Cap the tubes and visually inspect to ensure all are capped securely. Place tubes in the machine and program it according to the section 15.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*


Fig1. Tubes setup for Mixes and samples.

## 18. EGFR Mutation Detection Analysis

Before analyzing results for EGFR mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be examined. If the test is valid, then may proceed to EGFR mutation analysis.

Briefly, click on Analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold at 0.05.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**If all above are met, then the Test is valid, and results can be analyzed further for EGFR mutations.**

### **A) Test is valid only if:**

- 1) NTC is Negative in the FAM/Green channel with all Mixes, and
- 2) Positive Ctrl, Negative Ctrl and NTC are Positive in the VIC/Yellow channel with each Mix with a CT of 27-33, and
- 3) Positive control in the FAM/Green channel with Control Mix with a CT of 22-27, and
- 4) Positive control is Positive in the FAM/Green channel with all mixes with of 27-35,
- 5) Sample is Positive in the VIC/Yellow channel with Control Mix with a CT of 27-33,
- 6) Sample is Positive in the FAM/Green channel with Control Mix with a CT of 22-30,

Above steps for Test validation are summarized in Figure 2.


Fig 2. Test validation flowchart

## B) Data analysis for EGFR mutations

- 1) Select the patient's samples, which are positive in FAM/Green channel with EGFR specific mixes with a CT of 20-40, and document the CTs in the Table 4.
- 2) Calculate  $\Delta CT$  for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the  $\Delta CT$  in the Table 4 too.
- 4) If the  $\Delta CT$  is within the valid range as mentioned in the Table 4, the sample has the EGFR mutation and Positive.

Above steps for detection of EGFR mutations are summarized in Figure 3.

EGFR Mix	Valid CT	Sample CT	Valid $\Delta CT$	Sample $\Delta CT$	Result
EGFR Ctrl Mix	22-30		-		-
G719X Mix	20-40		$\leq 7.4$		
19Del Mix	20-40		$\leq 7.5$		
20Ins Mix	20-40		$\leq 6.9$		
S768I Mix	20-40		$\leq 10.3$		
T790M Mix	20-40		$\leq 8.2$		
L858R Mix	20-40		$\leq 10.2$		
L861Q Mix	20-40		$\leq 10$		

Table 4. Valid CT and  $\Delta CT$  ranges for EGFR Mixes.

As an example, If CT of a sample is 25.3 with Control Mix, 38.4 with G719X Mix, 33.2 with 19Del Mix, 36.8 with 20Ins Mix, 40 with S768I Mix, 36.9 for T790M mix, 40 with L858R Mix and 40 with L861Q Mix, then  $\Delta CT$  is 13.1 for G719X (38.4-25.3), 7.0 for 19Del

## EGFR RQ (V1.3)

(32.3-25.3), 11.5 for 20Ins and 14.7 for S768I, 11.6 for T790M, 14.7 for L858R and 14.7 for L861Q. (Table. 5). According to table 5, delta CT for 19Del is in valid range, and patient is Positive for 19Del.

EGFR Mix	Valid CT	Sample CT	Valid $\Delta$ CT	Sample $\Delta$ CT	Result
EGFR Ctrl Mix	22-30	25.3	-	-	Valid
G719X Mix	20-40	38.4	$\leq 7.4$	13.1	Neg
19Del Mix	20-40	32.3	$\leq 7.5$	7.0	Positive
20Ins Mix	20-40	36.8	$\leq 6.9$	11.5	Neg
S768I Mix	20-40	40	$\leq 10.3$	14.7	Neg
T790M Mix	20-40	36.9	$\leq 8.2$	11.6	Neg
L858R Mix	20-40	40	$\leq 10.2$	14.7	Neg
L861Q Mix	20-40	40	$\leq 10$	14.7	Neg

Table 5. EGFR mutation analysis for a specific sample.

*Note: above  $\Delta$ CTs should be validated each run by using Positive and Negative Control which provided in the kit.*

Note that a normal sample may cross-react with some of the mutation specific mixes. However, the calculated  $\Delta$ CT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in the FAM/Green channel for EGFR mutation mixes, will only be considered positive only if  $\Delta$ CT is within the valid range.


Fig 3. Sample analysis flowchart for EGFR mutation detection

Results can also be interpreted as below:

- Sample is **Negative** for EGFR mutations, if it is negative in the FAM/Green channel with all mixes or has CT above 40.
- Sample is **Negative** for EGFR mutations if it is positive in the FAM/Green channel with CT of 20-40 for one or more EGFR mixes, but  $\Delta CT$  is higher than cutoff range (Table 4).

- Sample is **Positive** for EGFR mutations, if sample has a valid CT and  $\Delta$ CT is lower than cutoff range (Table 4) in the FAM/Green channel for one of the seven EGFR mixes. Briefly,
  - It is Positive for G719X if with G719X Mix has a  $CT \leq 40$  and  $\Delta CT \leq 7.4$ . Note that, this kit doesn't differentiate between G719A, G719S and G719C mutations.
  - It is Positive for 19Del if with 19Del Mix has a  $CT \leq 40$  and  $\Delta CT \leq 7.5$ . Note that this kit detects 24 deletion mutations (Table1) but not distinguish between them.
  - It is Positive for 20Ins if with 20Ins Mix has a  $CT \leq 40$  and  $\Delta CT \leq 6.9$ . Note that this kit detects 3 insertion mutations (Table1) but not distinguish between them.
  - It is Positive for S768I if with S768I Mix has a  $CT \leq 40$  and  $\Delta CT \leq 10.3$ .
  - It is Positive for T790M if with T790M Mix has a  $CT \leq 40$  and  $\Delta CT \leq 8.2$ .
  - It is Positive for L858R if with L858R Mix has a  $CT \leq 40$  and  $\Delta CT \leq 10.2$ .
  - It is Positive for L861Q if with L861Q Mix has a  $CT \leq 40$  and  $\Delta CT \leq 10$ .

Interpretation of results for EGFR mutation test are summarized in Table 6.

**Note that, if a sample is Negative with this kit in terms of EGFR mutations, the following conditions should be considered:**

- Sample is Negative only for EGFR mutations which was mentioned in Table1, and it could be Positive for other EGFR mutations.
- Sample is Positive for EGFR mutations in quantities below the sensitivity of the kit.

EGFR Mix	Sample CT	Sample $\Delta$ CT	Conclusion
G719X Mix	>40	-	Neg for G719A/ G719S/ G719C mutations
	20-40	>7.4	Neg for G719A/ G719S/ G719C mutations
	20-40	$\leq 7.4$	Pos for G719A/ G719S/ G719C mutations
19Del Mix	>40	-	Neg for Exon 19 Deletions
	20-40	>7.5	Neg for Exon 19 Deletions
	20-40	$\leq 7.5$	Pos for Exon 19 Deletions
20Ins Mix	>40	-	Neg for Exon 20 Insertions
	20-40	>6.9	Neg for Exon 20 Insertions
	20-40	$\leq 6.9$	Pos for Exon 20 Insertions
S768I Mix	>40	-	Neg for S768I mutation
	20-40	>10.3	Neg for S768I mutation
	20-40	$\leq 10.3$	Pos for S768I mutation
T790M Mix	>40	-	Neg for T790M mutation
	20-40	>8.2	Neg for T790M mutation
	20-40	$\leq 8.2$	Pos for T790M mutation
L858R Mix	>40	-	Neg for L858R mutation
	20-40	>10.2	Neg for L858R mutation
	20-40	$\leq 10.2$	Pos for L858R mutation
L861Q Mix	>40	-	Neg for L861Q mutation
	20-40	>10	Neg for L861Q mutation
	20-40	$\leq 10$	Pos for L861Q mutation

Table 6. Interpretation of results for EGFR mutation.

## 19. Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated EGFR DNA that can be identified in the background of wild-type DNA. Sensitivity of the assay is determined by dilution of plasmids containing target sequence in human genomic DNA

and is mentioned in Table 7. The sensitivity depends on the DNA quality and quantity. Maximum sensitivity is achieved when a sample with Control Mix has a CT of 22-27 in the FAM/Green channel.

Reaction	Ctrl Mix CT: 22-27
G719X Mix	2%
19Del Mix	1%
20Ins Mix	2%
S768I Mix	1%
T790M Mix	8%
L858R Mix	1%
L861Q Mix	1%

Table 7. Sensitivity of the EGFR mutation detection assay.

## 20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 21. Technical Support

For technical support, contact us via  
Phone: +98 993-6223241  
Email: [info@novingene.com](mailto:info@novingene.com)

## 22. Contact Information

**NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124



Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)

## 23. References

- Hu, S. ed., 2022. Novel Sensitizing Agents for Therapeutic Anti-EGFR Antibodies (Vol. 1). Academic Press.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Rotow, J., Cheema, P., Pisapia, P. and Troncone, G. 2022. Fast Facts: EGFR Exon 20 Insertion Mutations in NSCLC. Karger Medical and Scientific Publishers.
- Zhang, Y.L., Yuan, J.Q., Wang, K.F., Fu, X.H., Han, X.R., Threapleton, D., Yang, Z.Y., Mao, C. and Tang, J.L., 2016. The prevalence of EGFR mutation in patients with non-small cell lung cancer: a systematic review and meta-analysis. Oncotarget, 7(48), p.78985.

## 24. Symbols

 <b>RUO</b>	Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b>	Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b>	Catalogue number	 <b>SN</b>	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**