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## ETV6-RUNX1 RQ (TEL-AML1)

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# **ETV6-RUNX1 RQ (TEL-AML1) Kit Manual**

Autumn 2025, Version 1.2

For Real-Time PCR Detection and Quantitation of ETV6-RUNX1  
(TEL-AML1) Transcripts  
For Research Use Only

 24 (Cat# ETV6-RUNX1RQ24)

 48 (Cat# ETV6-RUNX1RQ48)

 96 (Cat# ETV6-RUNX1RQ96)


 NG-WI-ASL-29-102

**RUO**

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

## Table of Contents

1. Introduction.....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle.....	3
5. Kit Contents.....	4
6. Packaging models .....	4
7. Storage and Stability .....	5
8. Product Use Limitations.....	5
9. Additionally Required Materials .....	5
10. General Precautions.....	6
11. Specimen, Storage and Transport.....	6
12. Interfering Substances.....	7
13. RNA Isolation .....	7
14. cDNA Synthesis.....	7
15. PCR Protocol.....	8
16. Devices and software .....	8
17. Programming RotorGene.....	8

18.	Programming of StepOne.....	10
19.	Programming Other Machines.....	10
20.	Data Analysis: RotorGene.....	10
21.	Data Analysis: StepOne .....	13
22.	ETV6-RUNX1% calculation .....	16
23.	Analytical Sensitivity.....	16
24.	Disposal Method.....	17
25.	Technical Support .....	17
26.	Contact Information.....	17
27.	References .....	17
28.	Symbols .....	18

## 1. Introduction

ETV6-RUNX1 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying ETV6-RUNX1 transcripts. All required reagents are included in the PCR Mix provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

**Important Note:** *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

## 2. Intended Use

ETV6-RUNX1 (TEL-AML1) RQ kit is intended for the Quantitation of ETV6-RUNX1 (TEL-AML1) transcripts. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## 3. Background Information

ETV6-RUNX1 (TEL-AML1) is an abnormality resulted from 21;12 translocation. RUNX1 gene, encodes a protein involved in transcriptional control of hematopoiesis. However, as a result of this translocation, it is repressed by ETV6-RUNX1 fusion protein. This alteration occurs in approximately 25% of childhood ALL diagnosed between the ages of 2-10 years.

ETV6-RUNX1 (TEL-AML1) RQ provides a ready-to-use system for detection and quantitation of ETV6-RUNX1 transcripts as well as calculation of ETV6-RUNX1 percentage.

## 4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR

facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual or a flash card and the following reagents:

Label	Content	Quantity
<b>ETV6-RNX1 Mix*</b>	Mix for ETV6-RUNX1	480 µl
<b>ABL Mix*</b>	Mix for ABL	480 µl
<b>ER1</b>	ETV6-RUNX1 Standard 1: 100,000 copies/µl	150 µl
<b>ER2</b>	ETV6-RUNX1 Standard 2: 10,000 copies/µl	150 µl
<b>ER3</b>	ETV6-RUNX1 Standard 3: 1,000 copies/µl	150 µl
<b>ER4</b>	ETV6-RUNX1 Standard 4: 100 copies/µl	150 µl
<b>ER5</b>	ETV6-RUNX1 Standard 5: 10 copies/µl	150 µl
<b>A1</b>	ABL Standard 1: 100,000 copies/µl	150 µl
<b>A2</b>	ABL Standard 2: 10,000 copies/µl	150 µl
<b>A3</b>	ABL Standard 3: 1,000 copies/µl	150 µl
<b>A4</b>	ABL Standard 4: 100 copies/µl	150 µl
<b>Water</b>	PCR Grade Water	200 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- **Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## 11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 48 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used. For optimum results, a sample should include 5 to 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

## 12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, and Mannheim, Germany).
- TRIzol isolation reagent (Cat. no. 15596026, Invitrogene/Thermo Fisher, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime/Thermo Fisher, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

## 14. cDNA Synthesis

2-4ug of total RNA is required and should be reverse transcribed to cDNA using random hexamers or specific primers. Different kits are available in the market for this purpose.

These specific primers are available upon request. If the CT observed for ABL using Random Hexamers is above 27, it is recommended to use specific primers for decreasing CT for ABL and ETV6-RUNX1 gene and increasing sensitivity of the test.

Dissolve RNA pellet in 20ul of the Primers solution and then add other reagents required for cDNA synthesis.

Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20ul of cDNA add 30ul of nuclease free water.

## 15. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

Each sample should be examined for both in ETV6-RUNX1 fusion gene and for ABL control gene. So, two sets of reactions are required. In ETV6-RUNX1 set, consider 1 tube for each sample as well as 6 tubes for the 5 standards (ER1 to ER5) and Negative control or NTC. In ABL set, consider 1 tube for each sample and 5 tubes for the 4 standards (A1 to A4) and Negative control or NTC. Place required number of tubes on cold block.

**Pipette 20µl of ETV6-RNX1 Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA or standard or control to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using RotorGene attach the locking ring.*

## 16. Devices and software

ETV6-RUNX1 RQ kit is designed to work with the Rotor-Gene, StepOne and MIC.

## 17. Programming RotorGene


*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the ETV6-RUNX1 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); ETV6-RUNX1 0.1 is for strip tubes and ETV6-RUNX1 0.2 is for 0.2ml tubes. Program starts.

## ETV6-RUNX1 RQ (V1.2)

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Select a tube number containing ETV6-RNX1 Mix for the Green channel and tube with ABL Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file.

Edit sample names on both ETV6-RUNX1 and ABL pages. Remember that tubes containing ETV6-RNX1 Mix should only be named in ETV6-RUNX1 page and tubes containing ABL Mix should only be named in ABL page.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been

entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

## 18. Programming of StepOne

Open the StepOne software (V 2.\*). On the Set Up menu, click on Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. One negative control, five standards for ETV6-RUNX1, four standards for ABL and a few samples are defined. You may change plate set up using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished click on Start Run and save the experiment. Instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60C for FAM and VIC dyes. Both of ETV6-RUNX1 Mix and ABL Mix contain ROX with final concentration of 300nM in reaction.

## 20. Data Analysis: RotorGene


Before analyzing results, make sure in the sample menu all, the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be

defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC", respectively.


Analyze the data according to RotorGene manual. Perform quantitative analysis for both **ETV6-RUNX1 (the Green channel)** and **ABL (the Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on "Cycling A. Green".

Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1. Figures 1 and 2 represent typical graphs for the RotorGene machine.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 1.** Typical ETV6-RUNX1 graph in Green channel for Rotor-Gene

**Fig 2.** Typical ABL graph in Green channel for Rotor-Gene


Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the Green/ETV6-RUNX1 and the Yellow/ABL channels with sigmoid graphs and CT of 20-40 for the Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the Green/ETV6-RUNX1 channel while it is positive in the Yellow/ABL channel with a sigmoid graph and CT of 20-27 with titer of above 20,000 copies/ul.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green/ETV6-RUNX1 and Yellow/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green/ETV6-RUNX1 channel while it is positive in the Yellow/ABL channel with sigmoid graph and CT of above 27. Improper extraction or low RNA input could cause that.


*Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used. CTs higher than 27 reduce the sensitivity of the test and may result in **false negative** reports.*

## 21. Data Analysis: StepOne

Analyze the data according to StepOne manual. Briefly, click on "Analyze" and set the threshold for both **ETV6-RUNX1/FAM** and **ABL/VIC** on 0.1. Figures 3 and 4 represent typical graphs for the StepOne machine

**Fig 3.** Typical ETV6-RUNX1 graph in FAM channel for StepOne

**Fig 4.** Typical ABL graph in VIC channel for StepOne

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the FAM/ETV6-RUNX1 and the VIC/ABL channels with sigmoid graphs and CT of 20-40 for the FAM and 20-30 for the VIC.
- A sample is **Negative** if it is negative in the FAM/ETV6-RUNX1 channel while it is positive in the VIC/ABL channel with a sigmoid graph and CT of 20-27 with titer of above 20,000 copies/ul.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/ETV6-RUNX1 and the VIC/ABL channels. Improper extraction or test setup could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM/ETV6-RUNX1 channel while it is positive in the VIC/ABL channel with sigmoid graph and CT of above 27. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used. CTs higher than 27 reduce the sensitivity of the test and may result in **false negative** reports.*

## 22. ETV6-RUNX1% calculation

To assess the response to therapy, ETV6-RUNX1% value for each patient can be calculated. This kit uses NCN method for this purpose (*Beillard E. 2003, Leukemia 17:2474*). In this method ETV6-RUNX1% value is the ETV6-RUNX1 transcripts (titer) normalized by the ABL transcripts (titer) and then multiplied by 100%.

Usually, titer of ABL gene transcripts is higher than for ETV6-RUNX1 fusion gene. Therefore, results of above calculation would be a number less than 100%. This value may fall even below 0.001% in optimum cases. However, at diagnosis and in case of relapse, ETV6-RUNX1 transcripts may surpass ABL gene transcripts. As a result, ratio of ETV6-RUNX1/ABL would be above 100%. Such results have been reported previously too. For example, please note table 12 of the paper published by J. Gabert *et al* (2003, *Leukemia 17:2318*) for B-ALL patients. He reported a range of 20% to 1520% with average of 471% for peripheral blood and range of 46% to 3270% with average of 445% for bone marrow samples. Therefore, a ratio above 100% is not unexpected and simply means higher transcripts of the fusion gene compared to the control gene.

## 23. Analytical Sensitivity

Analytical sensitivity of the kit has been determined by examining dilution series of cloned cDNA and was estimated about 2 copeis/ul.

## 24. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 25. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: [info@novingene.com](mailto:info@novingene.com)

## 26. Contact Information

### NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124




Website: [www.novingene.com](http://www.novingene.com)

## 27. References

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Miles, G., Rajagopal, G., Strair, R. and Sabaawy, H.E., 2010. Identification and Modeling of TEL-AML1 (ETV6-RUNX1) Molecular Signature In Acute Lymphoblastic Leukemia. Blood, 116(21), p.2500.
- Zapotocky, M., Starkova, J., Hrusak, O., Horak, J., Trka, J., Zuna, J., Madzo, J., Krejci, O., Zemanova, Z., Kalinova, M. and Muzikova, K., 2011. TEL/AML1 (ETV6/RUNX1. BLOOD, 117(1).

## 28. Symbols

 <b>RUO</b>	Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b>	Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b>	Catalogue number	 <b>SN</b>	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**