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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Chen, B.S., Lee, H.C., Lee, K.M., Gong, Y.N. and Shih, S.R., 2020. Enterovirus and encephalitis. *Frontiers in microbiology*, 11, p.261.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.
- Sawyer, M.H., 2002, January. Enterovirus infections: diagnosis and treatment. In *Seminars in pediatric infectious diseases* (Vol. 13, No. 1, pp. 40-47). WB Saunders.
- Solomon, T., Lewthwaite, P., Perera, D., Cardoso, M.J., McMinn, P. and Ooi, M.H., 2010. Virology, epidemiology, pathogenesis, and control of enterovirus 71. *The Lancet infectious diseases*, 10(11), pp.778-790.

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Enterovirus RQ Kit Manual

Autumn 2025, Version 1.0

For Real-Time RT-PCR Detection of Enterovirus RNA
For Research Use Only

 24 (Cat# ENTVRQ24)

 48 (Cat# ENTVRQ48)

 96 (Cat# ENTVRQ96)

 NG-WI-ASL-70-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Internal Control (IC)	6
13. RNA Isolation	7
14. RT-PCR Protocol	7
15. Devices and software	8
16. Programming Rotor-Gene	8
17. Programming StepOne	9
18. Programming Other Machines	10

19. Data Analysis: Rotor-Gene	10
20. Data Analysis: StepOne	12
21. Sensitivity	14
22. Disposal Method	14
23. Technical Support	14
24. Contact Information	14
25. References	15
26. Symbols	15

1. Introduction

Enterovirus RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting Enterovirus RNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

The Enterovirus RQ kit is intended for detecting RNA of Enterovirus. Detection is achieved using Real-Time RT-PCR system with StepOne, Rotor-Gene and MIC machines.

3. Background Information

Enteroviruses are non-enveloped viruses with single-stranded RNA genomes. The genus Enterovirus is a member of the Picornaviridae family. Transmission primarily occurs via the fecal–oral route, although respiratory spread is also reported. While many infections remain asymptomatic or mild, enteroviruses can cause a broad clinical spectrum, including respiratory infection, aseptic meningitis, and acute flaccid paralysis. Because cell culture is time-consuming and has limited sensitivity, real-time reverse transcription PCR is recommended for rapid and accurate detection.

Please note that Enteroviruses and Rhinoviruses both belong to the Picornaviridae family and sequence similarity may lead to non-specific reactions. Therefore, the Enterovirus RQ kit may, in

rare cases (particularly with respiratory samples) show cross-reactivity with certain rhinoviruses.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual a flash card with templates and the following reagents:

Label	Content	Quantity
Entero Mix	PCR mix*	360 µl
Entero Pos Ctrl	Enterovirus Positive Control	150 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes,

and c) Reaction preparation area for addition of extracted RNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

The appropriate clinical sample for the detection of Enterovirus varies depending on the site of infection and patient's symptoms. Commonly used specimens include cerebrospinal fluid (CSF), respiratory sample and stool.

Note: With respiratory samples such as nasopharyngeal swabs, in rare cases the assay may show cross-reactivity with certain rhinoviruses.

Samples should be shipped at +4°C and can be stored at +4°C for 48 hours or aliquoted and stored at -20°C for up to a few weeks.

12. Internal Control (IC)

To assess the possibility of RNA extraction failure and PCR inhibition, and to prevent false negative results, the Enterovirus RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the Entero Mix. To monitor both RNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer.

For instance, if the extracted RNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to Entero Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC reagent should be added to each reaction. For example, for 10 PCR reactions, 15ul of the internal control should be added to 150ul of Entero Mix before it is added to the tubes.

In a successful RNA extraction and PCR test, the Internal control should generate a CT of 29-34 in the Yellow Channel on Rotor-Gene.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat# 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52906, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany)

14. RT-PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold

block. Consider one tube for each sample plus one for positive sample and one for water.

If the IC is introduced during the extraction process, pipette 15µl of the [Entero Mix](#) to each PCR tube.

If the IC is added to the Entero Mix, add 15ul of the [prepared mix](#) (as described in section 13) to each PCR tube.

Aliquot 15µl of [Entero Mix](#) directly to each PCR tube. Then add 10ul of extracted RNA, [Positive Control](#) or water.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

15. Devices and software


Enterovirus RQ kit is designed to work with StepOne, Rotor-Gene and MIC.

16. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Enterovirus template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Enterovirus 0.1 is for strip tubes and Enterovirus 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for both channels (note that Tube number 1 should contain Entero Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

17. Programming of the StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. Positive and negative control, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment to the desired location. The instrument will start shortly.

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Entro Mix contain ROX with a final concentration of 300nM in reaction.

19. Data Analysis: Rotor-Gene

To analyze data briefly, click on analysis menu and then under Quantitation tab double-click on cycling A. Green. Manually put threshold at 0.05. Repeat the above for the Yellow channel and put threshold at 0.1. Figures 1 and 2 represent typical graphs for the Rotor-Gene. To interpret the results, please note that:

Signal in the **Green** channel is due to **Enterovirus**, and a signal in the **Yellow** channel is due to **IC**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical Control graph in Green channel for Rotor-Gene


Fig 2. Typical Internal Control graph in Yellow channel for Rotor-Gene

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green/Enterovirus channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the Green/Enterovirus channel while it is positive in the Yellow/IC channel with a sigmoid graph and CT of 29-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green/ Enterovirus and Yellow/IC channels.

The interpretation of results is summarized in the bellow Table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

20. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for Enterovirus/FAM and IC/VIC at 0.1. Figures 3 and 4 represent typical graphs for the StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical Controls graph in FAM channel for StepOne


Fig 4. Typical Internal Control graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/ Enterovirus channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the FAM/ Enterovirus channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 27-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/ Enterovirus and VIC/IC channels.

The interpretation of results is summarized in the table on page 12.

21. Sensitivity

The analytical detection limit of the kit was assessed with cloned target, of Enterovirus genome and showed a limit of detection equal to 21 copies/ μ l.

22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

23. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

24. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124




Email: info@novingene.com

Website: www.novingene.com

25. References

- Chen, B.S., Lee, H.C., Lee, K.M., Gong, Y.N. and Shih, S.R., 2020. Enterovirus and encephalitis. *Frontiers in microbiology*, 11, p.261.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.
- Sawyer, M.H., 2002, January. Enterovirus infections: diagnosis and treatment. In *Seminars in pediatric infectious diseases* (Vol. 13, No. 1, pp. 40-47). WB Saunders.
- Solomon, T., Lewthwaite, P., Perera, D., Cardosa, M.J., McMinn, P. and Ooi, M.H., 2010. Virology, epidemiology, pathogenesis, and control of enterovirus 71. *The Lancet infectious diseases*, 10(11), pp.778-790.

26. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com