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ایمیل: info@novingene.com

وبسایت: www.novingene.ir

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- Hotoleanu, C., 2017. Genetic risk factors in venous thromboembolism. Thrombosis and Embolism: from Research to Clinical Practice: Volume 1, pp.253-272.
- Kaushansky, K., 2021. Williams hematology. New York: Mcgraw Hill.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3., 2004: 190-212.

- Varga, E.A. and Moll, S., 2004. Prothrombin 20210 mutation (factor II mutation). Circulation, 110(3), pp.e15-e18.

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# Factor II RQ Kit Manual

Autumn 2025, Version 3.1

For Real-Time PCR Detection of Factor II (G20210A) mutation  
For Research Use Only

 24 (Cat# F2RQ24)

 48 (Cat# F2RQ48)

 96 (Cat# F2RQ96)

 NG-WI-ASL-16-301

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

|   |   |
|---|---|
| 1. Introduction .....                     | 3 |
| 2. Intended Use .....                     | 3 |
| 3. Background Information .....           | 3 |
| 4. Test Principle .....                   | 3 |
| 5. Kit Contents .....                     | 4 |
| 6. Packaging models .....                 | 4 |
| 7. Storage and Stability .....            | 4 |
| 8. Product Use Limitations .....          | 4 |
| 9. Additionally Required Materials .....  | 5 |
| 10. General Precautions .....             | 5 |
| 11. Specimen, Storage and Transport ..... | 6 |
| 12. Interfering Substances .....          | 6 |
| 13. Internal Control .....                | 6 |
| 14. DNA Isolation .....                   | 6 |
| 15. PCR Protocol .....                    | 7 |
| 16. Devices and software .....            | 7 |
| 17. Programming Rotor-Gene .....          | 7 |

|                                     |    |
|-------------------------------------|----|
| 18. Programming StepOne.....        | 8  |
| 19. Programming Other Machines..... | 9  |
| 20. Data Analysis: Rotor-Gene.....  | 9  |
| 21. Data Analysis: StepOne .....    | 11 |
| 22. Disposal Method.....            | 13 |
| 23. Technical Support .....         | 13 |
| 24. Contact Information .....       | 14 |
| 25. References.....                 | 14 |
| 26. Symbols .....                   | 14 |

## **1. Introduction**

Factor II RQ kit provides a ready-to-use Real-Time PCR test for the detection of F2 using Real-Time PCR (based on polymerase chain reaction). All required reagents are included in the PCR Mix provided in the kit.

This kit is for Research Use Only!

## **2. Intended Use**

Factor II RQ kit is intended for detection of prothrombin G20210A (Factor II) mutation in human DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## **3. Background Information**

Prothrombin or factor II is the precursor for thrombin in the coagulation cascade. G20210A mutation in prothrombin gene is the second most prevalent genetic risk factor for venous thromboembolism (VTE) after Factor V Leiden. This mutation increases the risk of venous thrombosis about 2 to 7 times. It has also been linked to increased risk of pregnancy loss. Prevalence of this mutation ranges between 0.5 to 5% depending on the race and ethnic background.

The diagnosis of factor II G20210A mutation requires DNA analysis of factor II gene.

## **4. Test Principle**

The target sequence is detected using PCR, where primers specific to the target sequence amplify it. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-

amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

| Label         | Content                       | Quantity    |
|---------------|-------------------------------|-------------|
| F2 RQ Mix     | PCR Mix*                      | 480 $\mu$ l |
| Thrmb MM Ctrl | Homozygous Positive Control   | 100 $\mu$ l |
| Thrmb WM Ctrl | Heterozygous Positive Control | 100 $\mu$ l |
| Thrmb WW Ctrl | Negative Control              | 100 $\mu$ l |
| Water         | PCR Grade Water               | 200 $\mu$ l |

\*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.

- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

## **11. Specimen, Storage and Transport**

Whole blood (0.5ml) is the preferred sample. Whole blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for few days). For longer terms, sample should be aliquoted and stored at -20°C which is stable for up to a few months.

## **12. Interfering Substances**

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## **13. Internal Control**

Since each person carries wild type, mutant or both alleles, it serves as both the target and the Internal Control for the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, then the test should be repeated.

## **14. DNA Isolation**

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

## 15. PCR Protocol

Thaw the reagents on crushed ice completely, followed by a brief mixing and a quick spin. Place the required number of microtubes on the cold block. Consider one tube for each sample plus one for each control and one for the NTC.

**Pipette 20µl of F2 RQ Mix directly to each tube followed by adding 5ul of controls or sample DNA.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring too.*

*Please note that genotype of the samples can be called only if all the four controls provided with kit are used in each test!*

## 16. Devices and software


F2 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

## 17. Programming Rotor-Gene

*- Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the F2 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); F2 0.1 is for strip tubes and F2 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file. Edit sample names.

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on "Plate Setup". Controls and a few samples are defined. You may change the plate setup using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment on the desired location. Instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

| Step | Temperature and time | Cycles |
|------|----------------------|--------|
| 1    | <b>95°C x 3 min</b>  | 1      |
| 2    | <b>95°C x 15 sec</b> | 45     |
|      | <b>58°C x 40 sec</b> |        |

Fluorescence should be collected at 58°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

## 20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to Rotor-Gene manual. Briefly, click on Analysis menu select "Other" and then "Scatter Graph Analysis". Using Ctrl button, mark both Green and Yellow channels and then click on "Show".

Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test!


Also, make sure that in the Setup section, the alleles MM, WM, and WW are defined.

The "Scatter Graph Window" is now shown in the lower left window. Each dot represents one sample. The vertical axis shows the Green fluorescence and horizontal axis is for the Yellow fluorescence. Note that M allele is detected in the Green channel and W allele in the Yellow channel.

Therefore, samples are located in four regions of the graph. MM


samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.

To determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 1).

**Fig 1.** Defining genotypes of different controls of the test on Rotor-Gene

Now turn on all the samples and the results are displayed in the tab “Scatter Analysis Results” (Figure 2).

**Fig 2.** Typical scatter graph for controls and samples in Rotor-Gene

**Note!** If location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also, if a sample is located in between regions, it should be re-examined!

## 21. Data Analysis: StepOne


Analyze data according to StepOne manual. Briefly, click on Analyze and select “Allelic Discrimination” tab. Genotype of the samples are automatically determined by the software based on comparison of their fluorescence with the controls.

Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test!

please ensure that the alleles MM, WM, and WW are defined in the Setup section.

Figure 3 represents the typical graph. Each dot represents one sample. All samples are color coded based on their genotype as blue for MM, Red for WW, Green for WM, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "×". Note that, the vertical axis shows the FAM fluorescence and horizontal axis is for the VIC fluorescence. Accordingly, M allele is detected in FAM channel and W allele in VIC channel. Therefore, samples are located in four regions of the graph. MM samples with high FAM and low VIC fluorescence gather in upper left; WW samples with low FAM fluorescence and high VIC gather in lower right; WM samples with high FAM and high VIC fluorescence gather in upper right; and finally, NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence. Genotypes of samples are automatically determined by the software based on comparison with the controls.

**Fig 3.** Typical graph for controls and samples on StepOne

**Note!** If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid and test should be repeated!

## 22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 23. Technical Support

For technical support, contact us via  
Phone +98-9936223241  
Email: [info@novingene.com](mailto:info@novingene.com)

## 24. Contact Information

### NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124




Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)

## 25. References

- Hotoleanu, C., 2017. Genetic risk factors in venous thromboembolism. Thrombosis and Embolism: from Research to Clinical Practice: Volume 1, pp.253-272.
- Kaushansky, K., 2021. Williams hematology. New York: Mcgraw Hill.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3., 2004: 190-212.
- Varga, E.A. and Moll, S., 2004. Prothrombin 20210 mutation (factor II mutation). Circulation, 110(3), pp.e15-e18.

## 26. Symbols

|  |                   |   |   |
|--|-------------------|---|---|
|  <b>RUO</b> | Research use only |  <b>Manufacturer</b>                           |  <b>Consult instructions for use</b> |
|  <b>LOT</b> | Lot number        |  <b>Content sufficient for &lt;n&gt; tests</b> |  <b>Use-by date</b>                  |
|  <b>REF</b> | Catalogue number  |  <b>SN</b>                                     |  <b>Temperature limit</b>            |

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**