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ایمیل: info@novingene.com

وبسایت: www.novingene.ir

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- Hotoleanu, C., 2017. Genetic risk factors in venous thromboembolism. *Thrombosis and Embolism: from Research to Clinical Practice: Volume 1*, pp.253-272.
- Kaushansky, K., 2021. *Williams hematology*. New York: Mcgraw Hill.

- Kujovich, J.L., 2011. Factor v Leiden thrombophilia. Genetics in medicine, 13 (1), pp.1-16.
- Mackay, Ian M., 2004 "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3., 190-212.

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Factor V RQ Kit Manual

Autumn 2025, Version 3.1

For Real-Time PCR Detection of Factor V Leiden (G1691A)
For Research Use Only

 24 (Cat# F5RQ24)

 48 (Cat# F5RQ48)

 96 (Cat# F5RQ96)

 NG-WI-ASL-18-301

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.

Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents.....	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations.....	5
9. Additionally Required Materials	5
10. General Precautions.....	5
11. Specimen, Storage and Transport.....	6
12. Interfering Substances	6
13. Internal Control.....	6
14. DNA Isolation	7
15. PCR Protocol	7
16. Devices and software	7
17. Programming Rotor-Gene	7
18. Programming StepOne.....	8

19. Programming Other Machines.....	9
20. Data Analysis: Rotor-Gene.....	9
21. Data Analysis: StepOne	11
22. Disposal Method.....	13
23. Technical Support	13
24. Contact Information	14
25. References.....	14
26. Symbols	14

1. Introduction

Factor V Leiden RQ kit provides a ready-to-use Real-Time PCR test for detection of F5 using Real-Time PCR (based on polymerase chain reaction) with Rotor-Gene or StepOne machines. All required reagents are included in the PCR Mix provided in the kit.

This kit is for Research Use Only!

2. Intended Use

Factor V Leiden RQ kit is intended for detection of Factor V Leiden (G1691A) mutation in human DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Factor V Leiden is an inherited disorder of blood clotting and refers to specific G1691A substitution in Factor V gene (R506Q). It was first reported in 1994.

3 to 8 percent of Caucasian population carries one copy of this mutation (heterozygous) and about 1 in 5000 carry two copy of it (homozygous). Heterozygosity for factor V Leiden is associated with 4-to-8-fold increased risk of venous thromboembolism (VTE) including deep venous thrombosis (DVT) and pulmonary emboly. Factor V Leiden is also associated with a two to threefold increase in relative risk of pregnancy loss, and possibly other pregnancy complications such as preeclampsia, fetal growth retardation and placental abruption.

The diagnosis of factor V Leiden requires the APC resistance assay as a coagulation screening test or sequence analysis of factor V gene.

4. Test Principle

The target sequence is detected using PCR, where primers specific to the target sequence amplify it. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
F5Q Mix	PCR Mix*	480 μ l
Thrmb MM Ctrl	Homozygous Positive Control	100 μ l
Thrmb WM Ctrl	Heterozygous Positive Control	100 μ l
Thrmb WW Ctrl	Negative Control	100 μ l
Water	PCR Grade Water	200 μ l

*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes,

and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw the kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Whole blood (0.5ml) is the preferred sample. Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for few days). For longer terms, sample should be aliquoted and stored at -20°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control

Each sample carries wild type, mutant or both alleles of factor V. These alleles serve as both the target and the Internal Control of the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, the reaction is inhibited, and the test should be repeated.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

15. PCR Protocol

Thaw the reagents on crushed ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus one for each control and one for the NTC.

Pipette 20µl of F5Q Mix directly to each tube followed by adding 5µl of controls or sample DNA.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring too.

Please note that genotype of the samples can be called only if all the four controls provided with kit are used in each test.

16. Devices and software

F5 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.


17. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring to the rotor!

F5 RQ (V3.1)

Open the F5 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); F5 0.1 is for strip tubes and F5 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file. Edit sample names.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click on Template (provided in the flash card, or accessible by kit QR

code). Click on Plate Setup. Controls and a few samples are defined. You may change the plate Setup using right-click options (copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished, click on Start Run and save the experiment on the desired location. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	58°C x 40 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

20. Data Analysis: Rotor-Gene

Analyze data according to Rotor-Gene manual. Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Green and Yellow channels and then click on “Show”.

Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test!

Also, make sure that in the Setup section, the alleles MM, WM, and WW are defined.

The “Scatter Graph Window” is now shown in the lower left window. Each dot represents one sample. The vertical axis shows the Green fluorescence and horizontal axis is for the Yellow

fluorescence. Note that M allele is detected in the Green channel and W allele in Yellow channel.

Therefore, samples are located in four regions of the graph. MM samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.

To determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 1).


Fig1. Defining genotypes of different controls of the test on Rotor-Gene

Now turn on all the samples and the results are displayed in the tab “Scatter Analysis Results” (Figure 2).


Fig 2. Typical scatter graph for controls and samples on Rotor-Gene

Note! If the location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also, if a sample is located in between regions, it should be re-examined.

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and select “Allelic Discrimination” tab. Genotype of the samples are automatically determined by the software based on

Please note that the device software can only detect the genotype

of the samples if all three kit controls and water or a DNA-free control have been used in the test!

please ensure that the alleles MM, WM, and WW are defined in the Setup section.

Figure 3 represents the typical graph. Each dot represents one sample. All samples are color coded based on their genotype as blue for MM, Red for WW, Green for WM, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "×".

Note that the vertical axis shows the FAM fluorescence and horizontal axis is for the VIC fluorescence. Accordingly, M allele is detected in FAM channel and W allele in the VIC channel. Therefore, samples are located in four regions of the graph. MM samples with high FAM and low VIC fluorescence gather in upper left; WW samples with low FAM fluorescence and high VIC gather in lower right; WM samples with high FAM and high VIC fluorescence gather in upper right; and finally, NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence. Genotypes of samples are automatically determined by the software based on comparison with the controls (Figure 3).


Fig 3. Typical graph for controls and samples on StepOne

Note! If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid and test should be repeated.

22. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

23. Technical Support

For technical support, contact us via
Phone +98-9936223241
Email: info@novingene.com

24. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124




Email: info@novingene.com

Website: www.novingene.com

25. References

- Hotoleanu, C., 2017. Genetic risk factors in venous thromboembolism. Thrombosis and Embolism: from Research to Clinical Practice: Volume 1, pp.253-272.
- Kaushansky, K., 2021. Williams hematology. New York: Mcgraw Hill.
- Kujovich, J.L., 2011. Factor v Leiden thrombophilia. Genetics in medicine, 13 (1), pp.1-16.
- Mackay, Ian M. "Real-time PCR in the microbiology., laboratory." Clinical microbiology and infection 10, no. 3 .2004: 190-212.

26. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com