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ايميل: info@novingene.com

وبسایت: www.novingene.com

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- Griffiths, C., Drews, S.J. and Marchant, D.J., 2017 'Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment', Clinical Microbiology Reviews, 30(1), pp. 277–319.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.

- Palese, P. and Kingsbury, D.W., 2011. Genetics of Influenza Viruses. Springer.
- Paules, C. and Subbarao, K., 2017. 'Influenza', The Lancet, 390(10095), pp. 697–708.
- Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022. Influenza. Lancet (London, England), 400(10353), 693–706.

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Flu-RSV RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time RT-PCR Detection of Influenza A, Influenza B
and RSV (RNA)

For use with Rotor-Gene and MIC

For Research Use Only

 24 (Cat# Flu-RSVRQ24)

 48 (Cat# Flu-RSVRQ48)

 96 (Cat# Flu-RSVRQ96)

 NG-WI-ASL-39-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.

Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	4
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Internal Control (IC)	6
13. RNA Isolation	6
14. RT-PCR Protocol	7
15. Devices and software	7
16. Programming of the Rotor-Gene	7
17. Programming Other Machines	9
18. Data Analysis: Rotor-Gene	9
19. Sensitivity	13
20. Disposal Method	13

21. Technical Support.....	13
22. Contact Information.....	14
23. References	13
24. Symbols.....	14

1. Introduction

Flu-RSV RQ kit provides a ready-to-use Real-Time RT-polymerase chain reaction (PCR) test designed for detecting Influenza A virus, Influenza B virus and RSV. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains different series of primers and probe for the detection of a housekeeping gene (RNase P) as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps. This kit is intended for Research Use Only!

2. Intended Use

Flu-RSV RQ kit provides a ready-to-use One-Step Multiplex Real-Time RT-PCR system for the detection of Influenza A virus, Influenza B virus and RSV with Rotor-Gene and MIC.

3. Background Information

Human influenza A and B are enveloped RNA viruses of the Orthomyxoviridae family. While both of these viruses cause seasonal flu epidemics, Influenza A virus is responsible for the majority of them and may also cause pandemics.

Respiratory Syncytial Virus (RSV) is also an enveloped RNA virus and belongs to Paramyxoviridae family. RSV mostly causes lower respiratory tract infections with cold like symptoms in healthy individuals, but it may cause also severe respiratory infections in some.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through

fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual a flash card and the following reagents:

Label	Content	Quantity
Flu-RSV Mix	RT-PCR Master mix	360 μ l*
Pos Ctrl	Positive Control	100 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.

- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

This kit does not require cDNA synthesis reagents!

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted RNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.**
- **Keep RT-PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.**

- Do not place PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

We recommend upper respiratory tract samples mainly nasopharyngeal swabs. Samples can be stored at 2-8°C for a few hours or at -20°C or lower for up to a few days.

12. Internal Control (IC)

To evaluate the possibility of RNA extraction failure and PCR inhibition and prevent false negative results, Flu-RSV Mix contains primers and probe of a housekeeping gene (RNase P) as an internal control. In a successful RNA extraction and RT-PCR test, IC should generate a CT of 20-35 in the Orange/Rox Channel.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

14. RT-PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold

block. Consider one tube for each sample plus one for positive sample and one for the negative control/water.

Aliquot 15µl of Flu-RSV Mix directly to each PCR tube. Then add 5ul of extracted RNA, Positive Control or water.

Cap the tubes and visually inspect to ensure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring too.

15. Devices and software

Flu-RSV RQ kit is designed to work with Rotor-Gene and MIC.

16. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Flu-RSV template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Flu-RSV 0.1 is for strip tubes and Flu-RSV 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Flu-RSV Mix).

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	5FI	10FI	1	10
Orange	1	5FI	10FI	1	10
Red	1	10FI	15FI	1	10
Yellow	1	10FI	15FI	1	10

Click on Start button (Green button on the top menu). On the pop-up window click start again and save program on desired location. You may also set the machine according to the following table.

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	40
	60°C x 60 sec	

Fluorescence should be collected at 60°C via Green, Yellow, Orange and Red channels.

17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table (please note that following program is not applicable to Rotor-Gene):

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	40
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

Please note that Flu-RSV Mix does not contain ROX dye as a normalizer!

18. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. To analyze data briefly, click on the Analysis menu and then under Quantitation tab double-click on Cycling the A. Green. Manually put threshold at 0.1. Repeat the above for the Yellow, the Red and put threshold at 0.05 for the Orange Channels.

Figures 1, 2, 3 and 4 represent typical graphs for the Rotor-Gene. To interpret the results, please note that:

- Signal in the **Green** channel is due to **Influenza A**, the **Yellow** channel due to **Influenza B**, the **Red** channel due to **RSV** and the **Orange** channel due to **IC or RNase P**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** for Influenza A if it is positive in the **Green** channel with sigmoid graphs and CT of less than 35 and also positive in the **Orange** channel with a CT of 20-35.
- A sample is **Positive** for Influenza B if it is positive in the **Yellow** channel with sigmoid graphs and CT of less than 35 and also positive in the **Orange** channel with a CT of 20-35.
- A sample is **Positive** for RSV if it is positive in the **Red** channel with sigmoid graphs and CT of less than 35 and also positive in the **Orange** channel with a CT of 20-35.
- A sample is **Negative** for Influenza A, B and RSV if it is negative in the Green and Yellow and Red channels while it is positive in Orange channel with a sigmoid graph and CT of 20-35.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in all four channels.

Interpretation of results is summarized in the following table.

Green	Yellow	Red	Orange	Result
+	-	-	+	Positive for Flu A
-	+	-	+	Positive for Flu B
-	-	+	+	Positive for RSV
-	-	-	+	Negative
-	-	-	-	Inconclusive


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene


Fig 3. Typical Controls graph in Red channel for Rotor-Gene


Fig 4. Typical Controls graph in Orange channel for Rotor-Gene

19. Sensitivity

The analytical detection limit of the kit was assessed with a cloned target of Influenza A, B and RSV genome and showed a limit of detection equal to 50 copies/ μ l for Influenza A, 250 copies/ μ l for Influenza B and 150 copies/ μ l for RSV.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via
Phone: +98 993-6223241
Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124




Email: info@novingene.com

Website: www.novingene.com

23. References

- Griffiths, C., Drews, S.J. and Marchant, D.J., 2017 'Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment', Clinical Microbiology Reviews, 30(1), pp. 277–319.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Palese, P. and Kingsbury, D.W., 2011. Genetics of Influenza Viruses. Springer.
- Paules, C. and Subbarao, K., 2017. 'Influenza', The Lancet, 390 (10095), pp. 697–708.
- Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022. Influenza. Lancet (London, England), 400 (10353), 693–706.

24. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com