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- High Pure Viral Nucleic Acid Kit (Cat# 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp UltraSens® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany)

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Auto-Gain Optimisation Setup [X]

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 Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

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Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10Fl	15Fl	1	10
Red	1	10Fl	15Fl	1	10
Yellow	1	10Fl	15Fl	1	10

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HCV Genotype RG Kit Manual

Autumn 2025, Version 2.1

For Genotyping of Hepatitis C Virus (HCV)
With Real-Time RT-PCR
For use with Rotor-Gene
For Research Use Only

 24 (Cat# HCVGenRG24)

 48 (Cat# HCVGenRG48)

 96 (Cat# HCVGenRG96)

 NG-WI-ASL-03-201

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents:

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models.....	4
7. Storage and Stability.....	4
8. Product Use Limitations	5
9. Additionally Required Materials.....	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. RNA Isolation	7
14. RT-PCR Protocol	7
15. Devices and software.....	7
16. Programming Rotor-Gene.....	7
17. Programming Other Machines	8
18. Data Analysis.....	9

19. Sensitivity.....	12
20. Disposal Method	12
21. Technical Support	12
22. Contact Information.....	13
23. References	13
24. Symbols	13

1. Introduction

HCV Genotype RG kit provides a ready-to-use Real-Time RT-polymerase chain reaction (PCR) test designed for genotyping of Hepatitis C virus RNA using Real-Time PCR (based on polymerase chain reaction) with Rotor-Gene. All required reagents are included in the PCR Mix provided in the kit.

This kit is intended for Research Use Only!

2. Intended Use

HCV Genotype RG kit is intended for the genotyping of HCV type 1, 2 and 3 using Real-Time RT-PCR. Detection is achieved using Real-Time RT-PCR and is compatible with Rotor-Gene and MIC machines.

3. Background Information

Hepatitis C virus (HCV) is the major cause of cirrhosis and hepatocellular carcinoma. So far 6 major genotypes and more than 80 subtypes have been identified.

Based on current clinical trials, pre-treatment viral load and HCV genotype are the two best indicators for medication choice as well as patient response to therapy. Therefore, once a patient tested positive for HCV, both viral load and genotype should be determined. Moreover, since HCV genotypes have significant differences in their global distribution and prevalence, genotyping can be a useful tool for epidemiologic studies.

The gold standard method for HCV genotyping is sequencing and phylogenetic studies. Since this method is not feasible for routine practice in a clinical laboratory, other methods have been used for this purpose which most of them are PCR-based.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card, and the following reagents:

Label	Content	Quantity
HCV Gen Mix	RT-PCR Mix*	300 µl
HCV-1	HCV Type 1 control	150 µl
HCV-2	HCV Type 2 control	150 µl
HCV-3	HCV Type 3 control	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- RNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

Note! This kit does not require cDNA synthesis reagents. They are already included in the Mix.

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the HCV Gen Mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted RNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice after.**
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Proper samples to test for HCV Gen are whole blood and peripheral blood, which should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate plasma for HCV detection and genotyping. Whole blood or plasma should be shipped at +4°C. Upon receipt, plasma should be separated from whole blood, then stored at +4°C for a few days or aliquoted and stored at -20°C for up to a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using either of the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany)

14. RT-PCR Protocol

Thaw kit reagents on crushed ice completely, followed by a brief mixing and a quick spin. **Pipette 12.5ul of HCV Gen Mix RG directly to each tube followed by adding 12.5ul of isolated RNA or Positive Ctrl or water.** Cap the tubes and visually inspect to make sure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring too.

15. Devices and software


HCV Gen RG kit is designed to work with Rotor-Gene and MIC.

16. Programming Rotor-Gene

Before you start the machine make sure you have attached the locking ring to the rotor!

Open the HCV Gen template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HCV Gen 0.1 is for strip tubes and HCV Gen 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HCV Gen Mix). Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10FI	15FI	1	10
Red	1	10FI	15FI	1	10
Yellow	1	10FI	15FI	1	10

17. Programming Other Machines

HCV Gen RG kit is intended for use with the Rotor-Gene machine. If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	10
	63°C x 15 sec	
	72°C x 15 sec	
4	95°C x 15 sec	45
	62°C x 15 sec	
	72°C x 15 sec	

Fluorescence should be collected at 62°C for FAM, VIC and Cy5 dyes. The Mix contains ROX with the final concentration of 30nM in reaction.

18. Data Analysis

Before analyzing results, make sure that in the sample menu, Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Analyze the data according to the Rotor-Gene manual. Each of HCV genotypes generates a signal in a different channel: the **Green channel: Type 1**, **Red channel: Type 2** and the **Yellow channel: Type 3**.

Briefly, click on Analysis menu and then under Quantitation tab double click on "Cycling A. Green". Close the pop-up window for Automatic Threshold and set threshold manually on 0.1 for the Green channel and 0.2 for the Yellow and Red channels. Figures 1, 2 and 3 represent typical graphs for the Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** for **HCV Genotype 1** when it is positive in the **Green channel**.
- A sample is **Positive** for **HCV Genotype 2** when it is positive in the **Red channel**.
- A sample is **Positive** for **HCV Genotype 3** when it is positive in the **Yellow channel**.

If a sample is Negative in all three the Green, Yellow and Red channels, there are following possibilities:

- Sample is Negative for HCV.
- Sample is Positive for HCV genotypes 1, 2 or 3 but the viral load is lower than the kit sensitivity.
- Sample is Positive for HCV genotypes 4, 5 or 6.

Note! according to published studies, type 1 is the most prevalent type in Iran with type 3 as the second most prevalent type. Type 2 has also been reported. Types 4, 5 and 6 have not been reported from Iran yet.

The Interpretation of results is summarized in the following Table.

FAM	HEX	CY5	Result
+	-	-	Pos: HCV Type 1
-	+	-	Pos: HCV Type 3
-	-	+	Pos: HCV Type 2
-	-	-	Neg for HCV or Pos for other types


Fig 1. Typical HCV graph in Green channel for Rotor-Gene


Fig 2. Typical HCV graph in Yellow channel for Rotor-Gene


Fig 3. Typical HCV graph in Red channel for Rotor-Gene

19. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned HCV genotype 1 sequence and showed a limit of detection equal to 6 IU/ μ l.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via
 Phone: +98 993-6223241
 Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393

+98-990 11813124





Email: info@novingene.com

Website: www.novingene.com

23. References

- Hnatyszyn, H.J., 2005. Chronic hepatitis C and genotyping: the clinical significance of determining HCV genotypes. Antiviral therapy, 10(1), pp.1-11.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Martinez, M.A. and Franco, S., 2020. Therapy implications of hepatitis C virus genetic diversity. Viruses, 13(1), p.41.
- Tan, S.-L., 2006. Hepatitis C viruses: Genomes and molecular biology. Wymondham, Norfolk, U.K.: Horizon Bioscience.

24. Symbols

 Research use only	 Manufacturer	 Consult instructions for use
 Lot number	 Content sufficient for <n> tests	 Use-by date
 Catalogue number	 Serial number	 Temperature limit

For more information and resources please visit our website:
www.novingene.com