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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Martinez, M.A. and Franco, S., 2020. Therapy implications of hepatitis C virus genetic diversity. Viruses, 13(1), p.41.
- Pham, T.N., MacParland, S.A., Mulrooney, P.M., Cooksley, H., Naoumov, N.V. and Michalak, T.I., 2004. Hepatitis C virus persistence after spontaneous or treatment-induced resolution of hepatitis C. Journal of virology, 78(11), pp.5867-5874.
- Tan, S.L., 2006. Hepatitis C viruses: genomes and molecular biology.

- Watzinger, F., Ebner, K. and Lion, T., 2006. Detection and monitoring of virus infections by real-time PCR. Molecular aspects of medicine, 27(2-3), pp.254-298.

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HCV RQ Kit Manual

Winter 2025, Version 6.0

For Real-Time RT-PCR Detection and Quantitation of
Hepatitis C virus (HCV) RNA
For Research Use Only

 24 (Cat# HCVRQ24)

 48 (Cat# HCVRQ48)

 96 (Cat# HCVRQ96)

 NG-WI-ASL-02-600

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	4
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. Internal Control (IC)	6
14. RNA Isolation	7
15. RT-PCR Protocol	7
16. Devices and software	8
17. Programming of the Rotor-Gene	8
18. Programming of StepOne	9
19. Programming Other Machines	10

20. Data Analysis: Rotor-Gene	10
21. Data Analysis: StepOne	12
22. Quantitation	14
23. Linear Range	15
24. Sensitivity.....	15
25. Disposal Method	15
26. Technical Support.....	15
27. Contact Information.....	16
28. References	16
29. Symbols	16

1. Introduction

HCV RQ kit provides a ready-to-use Real-Time reverse transcription polymerase chain reaction (RT-PCR) test designed for detecting and quantifying HCV RNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

HCV RQ kit is intended for detecting and quantifying HCV RNA. Detection is achieved using One-Step Real-Time RT-PCR and is compatible with Rotor-Gene, MIC or StepOne machines.

3. Background Information

Hepatitis C virus (HCV) is an RNA virus and an important human pathogen. It is estimated that about 180 million individuals are infected worldwide. It can cause acute and chronic hepatitis followed by liver cirrhosis and hepatocellular carcinoma.

Currently infection can be diagnosed by detection of HCV RNA through RT-PCR in patient sample, weeks to months before antibody response. In addition, monitoring viral load by the same method has become a clinical gold standard for evaluation of therapy success as well as disease progress.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR

facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HCV Mix	RT-PCR mix*	360 µl
HCV S1	Standard 1: 50,000 IU/ul	150 µl
HCV S2	Standard 2: 5,000 IU/ul	150 µl
HCV S3	Standard 3: 500 IU/ul	150 µl
HCV S4	Standard 4: 50 IU/ul	150 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.

- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- RNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

This kit does not require cDNA synthesis reagents!

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted RNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep RT-PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

We recommend EDTA or citrate plasma for HCV detection. Peripheral blood should be collected in sterile condition in proper and sterile tubes. Whole blood or plasma should be shipped at +4°C. Upon receipt plasma should be separated from whole blood and can be stored at +4°C for 48 hours or aliquoted and stored at -20°C for up to a few weeks.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of RNA extraction failure and PCR inhibition and prevent false negative results, the HCV RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the HCV Mix. To monitor both RNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample. The required volume

of IC is 10% of the elution buffer. For instance, if the extracted RNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample. **Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.**

If the IC is added to HCV Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of HCV Mix before it is added to the tubes.

In a successful RNA extraction and PCR test, the IC should generate a CT of 27-34 in the Yellow/VIC Channel.

14. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

To monitor RNA extraction process, internal control should be applied to the extraction process. For more details, please, refer to section 13 of this handbook.

15. RT-PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold

block. Consider one tube for each sample, plus one for each standard and one for the negative control.

If the IC is introduced during the extraction process, pipette 15µl of HCV Mix directly to each PCR tube.

If the IC is added to the HCV Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.

Then add 10ul of extracted RNA, **standard, or water to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring too.

16. Devices and software

HCV RQ kit is designed to work with Rotor-Gene, MIC and StepOne.

17. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the HCV template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HCV 0.1 is for strip tubes and HCV 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HCV Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition
☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10FI	15FI	1	10
Yellow	1	5FI	10FI	1	10

18. Programming of StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click HCV template file (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The HCV Mix contains ROX. Final concentration of ROX in the reaction is 300nM.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Analyze the data according to the Rotor-Gene Manual. A signal in the **Green** channel indicates **HCV** and a signal in the **Yellow channel** indicates **IC**.

Briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green. In the pop-up for Automatic Threshold, increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK or simply set the threshold on 0.1 for the Green and Yellow channels. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.


Fig 1. Typical HCV graph in Green channel for Rotor-Gene


Fig 2. Typical HCV graph in Yellow channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 27-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.

The interpretation of results is summarized in the following table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

21. Data Analysis: StepOne

Analyze data according to the StepOne Manual. Briefly, click on Analyze and set the threshold for the **HCV/FAM** at 0.1 and at 0.05 for the **IC/VIC**.

Figures 3 and 4 represent typical graphs for the StepOne machine.


Fig 3. Typical HCV graph in FAM channel for StepOne


Fig 4. Typical HCV graph in VIC channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/HCV channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in the FAM/HCV channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 27-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/HCV and VIC/IC channels.

The interpretation of results is summarized in the table on page 12.

22. Quantitation

The kit provides four quantitation standards with defined titers to generate a standard curve for the quantification of samples, viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently, using all four standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ul. To convert the result to copy/ml the following equation should be used:

$$\text{Result (IU/ml)} = \frac{\text{Result (IU/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the plasma volume used for RNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated RNA.

In the HCV RQ Kit, each IU equals 4 copy numbers. Thus, to report the results in copy/ml just, multiply the results of IU/ml by 4 to reach the results in copy/ml.

23. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 2,500,000 IU/ μ l to 5 IU/ μ l.

24. Sensitivity

The analytical detection limit of the kit was assessed with a cloned target and showed a limit of detection equal to 1.5 IU/ μ l.

25. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

26. Technical Support

For technical support, contact us via
Phone: +98 993-6223241
email: info@novingene.com

27. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Website: www.novingene.com

28. References

- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.
- Martinez, M.A. and Franco, S., 2020. Therapy implications of hepatitis C virus genetic diversity. *Viruses*, 13(1), p.41.
- Pham, T.N., MacParland, S.A., Mulrooney, P.M., Cooksley, H., Naoumov, N.V. and Michalak, T.I., 2004. Hepatitis C virus persistence after spontaneous or treatment-induced resolution of hepatitis C. *Journal of virology*, 78(11), pp.5867-5874.
- Tan, S.L., 2006. Hepatitis C viruses: genomes and molecular biology.
- Watzinger, F., Ebner, K. and Lion, T., 2006. Detection and monitoring of virus infections by real-time PCR. *Molecular aspects of medicine*, 27(2-3), pp.254-298.

29. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 -30°C $+10^{\circ}\text{C}$ Temperature limit

For more information and resources please visit our website; www.novingene.com