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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Flamand, L., Lautenschlager, I., Krueger, G.R.F. and Ablashi, D.V. (2014). Human herpesviruses HHV-6A, HHV-6B & HHV-7 : diagnosis and clinical management. Amsterdam: Elsevier Science.
- Agut, H., Bonnafous, P. and Gautheret-Dejean, A., 2015. Laboratory and clinical aspects of human herpesvirus 6 infections. Clinical microbiology reviews, 28(2), pp.313-335.
- Komaroff, A.L., Pellett, P.E. and Jacobson, S., 2020. Human herpesviruses 6a and 6b in brain diseases: Association versus causation. Clinical microbiology reviews, 34(1), pp.10-1128.

- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.

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# HHV6 RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection and Quantitation of  
HHV-6A and HHV-6B DNA  
For Research Use Only

 24 (Cat# HHV6RQ24)

 48 (Cat# HHV6RQ48)

 96 (Cat# HHV6RQ96)


 NG-WI-ASL-38-101

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	3
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	4
8. Product Use Limitations .....	4
9. Additionally Required Materials .....	5
10. General Precautions .....	5
11. Specimen, Storage and Transport .....	6
12. Interfering Substances .....	6
13. Internal Control (IC) .....	6
14. DNA Isolation .....	7
15. PCR Protocol .....	7
16. Devices and software .....	8
17. Programming Rotor-Gene .....	8
18. Programming StepOne .....	9
19. Programming Other Machines .....	10

20. Data Analysis: Rotor-Gene .....	13
21. Data Analysis: StepOne .....	13
22. Quantitation .....	15
23. Linear Range .....	16
24. Sensitivity.....	16
25. Disposal Method .....	16
26. Technical Support.....	16
27. Contact Information.....	17
28. References .....	17
29. Symbols.....	18

## 1. Introduction

HHV6 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying HHV-6A and HHV-6B DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

## 2. Intended Use

HHV6 RQ kit is intended for detecting and quantifying HHV-6A and HHV-6B viruses. Detection is achieved using Real-Time PCR.

## 3. Background Information

Human Herpesvirus 6 (HHV-6) is a double-stranded DNA virus and a member of the Herpesviridae family. It comprises two distinct herpesvirus species of HHV-6A and HHV-6B. Primary infection occurs mostly in early childhood followed by a life-long latent infection. While both viruses may get re-activated later, HHV-6B is more prone to it especially in post-transplant patients. Both viruses have also been related to nervous system complications including meningitis and encephalitis.

## 4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence

provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HHV6 Mix	PCR mix*	360 µl
HHV6 S1	Standard 1: 10,000 copies/µl	150 µl
HHV6 S2	Standard 2: 1,000 copies/µl	150 µl
HHV6 S3	Standard 3: 100 copies/µl	150 µl
HHV6 S4	Standard 4: 10 copies/µl	150 µl
Internal Ctrl	Internal Control	250ul
Water	PCR Grade Water	200 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.

- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

### 11. Specimen, Storage and Transport

Proper samples to test for HHV6 are peripheral blood and whole blood, which should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate plasma for CMV detection. Whole blood or plasma should be shipped at +4°C. Sample can be stored at +4°C for 48 hours or aliquoted and stored at -20°C for up to a few weeks.

### 12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

### 13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the HHV6 RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the HHV6 Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and the patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance,

if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to HHV6 Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 28-32 in the Yellow/VIC Channel.

#### 14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

#### 15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, four for standards and one for the negative control.

**If the IC is introduced during the extraction process, pipette 15ul of HHV6 Mix to each PCR tube.**

**If the IC is added to the [HHV6 Mix](#), add 15ul of the [prepared Mix](#) (as described in section 13) to each PCR tube.**

**Then add 10ul of [standard](#), isolated DNA or water to each tube.**

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## 16. Devices and software

HHV6 RQ kit is designed to work with Rotor-Gene, MIC and StepOne.

## 17. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the flash card provided in the kit and double-click on HHV6 0.1 or HHV6 0.2 template depending on used tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HHV6 Mix).

Click on Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

**Auto-Gain Optimisation Setup**

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to  degrees.

☒ Perform Optimisation Before 1st Acquisition  
☐ Perform Optimisation At 60 Degrees At Beginning Of Run

**Channel Settings :**

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10FI	15FI	1	10
Orange	1	5FI	10FI	1	10
Yellow	1	5FI	10FI	1	10

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the "Define Targets and Samples" menu. Please note that Herpes Mix does not contain ROX dye as a normalizer. When finished, click on "Start Run" and save the experiment to the desired location. The instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:


Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM, VIC and ROX dyes. *Please note that Herpes Mix does not contain ROX dye as a normalizer!*


## 20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Analyze the data according to the Rotor-Gene manual.


Perform quantitative analysis for the **HHV-6A (Green channel)**, the **HHV-6B (Orange channel)** and qualitative analysis for the **Internal Control (Yellow channel)**. Briefly, click on Analysis menu and then under Quantitation tab double-click on Cycling A. Green and set the threshold on 0.1. Repeat above for the Orange and Yellow channels. Figures 1,2 and 3 represent typical graphs for the Rotor-Gene machine. **Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 1.** Typical HHV-6A graph in Green channel for Rotor-Gene

**Fig 2.** Typical HHV-6B graph in Orange channel for Rotor-Gene

**Fig 3.** Typical IC graph in Yellow channel for Rotor-Gene

Consider the following points when analyzing:


- A sample is **Positive** for HHV-6A if it is positive in the Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Positive** for HHV-6B if it is positive in the Orange channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Orange are valid.
- A sample is **Negative** if it is negative in the Green and Orange channels while it is positive in the Yellow channel with a sigmoid graph and CT of 28-32.
- Results are **Invalid** and DNA isolation should be repeated if a sample is negative in the FAM and ROX with CT of higher than 35 in the VIC channels.

- Results are **Invalid** and the test should be repeated if a sample is negative in all the green, orange and yellow channels.


## 21. Data Analysis: StepOne

Analyze data according to StepOne Manual. Briefly, click on Analyze and set the threshold for the FAM at 1000 and at 500 for the ROX and the VIC. Figures 4, 5 and 6 represent typical graphs for the StepOne machine.


**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 4.** Typical HHV-6A graph in FAM channel for StepOne

**Fig 5.** Typical HHV-6B graph in ROX Channel for StepOne

**Fig 6.** Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** for HHV-6A if it is positive in the FAM channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Positive** for HHV-6B if it is positive in the ROX channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in the FAM and ROX channels while it is positive in the VIC channel with a sigmoid graph and CT of 28-32.
- Results are **Invalid** and DNA isolation should be repeated if a sample is negative in the FAM and ROX with CT of higher than 35 in VIC channels.
- Results are **Invalid** and the test should be repeated if a sample is negative in all the FAM, ROX and VIC channels.

The interpretation of the results is summarized in the following table.

	Green/FAM	Yellow/VIC	Orange/ROX	Result
1	+	- / +	-	HHV-6A Pos
2	-	- / +	+	HHV-6B Pos
3	-	+ (CT 28-32)	-	HHV6 Neg
4	-	+ (CT > 35)	-	Invalid
5	-	-	-	Invalid

## 22. Quantitation

The kit provides four quantitation standards with defined titers to generate a standard curve for the quantification of samples, viral load. Working with Rotor-Gene machine, the standard curve from

a previous run can also be imported for the quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently, using all four standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ $\mu$ l. To convert the result to copy/ml the following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result(copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

### **23. Linear Range**

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 100,000,000 copies/ $\mu$ l to 10 copies/ $\mu$ l.

### **24. Sensitivity**

The analytical detection limit of the kit was assessed with a cloned target and showed a limit of detection equal to 10 copy/ $\mu$ l.

### **25. Disposal Method**

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

### **26. Technical Support**

For technical support, contact us via

Phone: +98 993-6223241  
email: [info@novingene.com](mailto:info@novingene.com)

## **27. Contact Information**

### **NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124






Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)

## **28. References**

- Agut, H., Bonnafous, P. and Gautheret-Dejean, A., 2015. Laboratory and clinical aspects of human herpesvirus 6 infections. *Clinical microbiology reviews*, 28(2), pp.313-335.
- Flamand, L., Lautenschlager, I., Krueger, G. and Ablashi, D. eds., 2014. *Human Herpesviruses HHV-6A, HHV-6B and HHV-7: diagnosis and clinical management*. Elsevier.
- Komaroff, A.L., Pellett, P.E. and Jacobson, S., 2020. Human herpesviruses 6a and 6b in brain diseases: Association versus causation. *Clinical microbiology reviews*, 34(1), pp.10-1128.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.

## 29. Symbols

 <b>RUO</b>	Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b>	Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b>	Catalogue number	 <b>SN</b>	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**