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HPV HR RQ3 Kit Manual

Autumn 2025, Version 2.1

For Real-Time PCR Detection of 14 High Risk and 2 Low Risk
Human Papilloma Virus (HPV) Types

For use with StepOne

For Research use only

 24 (Cat# HPVHR3RQ24)

 48 (Cat# HPVHR3RQ48)

 96 (Cat# HPVHR3RQ96)

 NG-WI-ASL-36-201

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction.....	3
2. Intended Use	3
3. Background Information	3
4. Test Principle.....	3
5. Kit Contents.....	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations.....	4
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. DNA Isolation.....	6
13. Internal Control (IC)	6
14. PCR Protocol	7
15. Programming Rotor-Gene.....	8
16. Devices and software.....	9
17. Programming StepOne	9
18. Programming Other Machines	10

19. Data Analysis: Rotor-Gene	10
20. Data Analysis: StepOne	14
21. Sensitivity.....	18
22. Disposal Method	18
23. Technical Support.....	18
24. Contact Information.....	19
25. References	19
26. Symbols	19

1. Introduction

HPV HR RQ3 kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting 14 High Risk Types and 2 Low Risk Types of HPV virus. All required reagents are included in the PCR Mix provided in the kit. HPV HR1 Mix also contains a different series of primers and probes for the detection of a housekeeping gene such as Internal Control (IC). IC prevents false negative results due to failure in extraction.

This kit is intended for Research Use Only!

2. Intended Use

The HPV HR RQ3 kit is intended for detection of HPV DNA including 14 High Risk Types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 and 2 Low Risk Types of 6 and 11. Detection is accomplished by polymerase chain reaction and is applicable with Rotor-Gene or StepOne machines.

3. Background Information

Human Papilloma Virus (HPV) is a non-enveloped double stranded DNA virus from papillomaviridae family. There are about 150 different types of HPV. HPV infection is considered the main cause of cervical cancer and has been linked to more than 99% of all cases. According to data linking HPV infection to cervical cancer, HPV types have been divided into three groups of High-Risk, Low-Risk and No-Risk. Currently 14 types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 are considered High-Risk and 6 types considered Low-Risk including 6, 11, 42, 43, 44 and 53.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR

facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HPV HR1 Mix	PCR Master mix*	480 µl
HPV HR2 Mix	PCR Master mix*	480 µl
HPV HR3 Mix	PCR Master mix*	480 µl
HPV HR4 Mix	PCR Master mix*	480 µl
Pos Ctrl	Positive Control	100 µl
Neg Ctrl	Negative Control	100 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaws more than 3 times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.

- The user manual should strictly be followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Pap-smear or collection of cervicovaginal cells can be used for DNA extraction. Sample should be shipped at +4°C and upon receipt should be stored at +4°C for a few hours or aliquoted and stored at -20°C for up to a few weeks.

12. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

13. Internal Control (IC)

To examine DNA extraction quality as well as the presence of PCR inhibitors and to prevent false negative results, primers and probe for an *Internal Control* (a housekeeping gene) is included in HPV HR1 Mix. With a proper sample, successful DNA extraction and PCR, Internal control should generate a CT of 20-30. Refer to the analysis section for more details.

14. PCR Protocol

This kit contains 4 different PCR mixes, each for detection of 3 to 6 HPV types. Therefore, each sample should be examined with all 4 mixes.

First, thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place 4 series of tubes on cold block. Each series should include one tube for each sample plus three for positive control, negative control and no template control (NTC).

For example, for examining 5 samples, 20 tubes for samples and 12 tubes for controls are required. So, a total of 32 tubes are arranged in 4 series of 8 tubes.

Pipette 20µl of HPV HR1 Mix to each of PCR tubes in the first series.

Pipette 20µl of HPV HR2 Mix to each of PCR tubes in the second series.

Pipette 20µl of HPV HR3 Mix to each of PCR tubes in the third series.

Pipette 20µl of HPV HR4 Mix to each of PCR tubes in the forth series.

Then add 5ul of extracted DNA, Pos Ctrl, Neg ctrl or water to each tube.

consider the first tube in each series for the first sample, second tube in each series for the second sample and continue respectively for the remaining samples or for the controls.

See Diagram 1 for details.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.



Diagram 1. Tubes setup for Mixes and samples.

15. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the HPV HR RQ3 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HPV HR RQ3 0.1 is for strip tubes and HPV HR RQ3 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HPV HR1 Mix).

Click on the Start button (the Green button on the top menu). On the pop-up window, click Start again and save the program on desired location.

16. Devices and software

The HPV RQ3 kit is designed to work with Rotor-Gene and StepOne.

17. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. Controls and few samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample name

on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1 cycle
2	95°C x 15 sec	45 cycles
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC and ROX dyes.

Please note that HPV HR Mixes does not contain ROX dye as normalizer!

19. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as “Positive control”. Patient samples are defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze data according to manufacturer recommendations. Briefly, click on the Analysis menu and then under Quantitation tab double click on Cycling A. Green. Close the pop-up window for Automatic Threshold and set the threshold at 0.1.

Repeat the above for Cycling A. Yellow and Orange. Figures 1, 2 and 3 represent typical graphs for Rotor-Gene machine. To interpret the results, please, note the following table. It shows each of the 4 mixes, detects which of the HPV in each channel.

	Green	Yellow	Orange
HR1 Mix	16	6/11	Albumin
HR2 Mix	18	39	59
HR3 Mix	31/33	52/58	56/66
HR4 Mix	35/45	51	68

According to the above table, results can be interpreted as below:

Note that a sample is considered Positive only if it has a log phase and sigmoid graph with a CT less than 40. In the absence of log phase and a sigmoid graph and with a CT of above 40, sample cannot be considered as Positive.

Please note that, for a valid result, each sample should be positive in Orange channel with CT of 20-30. When a sample is Negative in Orange channel or has a CT above 30, results are not valid and test should be repeated.

Negative Result: A sample is **Negative** if it is **only** positive with **HR1 Mix** in **Orange** channel and remains negative for all other mixes-channels combinations.

Positive Result: If a sample is positive with HR1 Mix in **Orange** channel with CT of 20-30, then:

- It is **Positive for HPV16** if it is positive in **Green** channel with **HR1 Mix**.
- It is **Positive for HPV6 or 11** if it is positive in **Yellow** channel with **HR1 Mix**.
- It is **Positive for HPV18** if it is positive in **Green** channel with **HR2 Mix**.
- It is **Positive for HPV39** if it is positive in **Yellow** channel with **HR2 Mix**.

- It is **Positive for HPV59** if it is positive in **Orange** channel with **HR2 Mix**.
- It is **Positive for HPV31 or 33** if it is positive in **Green** channel with **HR3 Mix**.
- It is **Positive for HPV52 or 58** if it is positive in **Yellow** channel with **HR3 Mix**.
- It is **Positive for HPV56/66** if it is positive in **Orange** channel with **HR3 Mix**.
- It is **Positive for HPV35 or 45** if it is positive in **Green** channel with **HR4 Mix**.
- It is **Positive for HPV51** if it is positive in **Yellow** channel with **HR4 Mix**.
- It is **Positive for HPV68** if it is positive in **Orange** channel with **HR4 Mix**.

The above results are summarized in following Table 1:

	Green	Yellow	Orange	Result
HR1 Mix	+	-	IC	Pos: HPV16
HR1 Mix	-	+	IC	Pos: HPV6 or 11
HR2 Mix	+	-	-	Pos: HPV 18
HR2 Mix	-	+	-	Pos HPV39
HR2 Mix	-	-	+	Pos: HPV59
HR3 Mix	+	-	-	Pos: HPV31 or 33
HR3 Mix	-	+	-	Pos: HPV52 or 58
HR3 Mix	-	-	+	Pos: HPV56 or 66
HR4 Mix	+	-	-	Pos: HPV35 or 45
HR4 Mix	-	+	-	Pos: HPV51
HR4 Mix	-	-	+	Pos: HPV68


Fig 1. Typical HPV HR graph in Green channel for Rotor-Gene


Fig 2. Typical HPV HR graph in Yellow channel for Rotor-Gene


Fig 3. Typical HPV HR graph in Orange channel for Rotor-Gene

20. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for FAM and VIC at 1000 and at 500 for ROX. To interpret the results, please, note the following table. It shows each of the 4 mixes, detects which of the HPV in each channel.

	FAM	VIC	ROX
HR1 Mix	16	6/11	Albumin
HR2 Mix	18	39	59
HR3 Mix	31/33	52/58	56/66
HR4 Mix	35/45	51	68

According to the above table, results can be interpreted as below:

Note that a sample is considered Positive only if it has a log phase and sigmoid graph with a CT less than 40. In the

absence of log phase and a sigmoid graph and with a CT of above 40, sample cannot be considered as Positive.

Please note that, for a valid result, each sample should be positive in ROX channel with CT of 20-30. When a sample is Negative in ROX channel or has a CT above 30, results are not valid, and test should be repeated.

Negative Result: A sample is **Negative** if it is **only** positive with **HR1 Mix** in **ROX** channel and remains negative for all other mixes-channels combinations.

Positive Result: If a sample is positive with HR1 Mix in **Orange** channel with CT of 20-30, then:

- It is **Positive for HPV16** if it is positive in **FAM** channel with **HR1 Mix**.
- It is **Positive for HPV6 or 11** if it is positive in **VIC** channel with **HR1 Mix**.
- It is **Positive for HPV18** if it is positive in **FAM** channel with **HR2 Mix**.
- It is **Positive for HPV39** if it is positive in **VIC** channel with **HR2 Mix**.
- It is **Positive for HPV59** if it is positive in **ROX** channel with **HR2 Mix**.
- It is **Positive for HPV31 or 33** if it is positive in **FAM** channel with **HR3 Mix**.
- It is **Positive for HPV52 or 58** if it is positive in **VIC** channel with **HR3 Mix**.
- It is **Positive for HPV56/66** if it is positive in **ROX** channel with **HR3 Mix**.
- It is **Positive for HPV35 or 45** if it is positive in **FAM** channel with **HR4**.

- It is **Positive for HPV51** if it is positive in **VIC** channel with **HR4 Mix**.
- It is **Positive for HPV68** if it is positive in **ROX** channel with **HR4 Mix**.

The above results are summarized in following Table:

	FAM	VIC	ROX	Result
HR1 Mix	+	-	IC	Pos: HPV16
HR1 Mix	-	+	IC	Pos: HPV6 or 11
HR2 Mix	+	-	-	Pos: HPV 18
HR2 Mix	-	+	-	Pos HPV39
HR2 Mix	-	-	+	Pos: HPV59
HR3 Mix	+	-	-	Pos: HPV31 or 33
HR3 Mix	-	+	-	Pos: HPV52 or 58
HR3 Mix	-	-	+	Pos: HPV56 or 66
HR4 Mix	+	-	-	Pos: HPV35 or 45
HR4 Mix	-	+	-	Pos: HPV51
HR4 Mix	-	-	+	Pos: HPV68

Figures 4, 5 and 6 represent typical graphs for StepOne machine.


Fig 4. Typical HPV HR graph in FAM channel for StepOne


Fig 5. Typical HPV HR graph in VIC channel for StepOne


Fig 6. Typical HPV HR graph in ROX channel for StepOne

21. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned HPV 16 target and the detection limit was determined as 50 copies/ μ l.

22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

23. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

24. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393

+98-990 11813124


Email: info@novingene.com

Website: www.novingene.com

25. References

- Crosbie, E.J., Einstein, M.H., Franceschi, S. and Kitchener, H.C., 2013. Human papillomavirus and cervical cancer. The Lancet, 382(9895), pp.889-899.
- Harden, M.E. and Munger, K., 2017. Human papillomavirus molecular biology. Mutation Research/Reviews in Mutation Research, 772, pp.3-12.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212 .
- Ramírez-Fort, M.K., Khan, F., Rady, P.L. and Tyring, S.K. eds., 2014. Human Papillomavirus. Karger Medical and Scientific Publishers.

26. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com