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Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

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Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

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Kit Manual

Autumn 2025, Version 2.1

For Real-Time PCR Detection of 14 High Risk and 2 Low Risk
Human Papilloma Virus (HPV) Types
For Research Use Only

 24 (Cat# HPVHR4RQ24)

 48 (Cat# HPVHR4RQ48)

 96 (Cat# HPVHR4RQ96)

 NG-WI-ASL-37-201

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. DNA Isolation	6
13. Internal Control (IC)	6
14. PCR Protocol	7
15. Programming Rotor-Gene	8
16. Programming Other Machines	9
17. Data Analysis: Rotor-Gene	10
18. Sensitivity	15
19. Disposal Method	15

20. Technical Support.....	15
21. Contact Information.....	15
22. References	15
23. Symbols.....	16

1. Introduction

HPV HR RQ4 kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting of 14 High Risk Types and 2 Low Risk Types of HPV virus. All required reagents are included in the PCR Mix provided in the kit. HPV HR1 Mix also contains a different series of primers and probes for the detection of a housekeeping gene as Internal Control (IC). IC prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

2. Intended Use

HPV HR RQ4 kit is intended for the detecting of HPV DNA including 14 High Risk Types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 and 2 Low Risk Types of 6 and 11.

Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene and MIC machine.

3. Background Information

Human Papilloma Virus (HPV) is a non-enveloped double stranded DNA virus from papillomaviridae family. There are about 150 different types of HPV. HPV infection is considered the main cause of cervical cancer and has been linked to more than 99% of all cases. According to data linking HPV infection to cervical cancer, HPV types have been divided into three groups of High-Risk, Low-Risk and No-Risk. Currently 14 types of 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 are considered High-Risk and 6 types considered Low-Risk including 6, 11, 42, 43, 44 and 53.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HPV HR1 Mix	PCR Master mix*	480 µl
HPV HR2 Mix	PCR Master mix*	480 µl
HPV HR3 Mix	PCR Master mix*	480 µl
HPV HR4 Mix	PCR Master mix*	480 µl
Pos Ctrl	Positive Control	100 µl
Neg Ctrl	Negative Control	100 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaws more than 3 times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should strictly be followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent

preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Pap-smear or collection of cervicovaginal cells can be used for DNA extraction. Sample should be shipped at +4°C and upon receipt should be stored at +4°C for a 48 hours or aliquoted and stored at -20°C for up to a few weeks.

12. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

13. Internal Control (IC)

To examine DNA extraction quality as well as the presence of PCR inhibitors and to prevent false negative results, primers and probe for an *Internal Control* (a housekeeping gene) is included in HPV HR1 Mix. With a proper sample, successful DNA extraction and PCR, Internal control should generate a CT of 20-30. Refer to analysis section for more details.

14. PCR Protocol

This kit contains 4 different PCR mixes, each for detection of 4 of HPV types. Therefore, each sample should be examined with all 4 mixes.

First, thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place 4 series of tubes on a cold block. Each series should include one tube for each sample plus three for positive control, negative control and no template control (NTC).

For example, for examining 5 samples, 20 tubes for samples and 12 tubes for controls are required. So, total of 32 tubes arranged in 4 series of 8 tubes. Diagram 1 shows tube setup for five sample.

Pipette 20µl of HPV HR1 Mix to each of PCR tubes in the first series.

Pipette 20µl of HPV HR2 Mix to each of PCR tubes in the second series.

Pipette 20µl of HPV HR3 Mix to each of PCR tubes in the third series.

Pipette 20µl of HPV HR4 Mix to each of PCR tubes in the forth series.

Then add 5ul of extracted DNA, Pos Ctrl, Neg ctrl or water to each tube.

consider the first tube in each series for the first sample, second tube in each series for the second sample and continue respectively for the remaining samples or for the controls.


Diagram1. Tubes setup for Mixes and samples

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.


Note: If using the RotorGene attach the locking ring.

15. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the HPV HR RQ4 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HPV HR RQ4 0.1 is for strip tubes and HPV HR RQ4 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HPV HR1 Mix)

Click on Start button (the Green button on the top menu). On the pop-up window, click Start again and save the program on desired location.

16. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1 cycle
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

Please note that, HPV HR Mixes does not contain ROX dye as normalizer!

17. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to manufacturer recommendations. Briefly, click on Analysis menu and then under Quantitation tab double-click on the Cycling A. Green. Close the pop-up window for Automatic Threshold and set the threshold at 0.1.

Repeat the above for the Cycling A. Yellow and the Orange and the Red. Figures 1, 2, 3 and 4 represent typical graphs for the Rotor-Gene machine. To interpret the results, please, note the following table. It shows each of the 4 mixes, detects which of the HPV in each channel.

	Green	Yellow	Orange	Red
HR1 Mix	16	6/11	Albumin	33
HR2 Mix	18	39	59	56
HR3 Mix	31	52	66	58
HR4 Mix	35	51	68	45

According to the above table, results can be interpreted as below:

Note that a sample is considered Positive only if it has a log phase and sigmoid graph with a CT less than 40. In the absence of log phase and a sigmoid graph and with a CT of above 40, sample cannot be considered as Positive.

Please note that, for a valid result, each sample should be positive in the Orange channel with CT of 20-30. When a

sample is Negative in the Orange channel or has a CT above 30, results are not valid and test should be repeated.

Negative Result: A sample is **Negative** if it is **only** positive with **HR1 Mix** in the **Orange** channel and remains negative for all other mixes-channels combinations.

Positive Result: If a sample is positive with HR1 Mix in the **Orange** channel with CT of 20-30, then:

- It is **Positive for HPV16** if it is positive in the **Green** channel with **HR1 Mix**.
- It is **Positive for HPV6 or 11** if it is positive in the **Yellow** channel with **HR1 Mix**.
- It is **Positive for HPV33** if it is positive in the **Red** channel with **HR1 Mix**.
- It is **Positive for HPV18** if it is positive in the **Green** channel with **HR2 Mix**.
- It is **Positive for HPV39** if it is positive in the **Yellow** channel with **HR2 Mix**.
- It is **Positive for HPV59** if it is positive in the **Orange** channel with **HR2 Mix**.
- It is **Positive for HPV56** if it is positive in the **Red** channel with **HR2 Mix**.
- It is **Positive for HPV31** if it is positive in the **Green** channel with **HR3 Mix**.
- It is **Positive for HPV52** if it is positive in the **Yellow** channel with **HR3 Mix**.
- It is **Positive for HPV66** if it is positive in the **Orange** channel with **HR3 Mix**.
- It is **Positive for HPV58** if it is positive in the **Red** channel with **HR3 Mix**.

- It is **Positive for HPV35** if it is positive in the **Green** channel with **HR4 Mix**.
- It is **Positive for HPV51** if it is positive in the **Yellow** channel with **HR4 Mix**.
- It is **Positive for HPV68** if it is positive in the **Orange** channel with **HR4 Mix**.
- It is **Positive for HPV45** if it is positive in the **Red** channel with **HR4 Mix**.

Above results are summarized in the following Table:

	Green	Yellow	Orange	Red	Result
HR1 Mix	+	-	IC	-	Pos: HPV16
HR1 Mix	-	+	IC	-	Pos: HPV6 or 11
HR1 Mix	-	-	IC	+	Pos: HPV33
HR2 Mix	+	-	-	-	Pos: HPV18
HR2 Mix	-	+	-	-	Pos HPV39
HR2 Mix	-	-	+	-	Pos: HPV59
HR2 Mix	-	-	-	+	Pos: HPV56
HR3 Mix	+	-	-	-	Pos: HPV31
HR3 Mix	-	+	-	-	Pos: HPV52
HR3 Mix	-	-	+	-	Pos: HPV66
HR3 Mix	-	-	-	+	Pos: HPV58
HR4 Mix	+	-	-	-	Pos: HPV35
HR4 Mix	-	+	-	-	Pos: HPV51
HR4 Mix	-	-	+	-	Pos: HPV68
HR4 Mix	-	-	-	+	Pos: HPV45


Fig 1. Typical HPV HR graph in Green channel for Rotor-Gene


Fig 2. Typical HPV HR graph in Yellow channel for Rotor-Gene


Fig 3. Typical HPV HR graph in Orange channel for Rotor-Gene


Fig 4. Typical HPV HR graph in Red channel for Rotor-Gene

18. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned HPV 16 target and the detection limit was determined as 50 copeis/ μ l.

19. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

20. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

21. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393

+98-990 11813124

Email: info@novingene.com




Website: www.novingene.com

22. References

- Crosbie, E.J., Einstein, M.H., Franceschi, S. and Kitchener, H.C., 2013. Human papillomavirus and cervical cancer. The Lancet, 382(9895), pp.889-899.

- Harden, M.E. and Munger, K., 2017. Human papillomavirus molecular biology. Mutation Research/Reviews in Mutation Research, 772, pp.3-12.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212 .
- Ramírez-Fort, M.K., Khan, F., Rady, P.L. and Tying, S.K. eds., 2014. Human Papillomavirus. Karger Medical and Scientific Publishers.

23. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <math>\leq n</math> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com