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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Arvin, A., Campadelli-Fiume, G., Mocarski, E., Moore, P.S., Roizman, B., Whitley, R. and Yamanishi, K. eds., 2007. Human herpesviruses: biology, therapy, and immunoprophylaxis.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3),pp.190-212.
- Thomasini, R.L. ed., 2020. Human Herpesvirus Infection: Biological Features, Transmission, Symptoms, Diagnosis and Treatment. BoD-Books on Demand.

- Whitley, R.J., 2015. Herpes simplex virus infections of the central nervous system. CONTINUUM: Lifelong Learning in Neurology, 21(6), pp.1704-1713.
- Widener, R.W. and Whitley, R.J., 2014. Herpes simplex virus. In Handbook of clinical neurology (Vol. 123, pp. 251-263). Elsevier.

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# HSV Typ RQ Kit Manual

Autumn 2025, Version 3.0

For Real-Time PCR Detection and Differentiation of Herpes  
Simplex Type 1 and 2  
For Research Use Only

 24 (Cat# HSVTypRQ24)

 48 (Cat# HSVTypRQ48)

 96 (Cat# HSVTypRQ96)

 NG-WI-ASL-07-300

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	4
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	5
8. Product Use Limitations .....	5
9. Additionally Required Materials .....	5
10. General Precautions .....	6
11. Specimen, storage and transport .....	6
12. Internal control .....	6
13. DNA isolation .....	7
14. PCR Protocol .....	7
15. Devices and software .....	8
16. Programming Rotor-Gene .....	8
17. Programming StepOne .....	9
18. Programming Other Machines .....	10
19. Data Analysis: Rotor-Gene .....	10

20.	Data Analysis: StepOne .....	12
21.	Sensitivity.....	15
22.	Disposal Method .....	15
23.	Technical Support.....	15
24.	Contact Information.....	15
25.	References .....	16
26.	Symbols .....	16

## 1. Introduction

HSV Typ RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and differentiating Herpes Simplex Type 1 and 2 DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

## 2. Intended Use

HSV Typ RQ kit is intended for the detecting and typing of HSV-1 and HSV-2 DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## 3. Background Information

Human Herpes Virus type 1 and 2 (HHV-1, HHV-2) known as Herpes simplex virus 1 and 2 (HSV-1, HSV-2) are large, enveloped DNA viruses with a genome of about 150 kbp double stranded DNA. HSV is a ubiquitous agent and can cause a variety of infections. Although most HSV infections in adults are usually benign, yet HSV is the most common cause of viral encephalitis and genital herpes is one of the most prevalent sexually transmitted diseases (STDs).

Despite overlapping clinical manifestations, HSV-1 and HSV-2 infections are different regarding the severity and outcome, mortality, frequency of reactivation, response to medication, and drug resistance. Therefore, defining the type is clinically important and routinely recommended.

Classically HSV-1 and HSV-2 infections are associated with oral or genital diseases, respectively. However, the differentiation of HSV-1 from HSV-2 based on anatomical site of infection is far from absolute and both HSV-1 and -2 can be isolated from skin, mucosa, and visceral organs. Also, HSV-1 has become the predominant cause of genital ulcers in some developed countries. Therefore, HSV1 and HSV-2 detection and typing based on clinical symptoms is not reliable and requires laboratory tests.

Also, in case of viral encephalitis, prompt laboratory diagnosis of HSV infection is essential for patient management and possible initiation of antiviral therapy and reducing the mortality rate.

## 4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HSV Typ Mix	PCR mix*	360 µl
HSV-1 Ctrl	HSV-1 Positive Control	150 µl
HSV-2 Ctrl	HSV-2 Positive Control	150 µl
Internal Control	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should strictly be followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## 11. Specimen, storage and transport

CSF or swab of the genital lesions should be collected in sterile condition in proper and sterile tubes. Samples should be shipped at +4°C. Upon receipt sample can be stored at +4°C for 48 hours or stored at -20°C for up to a few weeks.

## 12. Internal control

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the HSV Type RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the HSV Typ Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of

IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to HSV Typ Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the, HSV Typ Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 26-34 in the Rox/Orange Channel.

### 13. DNA isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

To monitor DNA extraction process, internal control should be applied to the extraction process. For more details, please, refer to section 12 of this handbook.

### 14. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample plus two for Positive controls and one for the Negative control.

**If the IC is introduced during the extraction process, pipette 15µl of HSV Typ Mix to each PCR tube.**

**If the IC is added to the HSV Type Mix, add 15ul of the prepared mix (as described in section 12) to each PCR tube.**

**Then add 10ul of extracted DNA, Controls, or water to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## **15. Devices and software**

HSV Typ RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

## **16. Programming Rotor-Gene**


*Before you start the machine make sure you have attached the locking ring on the rotor!*

Open the HSV Typ template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HSV Typ 0.1 is for strip tubes and HSV Typ 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain HSV Typ Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

.

## 17. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, and two positive controls, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

*Note! ROX should not be selected as reference dye in this test.*

## 18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 40 sec	

Fluorescence should be collected at 60°C for FAM, VIC and ROX dyes.

*Note! ROX should not be selected as reference dye in this test.*


## 19. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.


Analyze the data according to Rotor-Gene manual. Perform qualitative analysis for **HSV-1** on the **Green channel**, **HSV-2** on the **Yellow channel** and for **Internal Control** on the **Orange channel**. Briefly, click on Analysis menu and then under Quantitation tab double click on Cycling A. Green. Set the threshold at 0.1. Repeat the above for Cycling A. Yellow and Orange. Figures 1, 2 and 3 represent typical graphs for the Rotor-Gene machine.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**


**Fig 1.** Typical HSV controls graph in Green channel for Rotor-Gene

**Fig 2.** Typical HSV Controls graph in Yellow channel for Rotor-Gene

## HSV Typ RQ (v3.0)


**Fig 3.** Typical HSV Controls graph in Orange channel for Rotor-Gene

Consider the following points when analyzing:


- A sample is **Positive** for **HSV-1** if it is positive in the Green channel with a CT of less than 40.
- A sample is **Positive** for **HSV-2** if it is positive in the Yellow channel with a CT of less than 40.
- A sample is **Negative** if it is negative in both the Green and Yellow channels while it is positive for IC in the Orange channel with CT of 26-34.
- Results are **Inconclusive** and the test should be repeated if sample is negative in all three channels of the Green, Yellow and Orange. The interpretation of results is summarized in the table on page 15.

## 20. Data Analysis: StepOne


Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for **HSV-1/FAM** at 1,000; for **HSV-2/VIC** at 500; and for **IC/ROX** at 200. Figures 4, 5 and 6 represent typical graphs for the StepOne machine.

**Fig 4.** Typical HSV Controls graph in FAM channel for StepOne

**Fig 5.** Typical HSV Controls graph in VIC channel for StepOne

**Fig 6.** Typical HSV Controls graph in ROX channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** for **HSV-1** if it is positive in the FAM channel with a CT of less than 40.
- A sample is **Positive** for **HSV-2** if it is positive in the VIC channel with a CT of less than 40.
- A sample is **Negative** if it is negative in both the FAM and VIC channels while it is positive for IC in the ROX channel with CT of 26-34.
- Results are **Inconclusive** and the test should be repeated if sample is negative in all three channels of the FAM (HSV-1), the VIC (HSV-2) and the ROX channel (IC).

Interpretation of results is summarized in the following table.

Green/FAM	Yellow/VIC	ROX	Result
+ (CT<40)	-	+/-	HSV Typ I
-	+ (CT<40)	+/-	HSV Typ II
-	-	+ (CT 26-34)	HSV Neg
-	-	-	Invalid
-	-	+ (CT>34)	Invalid

## 21. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the target, sequence and showed a limit of detection equal to 0.5 copy/ul.

## 22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 23. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: [info@novingene.com](mailto:info@novingene.com)

## 24. Contact Information

### NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393

+98-990 11813124




Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)

## 25. References

- Arvin, A.M., 2007. Human herpesviruses: biology, therapy, and immunoprophylaxis. Cambridge; New York: Cambridge University Press.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Thomasini, R.L., 2020. Human herpesvirus infection: Biological features, transmission, symptoms, diagnosis and treatment. London, United Kingdom: IntechOpen.
- Whitley, R.J., 2015. Herpes simplex virus infections of the central nervous system. CONTINUUM: Lifelong Learning in Neurology, 21(6), pp.1704-1713.
- Widener, R.W. and Whitley, R.J., 2014. Herpes simplex virus. In Handbook of clinical neurology (Vol. 123, pp. 251-263). Elsevier.

## 26. Symbols

 <b>RUO</b> Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b> Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b> Catalogue number	 <b>SN</b> Serial number	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**