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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Zhu, T., 2008. Human Retrovirus Protocols. Springer Science & Business Media.
- Kannian, P. and Green, P.L., 2010. Human T Lymphotropic Virus Type 1 (HTLV-1): Molecular Biology and Oncogenesis. Viruses, 2(9), pp.2037–2077.
- Umberto Bertazzoni, Ciminale, V. and Maria Grazia Romanelli., 2019. Molecular Pathology of HTLV-1. Frontiers Media SA.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.

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HTLV RQ Kit Manual

Winter 2025, Version 3.0

For Real-Time PCR Detection and Quantitation of HTLV-1 provirus
For Research Use Only

 24 (Cat# HTLVRQ24)

 48 (Cat# HTLVRQ48)

 96 (Cat# HTLVRQ96)

 NG-WI-ASL-13-300

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering substances	6
13. DNA Isolation	6
14. PCR Protocol	7
15. Devices and software	7
16. Programming Rotor-Gene	7
17. Programming StepOne	9
18. Programming Other Machines	9
19. Data Analysis: Rotor-Gene	10
20. Data Analysis: StepOne	12
21. Quantitation	14

22. Linear Range	14
23. Sensitivity.....	15
24. Disposal Method	15
25. Technical Support.....	15
26. Contact Information.....	15
27. References	15
28. Symbols.....	16

1. Introduction

HTLV RQ kit provides a ready-to-use Real-Time PCR test for detection and quantitation of HTLV-I provirus DNA using Real-Time PCR (based on polymerase chain reaction) with Rotor-Gene or StepOne machines. All required reagents are included in the PCR Mix provided in the kit. The other Mix contains a different series of primers and probes for the detection of a house keeping gene as Internal Control to prevent false negative results due to failure in extraction.

This kit is for Research Use Only!

2. Intended Use

HTLV-I RQ kit is intended for the detection and quantitation of HTLV-I provirus and cellular titer (based on Albumin gene). Detection is accomplished by polymerase chain reaction and is applicable with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Human T-cell Lymphotropic Virus I (HTLV-I) is the first human retrovirus discovered (1979). Two years later HTLV-II was also discovered. About 15 to 20 million people are infected with these viruses worldwide. HTLV-I is endemic in North of Iran. Transmission is mostly through sexual contact, breastfeeding, transfusion, transplant, and intravenous drug use. Within the body, infection spreads slowly through virological synapses from cell to cell and takes 20-50 years to develop symptoms. Infection with HTLV-I may lead to leukemia/lymphoma or myelopathy while HTLV-II causes mild neurologic disorders or chronic pulmonary infection. Disease progress and pathogenesis depend partly on proviral load and monitoring that could be beneficial for patient management.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
HTLV Mix	PCR Master mix*	480 µl
Albumin Mix	PCR Master mix*	480 µl
HTLV-Alb S1	Standard 1: 100,000 copies/µl	250 µl
HTLV-Alb S2	Standard 2: 10,000 copies/µl	250 µl
HTLV-Alb S3	Standard 3: 1,000 copies/µl	250 µl
HTLV-Alb S4	Standard 4: 100 copies/µl	250 µl
HTLV-Alb S5	Standard 5: 10 copies/µl	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should strictly be followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes; and

- c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
 - Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
 - Do not place PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

Proper samples to test for HTLV are peripheral blood, which should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate tubes. Whole blood should be shipped at +4°C and aliquoted upon receipt. One aliquot should be examined for CBC (see section 9), while other aliquots will be used for DNA extraction and can be stored at +4°C for a few days or at -20°C for up to a few weeks.

Note: Do not freeze the aliquot specified for CBC!

12. Interfering substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. DNA Isolation

DNA can be extracted directly from whole blood. Alternatively, DNA can be extracted from peripheral blood monocytes (PBMC) collected using ficol.

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

14. PCR Protocol

Each sample should be examined for detection and quantitation of provirus with HTLV-I mix and also with Albumin mix for cell titer. Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider two series of tubes, one for HTLV1 and the other for Albumin. Each series should include tubes for standards, samples and negative control.

Pipette 20ul of HTLV Mix directly to each tube of series one, followed by adding 5ul of each Standard (HTLV-Alb S1-S5) or sample DNA or control.

Pipette 20ul of Albumin Mix directly to each tube of series two, followed by adding 5ul of each Standard (HTLV-Alb S1-S5) or sample DNA or control.

Cap the tubes and visually inspect to ensure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

15. Devices and software

HTLV RQ kit is designed to work with Rotor-Gene, StepOne and MIC.


16. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the HTLV template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); HTLV 0.1 is for strip tubes and HTLV 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the following image.

Select a tube number containing p190 Mix for the Green channel and tube with Albumin Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Edit sample names on both HTLV and ABL pages. Remember that

tubes containing HTLV mix should only be named in HTLV page and tubes containing Albumin Mix should only be named in Albumin page.

Make sure in the “Type” column, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and no template control as “NTC” respectively.

17. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, four standards, and a few samples are defined. Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The HTLV and Albumin mixes contain ROX. The final concentration of ROX in reaction is 300nM.

19. Data Analysis: Rotor-Gene

Before analyzing the results, make sure that in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. **Samples should be named in two separate pages.** Tubes with HTLV mix should be named in “HTLV” page, and tubes with Albumin mix should be named in Albumin page. Define standards in each page relatively. Patient samples should be defined as “unknown” and control as “Negative Control” on both pages.

Analyze the data according to Rotor-Gene manual. Perform quantitative analysis for the **HTLV-I (Green channel)**, and for the **Albumin/cell titer (Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then, click on OK. Repeat the above for the “Cycling A. Yellow”. You may also simply set threshold at 0.1 for both channels.

Figures 1 and 2 represent typical graphs for Rotor-Gene machine.


Fig 1. Typical HTLV-1 graph in Green channel for Rotor-Gene


Fig 2. Typical Albumin graph in Yellow channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT if present, is not reliable.

Consider following points when analyzing:

- A sample is **positive** if it is positive in both the Green and Yellow channels with CTs less than 40. The proviral load can be determined according to section 21.
- A sample is **negative** if it is negative in the Green channel (for HTLV-I) while it is positive for Albumin in the Yellow channel with a CT within the Albumin standards.
- Results are **inconclusive** if a sample is negative in the Yellow channel (for Albumin) and the test should be repeated. Improper DNA extraction or error in test set up could be the cause of such a result.

The interpretation of the results is summarized in below table.

	Green/FAM	Yellow/VIC	Result
1	+	+ (CT 28-40)	HTLV Pos
2	-	+ (CT 28-34)	HTLV Neg
3	-	+ (CT>34)	Invalid
4	-	-	Invalid

20. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for the **HTLV/FAM** and **Albumin/VIC** at 0.1. Figures 3 and 4 represent typical graphs for StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT if present, is not reliable.


Fig 3. Typical HTLV-1 graph in FAM channel for StepOne


Fig 4. Typical Albumin graph in VIC channel for StepOne

Consider following points when analyzing:

- A sample is **positive** if it is positive in both FAM and VIC channels with CTs less than 40. The proviral load can be determined according to section 15.
- A sample is **negative** if it is negative in FAM channel (for HTLV-I) while it is positive for Albumin in VIC channel with a CT within the Albumin standards.
- Results are **inconclusive** if a sample is negative in VIC channel (for Albumin) and the test should be repeated. Improper DNA extraction or error in test set up could be the cause of such a result.

21. Quantitation

The kit provides one series of quantitation standards with defined titers for estimation of HTLV-I proviral titer and Albumin gene or cell titer.

If DNA has been extracted from PBMC, the HTLV proviral load is twice the ratio of HTLV/Albumin.

If DNA has been extracted from the whole blood, use the equation below to calculate the proviral load:

$$\text{Proviral Load} = \frac{\text{HTLV titre} \times 2 \times 100\%}{\text{Albumin titer} \times (\text{Lymphocyte \%} + \text{Monocyte \%})}$$

22. Linear Range

The linear range of the kit was assessed with dilution series of the cloned, target and showed to be linear in the range of 1,000,000 copy/μl to 10 copy/μl.

23. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 2 copies/ μ l for HTLV-I and 1 copy/ μ l for Albumin.

24. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

25. Technical Support

For technical support, contact us via phone +98 993-6223241 at email: info@novingene.com

26. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com


Website: www.novingene.com

27. References

- Zhu, T., 2008. Human Retrovirus Protocols. Springer Science & Business Media.
- Kannian, P. and Green, P.L., 2010. Human T Lymphotropic Virus Type 1 (HTLV-1): Molecular Biology and Oncogenesis. *Viruses*, 2(9), pp.2037–2077.
- Umberto Bertazzoni, Ciminale, V. and Maria Grazia Romanelli., 2019. Molecular Pathology of HTLV-1. Frontiers Media SA.

- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.

28. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com