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# Hemochromatosis RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of H63D, C282Y mutations  
For Research Use Only

 24 (Cat# HMCHRQ24)

 48 (Cat# HMCHRQ48)

 96 (Cat# HMCHRQ96)


 NG-WI-ASL-50-101

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

## Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	3
5. Kit Contents.....	3
6. Packaging models .....	4
7. Storage and Stability .....	4
8. Product Use Limitations.....	4
9. Additionally Required Materials .....	5
10. General Precautions.....	5
11. Specimen, Storage and Transport.....	6
12. Interfering Substances .....	6
13. Internal Control.....	6
14. DNA Isolation .....	6
15. PCR Protocol .....	7
16. Devices and software .....	7
17. Programming Rotor-Gene .....	7
18. Programming StepOne.....	9

19. Programming Other Machines.....	9
20. Data Analysis: Rotor-Gene.....	9
21. Data Analysis: StepOne .....	14
22. Disposal Method.....	19
23. Technical Support .....	19
24. Contact Information .....	19
25. References.....	19
26. Symbols .....	20

## 1. Introduction

Hemochromatosis RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting H63D and C282Y mutations. All required reagents are included in the PCR Mix provided in the kit.

This kit is for Research Use Only!

## 2. Intended Use

Hemochromatosis RQ kit is intended for detecting of H63D and C282Y mutations in human DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## 3. Background Information

Hereditary Hemochromatosis is a disorder resulted from iron overload in tissues. Some mutations on HFE gene, mostly H63D, C282Y and S65C can lead to this disease. Clinical features vary and may involve liver, heart, endocrine glands and skin.

## 4. Test Principle

The target sequence is detected using PCR, where primers specific to the target sequence amplify it. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following

reagents:

Label	Content	Quantity
282/63 M RQ Mix	PCR Mix* for 282/63 Mix	480 µl
282/63 W RQ Mix	PCR Mix* for Wild Type Mix	480 µl
282/63 MM Ctrl	Homozygous Positive Control	50 µl
282/63 WM Ctrl	Heterozygous Positive Control	50 µl
282/63 WW Ctrl	Negative Control	50 µl
Water	PCR Grade Water	200 µl

\*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipetters and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice after.**
- Keep PCR Mix at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## 11. Specimen, Storage and Transport

Whole blood (0.5ml) is the preferred sample. Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for a 48 hours). For longer terms, sample should be aliquoted and stored at -20°C which is stable for a few months.

## 12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 13. Internal Control

This kit detects two mutations H63D and C282Y. Since each person carries wild type, mutant or both alleles, it serves as both the target and the Internal Control for the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, then the test should be repeated.

## 14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

Extracted sample should contain 15 ng/ul DNA, although 1-50

ng/ul DNA will also work.

## 15. PCR Protocol

Thaw the reagents on crushed ice completely, followed by a brief mixing and a quick spin. This kit includes two mixes **282/63 M Mix** which detect mutant alleles and **282/63 W mix** which detect Wild alleles. Each sample should be examined with both 282/63 M Mix (For Mutant alleles) and 282/63 W Mix (for Wild Type alleles) in two separate set of reactions. In each set, consider 1 tube for each sample, plus four tubes for controls and NTC.

Place the required number of microtubes on a cold block.

**Pipette 20ul of 282/63 M Mix to the first series of tubes and 20ul of 282/63 W Mix to each tube of second series followed by adding 5ul of controls, or sample DNA and NTC.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring too.*

*Please note that genotype of the samples can be called only if all the four controls provided with kit are used in each test!*

## 16. Devices and software

Hemochromatosis RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

## 17. Programming Rotor-Gene

*- Before you start the machine, make sure you have attached the locking ring to the rotor!*


Open the Hemochromatosis template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Hemoch 0.1 is for strip tubes and Hemoch 0.2 is for 0.2ml tubes.

## Hemochromatosis RQ (V1.1)

Program starts.

*Note:* For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Select a tube number containing 282/63 M Mix for the Green channel and tube with 282/63 W Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Edit sample names on both M Mix and W Mix pages. Remember that tubes containing 282/63 M Mix should only be named in M Mix page and tubes containing 282/63 W Mix should only be named in W Mix page.

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on "Plate Setup". Controls and a few samples are defined. You may change the plate setup using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment on the desired location. Instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

## 20. Data Analysis: Rotor-Gene

Analyze data according to Rotor-Gene manual. Perform qualitative analysis for **H63D (the FAM channel)** and **C282Y (the VIC channel)**. Briefly, click on the "Analysis" menu and then under the "Quantitation" tab, double click on "Cycling A. Green". Close the pop-up for Automatic Threshold and set the threshold to 0.05. Repeat the above for "Cycling A. Yellow" both for 282/63 W Mix and 282/63 M Mix, refer to figures 1 to 4 for typical graphs.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

- 1) Document the CTs of sample for both 282/63 W Mix and 282/63 M Mix in the Green and Yellow channels.
- 2) Calculate  $\Delta CT$  through the following equation:

**For H63D:**  $\Delta CT_{(FAM)} = FAM\ CT_{M\ Mix} - FAM\ CT_{W\ Mix}$

**For C282Y:**  $\Delta CT_{(VIC)} = VIC\ CT_{M\ Mix} - VIC\ CT_{W\ Mix}$


**Note! If a sample remains negative in a channel, use CT of 45 for calculation.**

- 3) Interpretation  $\Delta CT$  of sample, with Table 1.


	$\Delta CT < -7$	$-2 < \Delta CT < +4$	$\Delta CT > +10$
H63D (Green)	Pos Homozygote MM	Heterozygote WM	Neg Homozygote WW
C282Y (Yellow)	Pos Homozygote MM	Heterozygote WM	Neg Homozygote WW

**Table 1.** Interpretation of results


# Hemochromatosis RQ (V1.1)


**Fig 1.** Typical Control graph with 282/63 Mix in Green channel for Rotor-Gene

**Fig 2.** Typical Control graph with 282/63 M Mix in Yellow channel for Rotor-Gene

**Fig 3.** Typical Control graph with 282/63 W Mix in Green channel for Rotor-Gene

**Fig 4.** Typical Control graph with 282/63 W Mix in Yellow channel for Rotor-Gene

## Hemochromatosis RQ (V1.1)

As an example, CTs of Controls are documented on Table 2 and  $\Delta$ CT calculated for each one. Document CTs of patient's sample on Table 2 to interpretate with Table 1.

	Channel	M Mix	W Mix	$\Delta$ CT	Interpretation
MM Ctrl	Green H63D	24.9	32.1	$24.9 - 32.1 = -7.2$	$\Delta$ CT < -7 H63D MM
	Yellow C282Y	26	33.8	$26 - 33.8 = -7.8$	$\Delta$ CT < -7 C282Y MM
WM Ctrl	Green H63D	26	23.4	$26 - 23.4 = 2.6$	$-2 < \Delta$ CT < +4 H63D WM
	Yellow C282Y	27.2	25.1	$27.2 - 25.1 = 2.1$	$-2 < \Delta$ CT < +4 C282Y WM
WW Ctrl	Green H63D	37.6	22	$37.6 - 22 = 15.6$	$\Delta$ CT > +10 H63D WW
	Yellow C282Y	38	23.7	$38 - 23.7 = 14.3$	$\Delta$ CT > +10 C282Y WW
Sample	Green H63D				
	Yellow C282Y				

**Table 2.** CT and  $\Delta$ CT for controls and samples

- If a sample  $\Delta$ CT with 282/63 W Mix and 282/63 M Mix in the **Green** channel is less than -7, sample is Positive and **Homozygote MM (H63D/H63D)**.
- If a sample  $\Delta$ CT with 282/63 W Mix and 282/63 M Mix in the **Green** channel is less than 4 and more than -2, sample is **Heterozygote WM (H63H/H63D)**.
- If a sample  $\Delta$ CT with 282/63 W Mix and 282/63 M Mix in the **Green** channel is more than 10, sample is **Homozygote Wild type or WW (H63H/H63H)**.

- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **Yellow** channel is less than -7, sample is Positive and **Homozygote MM (C282Y/ C282Y)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **Yellow** channel less than 4 and more than -2, sample is **Heterozygote WM (C282C/ C282Y)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **Yellow** channel is more than 10, sample is **Homozygote Wild type or WW (C282C/ C282C)**.
- If a sample is negative for both alleles, then results are **inconclusive**, and the test should be repeated. Improper extraction could cause that.

## 21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on "Analyze" and set the threshold at 0.05 for **H63D (the FAM channel)** and **C282Y (the VIC channel)** do this for both 282/63 M Mix and 282/63 W Mix. Refer to figures 5 to 8 for typical graphs.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

- 1) Document the CTs of sample for both 282/63 W Mix and 282/63 M Mix in the FAM and VIC channels.
- 2) Calculate  $\Delta CT$  through the following equation.

**For H63D:**  $\Delta CT (FAM) = FAM CT_{M Mix} - FAM CT_{W Mix}$

**For C282Y:**  $\Delta CT (VIC) = VIC CT_{M Mix} - VIC CT_{W Mix}$

**Note! If a sample remains negative in a channel, use CT of 45 for**

calculation.

3) Interpretation sample  $\Delta CT$  with Table 3.

	$\Delta CT < -7$	$-2 < \Delta CT < +4$	$\Delta CT > +10$
H63D (FAM)	Pos Homozygote MM	Heterozygote WM	Neg Homozygote WW
C282Y (VIC)	Pos Homozygote MM	Heterozygote WM	Neg Homozygote WW

**Table 3.** Interpretation of results

As an example, CTs of Controls are documented on Table 4 and  $\Delta CT$  calculated for each one. Document CTs of patient's sample on Table 4 to interpretate with Table 3.


	Channel	M Mix	W Mix	$\Delta CT$	Interpretation
MM Ctrl	FAM H63D	24.9	32.1	$24.9 - 32.1 = -7.2$	$\Delta CT < -7$ H63D MM
	VIC C282Y	26	33.8	$26 - 33.8 = -7.8$	$\Delta CT < -7$ C282Y MM
WM Ctrl	FAM H63D	26	23.4	$26 - 23.4 = 2.6$	$-2 < \Delta CT < +4$ H63D WM
	VIC C282Y	27.2	25.1	$27.2 - 25.1 = 2.1$	$-2 < \Delta CT < +4$ C282Y WM
WW Ctrl	FAM H63D	37.6	22	$37.6 - 22 = 15.6$	$\Delta CT > +10$ H63D WW
	VIC C282Y	38	23.7	$38 - 23.7 = 14.3$	$\Delta CT > +10$ C282Y WW
Sample	FAM H63D				
	VIC C282Y				

**Table 4.** CT and  $\Delta CT$  for controls and samples


## Hemochromatosis RQ (V1.1)

---


- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **FAM** channel is less than -7, sample is Positive and **Homozygote MM (H63D/H63D)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **FAM** channel is less than 4 and more than -2, sample is **Heterozygote WM (H63H/H63D)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **FAM** channel is more than 10, sample is **Homozygote Wild type or WW (H63H/H63H)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **VIC** channel is less than -7, sample is Positive and **Homozygote MM (C282Y/ C282Y)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **VIC** channel less than 4 and more than -2, sample is **Heterozygote WM (C282C/ C282Y)**.
- If a sample  $\Delta CT$  with 282/63 W Mix and 282/63 M Mix in the **VIC** channel is more than 10, sample is **Homozygote Wild type or WW (C282C/ C282C)**.
- If a sample is negative for both alleles, then results are **inconclusive**, and test should be repeated. Improper extraction could cause that.


**Fig 5.** Typical Control graph with 282/63 M Mix in FAM channel for StepOne

**Fig 6.** Typical Control graph with 282/63 M Mix in VIC channel for StepOne

**Fig 7.** Typical Control graph with 282/63 W Mix in FAM channel for StepOne

**Fig 8.** Typical Control graph with 282/63 W Mix in VIC channel for StepOne

## **22. Disposal Method**

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## **23. Technical Support**

For technical support, contact us via

Phone +98-9936223241

Email: [info@novingene.com](mailto:info@novingene.com)

## **24. Contact Information**

### **NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)





## **25. References**

- Barton, J.C. and Edwards, C.Q., 2000. Hemochromatosis: genetics, pathophysiology, diagnosis and treatment. Cambridge: Cambridge University Press.
- Brissot, P., Troadec, M.B., Bardou-Jacquet, E., Le Lan, C., Jouanolle, A.M., Deugnier, Y. and Loreal, O., 2008. Current approach to hemochromatosis. Blood reviews, 22(4), pp.195-210.
- Mackay, Ian M., 2004. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3: 190-212.

## Hemochromatosis RQ (V1.1)

- Pietrangelo, A., 2006. Hereditary hemochromatosis. Annu. Rev. Nutr., 26, pp.251-270.

### 26. Symbols

 Research use only	 Manufacturer	 Consult instructions for use
 Lot number	 Content sufficient for <n> tests	 Use-by date
 Catalogue number	 Serial number	 Temperature limit

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**