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Herpes RQ Kit Manual

Autumn 2025, Version 2.0

For Multiplex Real-Time PCR Detection of
CMV, HSV-1, HSV-2, EBV and VZV DNA

For use with Rotor-Gene

For Research Use Only

 24 (Cat# HerpesRQ24)

 48 (Cat# HerpesRQ48)

 96 (Cat# HerpesRQ96)

 NG-WI-ASL-34-200

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Pathogen Information.....	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability.....	4
8. Product Use Limitations	4
9. Additionally Required Materials.....	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substance.....	6
13. Internal Control (IC)	6
14 DNA Isolation	7
15 PCR Protocol	7
16. Devices and software.....	8
17. Programming of the Rotor-Gene.....	8
18. Programming Other Machines	9
19. Data Analysis: Rotor-Gene	10
20. Sensitivity.....	14
21. Disposal Method	14

22. Technical Support.....	14
23. Contact Information.....	14
24. References	15
25. Symbols.....	15

1. Introduction

Herpes RQ kit provides a ready-to-use Real-Time PCR test for detection of Cytomegalovirus (CMV), Herpes Simplex virus-1 (HSV-1), Herpes Simplex virus-2 (HSV-2), Epstein Barr virus (EBV) and Varicella Zoster virus (VZV) DNA.

Kit contains two PCR Mixes, with all required reagents included. One of the Mixes also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

2. Intended Use

Herpes RQ kit is intended for detecting CMV, HSV-1, HSV-2, EBV and VZV. Detection is achieved using Multiplex Real-Time PCR. This kit is designed for use with Rotor-Gene and MIC.

3. Background Information

CMV, HSV-1, HSV-2, EBV and VZV all are members of Herpesviridae family and have double-stranded DNA genome. These viruses are among common human pathogens and involved in a wide spectrum of diseases including CNS infections.

Herpes RQ kit, provides a fast, sensitive, and also cost-effective Multiplex Real-Time PCR system for detection of multiple viruses from Herpesviridae family (CMV, HSV-1, HSV-2, EBV and VZV) in variety of samples.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR

facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
CMV-EBV Mix	PCR mix* for CMV and EBV	360 µl
HSV-VZV Mix	PCR mix* for HSV-1, HSV-2 and VZV	360 µl
HH Pos Ctrl	Positive control	100 µl
Internal Ctrl	Internal control	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.

- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its content should not be used if there is any sign of pink or red color on the WarmMark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. **Additionally Required Materials**

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. **General Precautions**

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR microtubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Peripheral blood plasma, CSF or tissue should be collected in sterile condition in proper and sterile containers. Samples can be stored at +4°C for 48 hours or stored at -20°C for up to a few weeks.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition and to prevent false negative results, the Herpes RQ kit contains an internal control (IC). This IC should be used during the extraction process to monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e. before the addition of lysis buffer) as it loses its efficiency.

In a successful DNA extraction and PCR test, IC should generate a CT of 28-32 in Yellow channel with CMV-EBV Mix.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Each sample should be evaluated for DNA of CMV, HSV-1, HSV-2, EBV and VZV, in two different series of microtubes. In the first series, examining CMV and EBV DNA, consider one tube for each sample plus one for each controls. In second series, examining HSV-1, HSV-2 and VZV DNA, consider one tube for each sample, plus one for each control. Place required number of tubes on cold block.

In first series of microtubes, pipette 15µl of CMV-EBV Mix to each PCR microtube.

In second series of microtubes, pipette 15µl of HSV-VZV Mix to each PCR microtube.

Then add 10ul of extracted DNA, positive or negative controls to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: In Rotor-Gene machine, attach the locking ring.

16. Devices and software

Herpes RQ kit is designed to work with Rotor-Gene and MIC.

17. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Herpes template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Herpes 0.1 is for strip tubes and Herpes 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain CMV-EBV Mix).

Click on the Start button (Green button on the top menu). On the pop up window click Start again and save the program on the desired location.

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	15FI	20FI	1	10
Orange	1	2FI	5FI	1	10
Yellow	1	5FI	10FI	1	10

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 20 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for Green/FAM, Yellow/VIC and Orange/ROX dyes.

Please note that Herpes Mix does not contain ROX dye as a normalizer!

19. Data Analysis

Before analyzing results, make sure that in the sample menu, Patient samples are defined as "unknown", Positive Control for the "Positive control", and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to RotorGene Manual. Briefly, click on Analysis menu and then under Quantitation tab double click on Cycling A. Green. Set the threshold at 0.1. Repeat the above for the Yellow and the Orange channels and manually put threshold on 0.1. Figures 1, 2 and 3 represent typical graphs for the RotorGene.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene


Fig 3. Typical Controls graph in Orange channel for Rotor-Gene

To interpret the results, please note that:

- In **CMV-EBV Mix**, a signal in the **Green channel** is due to **EBV**, in the **Yellow channel** due to **IC** and signal in the **Orange channel** is due to **CMV**.
- In **HSV-VZV Mix**, a signal in the **Green channel** is due to **HSV-1**, in the **Yellow channel** due to **HSV-2** and in the **Orange channel** is due to **VZV**.

Table 1 summarizes the above explanation.

Mix	Green Channel	Yellow Channel	Orange Channel
CMV-EBV	EBV	IC	CMV
HSV-VZV	HSV-1	HSV-2	VZV

Table 1: The relationship between results with each mixes in different channels

Consider following points when analyzing:

- A sample is **Positive for EBV**, if it is positive in Green channel with CMV-EBV Mix and a CT of less than 40.
- A sample is **Positive for CMV**, if it is positive in Orange channel with CMV-EBV Mix and a CT of less than 40.
- A sample is **Positive for HSV-1** if it is positive in Green channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Positive for HSV-2** if it is positive in Yellow channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Positive for VZV** if it is positive in Orange channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Negative for CMV and EBV**, if it is negative in Green and Orange channels while it is positive in Yellow channel with a sigmoid graph and a CT of 28-32.
- A sample is **Negative for HSV-1, HSV-2 and VZV**, if it is negative in Green, Yellow and Orange channels with HSV-VZV mix while it is positive in the Yellow channel of CMV-EBV Mix with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in all three channels (Green, Yellow and Orange) and also if IC is negative in Yellow channel of CMV-EBV Mix.

Interpretation of results is summarized in table 2.

Mix	Green	Yellow	Orange	Result
CMV-EBV	+	+/-	-	Pos: EBV
CMV-EBV	-	+/-	+	Pos: CMV
HSV-VZV	+	-	-	Pos: HSV-1
HSV-VZV	-	+	-	Pos: HSV-2
HSV-VZV	-	-	+	Pos: VZV

Mix	Green	Yellow	Orange	Result
CMV-EBV	-	+	-	Negative
HSV-VZV	-	-	-	

Mix	Green	Yellow	Orange	Result
CMV-EBV	-	-	-	Invalid
HSV-VZV	-	-	-	

Table 2. Interpretation of results in Green, Yellow and Orange channels of Rotor-Gene

20. Sensitivity

The analytical detection limit of the kit was assessed with a cloned target and showed a sensitivity of 0.5 copy/μl for CMV, HSV-1 and HSV-2; and 1 copy/μl for EBV and VZV.

21. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

22. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

23. Contact Information

NovinGene ParsVira






Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393
 +98-990 11813124
 Email: info@novingene.com
 Website: www.novingene.com

24. References

- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Spear, P.G. and Longnecker, R., 2003. Herpesvirus Entry: an Update. Journal of Virology, 77(19), pp.10179–10185.
- Spencer, J.V., 2005. Herpes. Philadelphia: Chelsea House Publishers.
- Välimaa, H., Seppänen, M. and Hukkanen, V., 2013. Herpes simplex. Duodecim; Laaketieteellinen Aikauskirja, 129(1), pp.31-40.

25. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com