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JAK2 MQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection and Quantitation of
JAK2 V617F Mutation
For Research Use Only

 24 (Cat# JAK2MQ24)

 48 (Cat# JAK2MQ48)

 96 (Cat# JAK2MQ96)

 NG-WI-ASL-48-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Mutation Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability.....	5
8. Product Use Limitations	5
9. Additionally Required Materials.....	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. DNA Isolation	7
14. PCR Protocol	7
15. Programming RotorGene	8
16. Programming of StepOne	9
17. Programming Other Machines	9
18. Data Analysis: RotorGene and StepOne	10
19. Analytical Sensitivity	14

20. Disposal Method	15
21. Technical Support.....	15
22. Contact Information.....	15
23. References	15
24. Symbols.....	16

1. Introduction

JAK2 MQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying of JAK2 V617F mutation in exon 14 human genomic DNA. All required reagents are included in the PCR Mix provided in the kit. The other Mix also contains a different series of primers and probes for the detection of a Wild type gene serving as internal control to prevent false negative results.

This kit is intended for Research Use Only!

2. Intended Use

JAK2 MQ kit is intended for detecting and quantifying of JAK2 V617F mutation in exon 14 human genomic DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

JAK2 (Janus Kinase 2) is a Tyrosine Kinase located in cytoplasm with an essential role in signaling pathways for cytokines and growth factors. The acquired mutation G1849T replaces valine with phenylalanine (V617F). This substitution results in constitutively active JAK2 which leads to uncontrolled cell proliferation in the absence of growth factors. This mutation is found in the majority of BCR-ABL-negative myeloproliferative disorders (MPDs) and has become a main diagnostic test for polycythemia vera (PV), essential thrombocythemia (ET) and primary myelofibrosis (PMF).

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through

fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, quick guide and a flash card and the following reagents:

Label	Content	Quantity
MJ Mix*	Mix for JAK2 (V617F)	480 µl
WJ Mix*	Mix for Wild Type	480 µl
MJ1	Standard 1 V617F: 100,000 copies/µl	150 µl
MJ2	Standard 2 V617F: 10,000 copies/µl	150 µl
MJ3	Standard 3 V617F: 1,000 copies/µl	150 µl
MJ4	Standard 4 V617F: 100 copies/µl	150 µl
MJ5	Standard 5 V617F: 10 copies/µl	150 µl
WJ1	Standard 1 Wild Type: 100,000 copies/µl	150 µl
WJ2	Standard 2 Wild Type: 10,000 copies/µl	150 µl
WJ3	Standard 3 Wild Type: 1,000 copies/µl	150 µl
WJ4	Standard 4 Wild Type: 100 copies/µl	150 µl
WJ5	Standard 5 Wild Type: 10 copies/µl	150 µl
Pos Ctrl 2.5%	Positive Control 2.5%	50 µl
Neg Ctrl	Negative Control	50 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and store on **crushed ice while working.**
- Do not place 0.2 ml PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 48 hours).

DNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

To reach the sensitivity of 0.1%, extracted DNA should have a concentration of 10-50 ng/ μ l. Samples with lower DNA content, result in reduced sensitivity.

14. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin.

Each sample should be examined with both MJ Mix (for V617F) and with WJ Mix (for Wild type allele) in two separate set of reactions. In each set, consider 1 tube for each sample as well as 5 tubes for the 5 standards (MJ or WJ) and 3 tubes for Positive Control, Negative Control and NTC. Use MJ standards for reactions with MJ mix and WJ standards for reactions with WJ mix. Place required number of tubes on a cold block.

Pipette 20 μ l of MJ Mix to the first series of tubes and 20 μ l of WJ Mix to each tube of the second group. Continue by adding 5 μ l of sample DNA, standards, Positive or Negative Control to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.


15. Programming RotorGene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the JAK2 RQ template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); JAK2 MQ 0.1 is for strip tubes and JAK2 MQ 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation.

Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for Green channel (note that Tube number 1 should contain Jak2 Mix MJ).

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10FI	15FI	-10	10

Click on Start button (Green button on the top menu). On the pop-up window click Start again and save the run file.

Edit samples on both MJ and WJ pages. Remember that tubes containing MJ Mix should only be named in MJ page and tubes containing WJ Mix should only be named in WJ page.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

16. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. Two targets of MJ and WJ are define in FAM channel. For MJ reactions, 5 standards (MJ1-MJ5) along with Positive and Negative Control, NTC and few samples have been defined. Also, for WJ reactions, 5 standards (WJ1-WJ5) along with Positive and Negative Control, NTC and few samples have been defined. You may change plate setup using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished click on Start Run and save the experiment. The instrument will start shortly.

17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	50
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM. Both MJ Mix

and WJ Mix contain ROX with final concentration of 300nM in reaction.

18. Data Analysis: RotorGene and StepOne

Analyze the data according to machine manual. Perform quantitative analysis for both **MJ** and **WJ** in the **Green/FAM channel** with threshold set at 0.1.

Figures 1 and 2 represent typical graphs for the StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical MJ graph in FAM channel for StepOne


Fig 2. Typical WJ graph in FAM channel for StepOne

Consider the following points when analyzing:

Sample analysis with JAK2 MQ kit can be done in two ways, Quantitative and Qualitative methods.

A) Quantitative analysis (V617F% calculation)

Note! To consider a sample valid, total copy number of wild and mutant alleles must be equal to 5,000-30,000 copy/ul.

V617F% value can be calculated for each patient with the below formula.

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- A sample is **Positive** if V617F% is equal or higher than 0.1% and calculated titer could be reported.
- A sample is **Negative** if V617F% is equal or lower than 0.01% and calculated number should be ignored.
- A sample is **Inconclusive** if V617F% is between 0.02% - 0.09%. So it could be Positive or Negative and calculated number should be ignored.

JAK2 V617F%	Results
$\geq 0.1\%$	Positive for V617F
$\leq 0.01\%$	Negative for V617F
$0.02\% \leq \text{Sample} \leq 0.09\%$	Inconclusive

Note: Calculate V617F percentage for JAK2 Neg Ctrl and JAK2 Pos Ctrl. Results must be according to the below table.

Controls	JAK2 V617F%
2.5% Pos Ctrl	1.5% - 4%
Neg Ctrl	$\leq 0.02\%$

B) Qualitative Analysis

For qualitative detection of JAK2 V617F, in addition to sample, Positive and Negative Controls must be assessed with both MJ and WJ Mix.

First, calculate Cut-off as below:

$$\text{Cut-off} = \Delta CT_{\text{Neg Ctrl}} = CT_{\text{MJ Mix}} - CT_{\text{WJ Mix}}$$

Note! If CT of Negative Control is higher than 45 or it remains negative in Green Channel, consider CT of 45 for it.

- Calculate delta CT of MJ and WJ Mix for each patient's sample according to the following equation.

$$\Delta CT_{\text{Sample}} = CT_{\text{MJ Mix}} - CT_{\text{WJ Mix}}$$

- Determine Inconclusive CT range by following formula.

$$\text{Cut-off} < \text{Inconclusive CT range} < \text{Cut-off} - 3.3$$

Consider following points when analyzing:

- A sample is **Positive** for V617F, If sample ΔCT is less than Inconclusive CT range.
- A sample is **Negative** for V617F, If sample ΔCT is higher than Cut-off.
- A sample is **Inconclusive**, If sample ΔCT falls within Inconclusive CT range.

JAK2 MQ (v1.1)

Enter the results in the table below.

Sample	CT _{MJ}	CT _{WJ}	Δ CT	Inconclusive CT	Evaluation
Pos Ctrl					
Neg Ctrl					
Sample					
Sample					
Sample					

As an example, If CT of Neg Ctrl with WJ is 24.8 in the Green/FAM channel and 37.2 with MJ Mix, the Δ CT or Cut-off is 12.4. and Inconclusive CT range is (12.4-3.3=9.1) 9.1 to 12.4. results are summarized in below table.

Sample	CT _{MJ}	CT _{WJ}	Δ CT	Inconclusive CT	Evaluation
Neg Ctrl	37.2	24.8	12.4	9.1 - 12.4	Negative
Pos Ctrl	29.2	23.4	5.8	$5.8 < 9.1$	Positive
Sample	41.5	28.5	13	$13 > 9.1$	Negative
Sample	30.7	29.0	1.7	$1.7 < 9.1$	Positive
Sample	39.5	28.7	10.8	$9.1 < 10.8 < 12.4$	Inconclusive

19. Analytical Sensitivity

Analytical sensitivity of the kit has been estimated about 0.1% for JAK2 V617F. To achieve this sensitivity, sample must include 5,000-30,000 copy/ μ l mutant and wild type alleles.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com


Website: www.novingene.com

23. References

- Lay, M., Mariappan, R., Gotlib, J., Dietz, L., Sebastian, S., Schrijver, I. and Zehnder, J.L., 2006. Detection of the JAK2 V617F mutation by LightCycler PCR and probe dissociation analysis. *The Journal of Molecular Diagnostics*, 8(3), pp.330-334.
- Mackay IM. Real-time PCR in the microbiology laboratory. *Clin. Microbiol. Infect.* 2004; 10 (3): 190 – 212.
- Yoo, J.H., Park, T.S., Maeng, H.Y., Sun, Y.K., Kim, Y.A., Kie, J.H., Cho, E.H., Song, J., Lee, K.A., Suh, B. and Choi, J.R., 2009. JAK2 V617F/C618R mutation in a patient with

polycythemia vera: a case study and review of the literature. Cancer genetics and cytogenetics, 189(1), pp.43-47.

24. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for $\leq n$ tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com