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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Ahsan, N., 2007. Polyomaviruses and Human Diseases. Springer Science & Business Media.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Pinto, M. and Dobson, S., 2014. BK and JC virus: a review. Journal of Infection, 68, pp.S2-S8.
- Prezioso, C., Pietropaolo, V., Moens, U. and Ciotti, M., 2023. JC polyomavirus: a short review of its biology, its association with progressive multifocal leukoencephalopathy, and the diagnostic value of different methods to manifest its activity or presence. Expert Review of Molecular Diagnostics, 23(2), pp.143-157.

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JC RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection and Quantitation of JC Virus DNA
For Research Use Only

 24 (Cat# JCRQ24)

 48 (Cat# JCRQ48)

 96 (Cat# JCRQ96)

 NG-WI-ASL-26-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction 3

2. Intended Use..... 3

3. Background Information 3

4. Test Principle 4

5. Kit Contents 4

6. Packaging models 4

7. Storage and Stability 5

8. Product Use Limitations 5

9. Additionally Required Materials 5

10. General Precautions 6

11. Specimen, Storage and Transport..... 6

12. Interfering Substances 6

13. Internal Control (IC)..... 7

14. DNA Isolation 7

15. PCR Protocol 8

16. Devices and software 8

17. Programming Rotor-Gene 8

18. Programming StepOne 9

19. Programming Other Machines	10
20. Data Analysis: Rotor-Gene.....	10
21. Data Analysis: StepOne	13
22. Quantitation.....	14
23. Linear Range	15
24. Sensitivity	15
25. Disposal Method	15
26. Technical Support	16
27. Contact Information	16
28. References.....	16
29. Symbols	17

1. Introduction

JC RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying JC virus DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

JC RQ kit is intended for detecting and quantifying JC virus DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

JC virus is a member of polyomaviruses genus with a double stranded DNA genome. BK virus belongs to this genus as well. Early infection occurs in early childhood and without clinical symptoms. Then JC virus establishes a lifelong infection usually in kidney, central nervous system or lymphocytes. Usually about 80% of adult population are infected with this virus without noticing it. However, JCV may lead to significant consequences in patients with immune deficiency or very weak immune system such as people with HIV/AIDs, lymphoma, multiple sclerosis (MS) and transplant patients. There, JCV reactivation may result on a serious brain infection called progressive multifocal leukoencephalopathy

(PML). Therefore, it is crucial to monitor JCV titer in above mentioned target group.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
JC Mix	PCR mix*	360 µl
JC S1	Standard 1: 100,000 copy/µl	150 µl
JC S2	Standard 2: 10,000 copy/µl	150 µl
JC S3	Standard 3: 1,000 copy/µl	150 µl
JC S4	Standard 4: 100 copy/µl	150 µl
JC IC	Internal Control	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Peripheral blood, cerebrospinal fluid or urine should be collected in sterile condition in proper and sterile tubes. Samples can be shipped and stored at +4°C for couple of days or at -20°C for longer periods. Upon receipt plasma should be separated from whole blood and then stored as above mentioned

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the JC RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the JC Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to JC Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, Internal control should generate a CT of 27-32 in the Yellow Channel on Rotor-Gene and a CT of 28-34 in the VIC channel on StepOne.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using either of:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)

- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, four for standards and one for the negative control.

If the IC is introduced during the extraction process, pipette 15µl of JC Mix to each PCR tube.

If the IC is added to the JC Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.

Then add 10ul of extracted DNA, standard, or water to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software


JC RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the JC template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); JC 0.1 is for strip tubes and JC 0.2 is for 0.2ml tubes.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain JC Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The JC Mix contains ROX with the final concentration of 300nM in the reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze the data according to the Rotor-Gene manual. A signal in the **Green** channel indicates **JC** and a signal in the **Yellow channel** indicates **IC**.

Briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green. In the pop-up for Automatic Threshold, increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK or simply set the threshold on 0.1 for the Green and Yellow channels.

Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 27-32.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.

The interpretation of the results is summarized in the following Table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive


Fig 1. Typical JC graph in Green channel for Rotor-Gene


Fig2. Typical IC graph in Yellow channel for Rotor-Gene

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for the **JC/FAM** at 0.1 and at 0.05 for the **JC/VIC**.

Figures 3 and 4 represent typical graphs for StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical JC graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in the FAM channel while it is positive in the VIC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM and VIC channels.

22. Quantitation

The kit provides four quantitation standards with defined titers to generate a standard curve for the quantification of samples, viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in

the current run. Apparently, using all four standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ μ l. To convert the result to copy/ml the following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result (copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

23. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target, and showed to be linear in the range of 100,000,000 copies/ μ l to 50 copies/ μ l.

24. Sensitivity

The analytical detection limit of the kit was assessed with a dilution series of the cloned target and showed a limit of detection equal to 5 copies/ μ l.

25. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

26. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

email: info@novingene.com

27. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124





Email: info@novingene.com

Website: www.novingene.com

28. References

- Ahsan, N., 2007, Polyomaviruses and Human Diseases. Springer Science & Business Media.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.
- Pinto, M. and Dobson, S., 2014. BK and JC virus: a review. *Journal of Infection*, 68, pp.S2-S8.
- Prezioso, C., Pietropaolo, V., Moens, U. and Ciotti, M., 2023. JC polyomavirus: a short review of its biology, its association with progressive multifocal leukoencephalopathy, and the diagnostic value of different methods to manifest its activity or presence. *Expert Review of Molecular Diagnostics*, 23(2), pp.143-157.

29. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com