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Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10Fl	15Fl	1	10
Yellow	1	10Fl	15Fl	1	10

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KRAS RQ Kit Manual

Autumn 2025, Version 3.2

For Real-Time PCR Detection of KRAS Mutations
For use with Rotor-Gene
For Research Use Only

 24 (Cat# KRASRQ24)

 48 (Cat# KRASRQ48)

 NG-WI-ASL-33-302

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.

Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Items	5
10. General Precaution	5
11. Specimen	6
12. DNA Extraction	6
13. Protocol: DNA Sample Assessment	6
14. Devices and software	7
15. Programming Rotor-Gene	7
16. Analysis: Sample Assessment	9
17. Protocol: Detection of KRAS mutations	10
18. KRAS Mutation Detection Analysis	12

19. Analytical Sensitivity	19
20. Disposal Method	20
21. Technical Support.....	20
22. Contact Information.....	20
23. References	21
24. Symbols	21

1. Introduction

KRAS RQ Kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting of the mutations in codons 12 and 13 of KRAS in human genomic DNA. All required reagents are included in the PCR Mixes provided in the kit. The Control Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC).

This kit is intended for Research Use Only!

2. Intended Use

KRAS RQ kit is intended for detecting of mutations in codons 12 and 13 of KRAS in human genomic DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene machine.

3. Background Information

Colorectal cancer (CRC) is among the most prevalent cancers worldwide. Treatment with monoclonal antibodies against Epidermal Growth Factor Receptor (Anti-EGFR) has shown to be effective for CRC patients. However, Anti-EGFR would be ineffective in presence of some KRAS mutations including mutations in codons 12 and 13. This is the same for treating non small cell lung cancer (NSCLC) with Anti-EGFR. Therefore, determining the KRAS mutation status is essential for these patients. It should be noted that, 95% and 88% of KRAS mutations in CRC and NSCLC patients respectively are reported in codons 12 and 13 including G12A, G12C, G12D, G12R, G12S, G12V and G13D.

KRAS RQ Kit provides ready to use mixes for detection of these 7 mutations based on Real-Time PCR technology.

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and following reagents:

Label	Content	Quantity	Volume
Ctrl Mix	PCR Mix for quality control	2	480 μ l
G12A Mix	PCR Mix to check G12A mutation	1	480 μ l
G12C Mix	PCR Mix to check G12C mutation	1	480 μ l
G12D Mix	PCR Mix to check G12D mutation	1	480 μ l
G12R Mix	PCR Mix to check G12R mutation	1	480 μ l
G12S Mix	PCR Mix to check G12S mutation	1	480 μ l
G12V Mix	PCR Mix to check G12V mutation	1	480 μ l
G13D Mix	PCR Mix to check G13D mutation	1	480 μ l
KRAS Pos	Positive control	1	250 μ l
KRAS Neg	Negative control	1	250 μ l
Water	PCR Grade Water	1	200 μ l

* 1 and 2 tubes for 24 and 48 reaction kits.

6. Packaging models

The kit is available in 24 and 48 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precaution

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation

where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen

Sections of paraffin block are the proper sample for patients with colorectal cancer and for NSCLC patients, core needle biopsy or aspiration fine needle sample may be used.

12. DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 50 microliters. We recommend extracting 100 microliters or higher volume of DNA, enough to repeat the test if necessary.

Note: Extracted sample should contain 10-50 ng/ul DNA.

13. Protocol: DNA Sample Assessment

Before examining a sample for KRAS mutations, quality of DNA should be assessed. If the results are within the desired range,

KRAS RQ (V3.2)

then the second test for detecting KRAS mutations will be performed. To qualify DNA extraction, follow steps below.

Note! Check the volume of extracted DNA sample and make sure it is more than 50 microliters and preferably about 100 microliters. Volumes less than 50 microliters are not enough to proceed.

First, thaw the **Ctrl Mix** on ice completely, mix by flickering followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for negative controls and water.

Pipette 20µl of Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and water to each tube.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

14. Devices and software

KRAS RQ kit is designed to work with Rotor-Gene.


15. Programming Rotor-Gene

Before you start the machine, make sure the locking ring is attached to the rotor!

Open the KRAS template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); KRAS 0.1 is for strip tubes and KRAS 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Ctrl Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

You can also use the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 20 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Mix contain ROX with final concentration of 300nM in reaction. Make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as

"unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

16. Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. In the first step, results of positive and negative controls will be evaluated.

Perform quantitative analysis for both **DNA sample (the Green channel)** and **Internal control (the Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on "Cycling A. Green". Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

A) Test is valid only if:

- 1) Water sample is **Negative** in the Green channel, and it is **Positive** in Yellow channel with CT of 28-31.
- 2) Positive control/Negative control is **Positive** in the Green channel with CT of 26-29. And it is **Positive** in the Yellow channel with CT of 28-31.

Expected results for a valid test are summarized in the Table 1.

Sample	Yellow	Green
NTC	Pos CT: 28-31	Neg
Pos/Neg Ctrl	Pos CT: 28-31	Pos CT: 26-29

Table 1. Expected results for Sample Assessment Test validity.

If above are met, then the Test results are valid and may proceed to DNA sample assessment.

B) DNA sample Quality Analysis:

- **DNA sample is qualified only if** sample is positive in the Green channel with CT of 22-30 and also positive in the Yellow channel with CT of 28-31. This sample can be further examined for KRAS mutations.

- **Sample should be diluted with water**, if a sample is positive in the Green channel with CT of less than 22 and also positive in the Yellow channel with CT of 28-31. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the KRAS mutation test.

- **Result is invalid**, if a sample in the Green channel has CT above 30. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 2.

	FAM/Green	VIC/Yellow	Result
1	+ CT: 22-30	+ CT: 27-33	Valid
2	+ CT<22	+ CT: 27-33	Sample dilution
3	+ CT>30	+ CT: 27-33	Invalid
4	+ CT: 22-30	+ CT>33	Invalid

Table 2. DNA Quality Analysis

17. Protocol: Detection of KRAS mutations

To detect KRAS mutations, each sample should be examined with eight mixes, one for DNA quantitation and one for each of seven mutations, each in a separate microtube. Therefore, to examine only one sample, 32 microtubes are required; 8 for the

sample, 8 for the Positive control, 8 for the Negative control and 8 for water (NTC). Respectively for each extra samples 8 microtubes will be added. Therefore for 2 samples, 40 microtubes and for 3 samples 48 microtubes are required. Figure1 shows tube setup for 3 samples. Place required number of tubes on cold block organized in series of 8 each.

Pipette 20 μ l of Ctrl Mix to each tube in the first series.

Pipette 20 μ l of G12A Mix to each tube in the 2nd series.

Pipette 20 μ l of G12C Mix to each tube in the 3rd series.

Pipette 20 μ l of G12D Mix to each tube in the 4th series.

Pipette 20 μ l of G12R Mix to each tube in the 5th series.

Pipette 20 μ l of G12S Mix to each tube in the 6th series.

Pipette 20 μ l of G12V Mix to each tube in the 7th series.

Pipette 20 μ l of G13D Mix to each tube in the 8th series.

Continue by adding 5 μ l of extracted DNA, Positive control, Negative Control and water to each tube. Consider the first, second and third tube in each series for positive and Negative control and water. Then next would be for samples.

See Fig 1 for details.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine and program it according to the section 15.

Note: If using Rotor-Gene attach the locking ring.


Fig 1. Tubes setup for Mixes and samples.

18. KRAS Mutation Detection Analysis

Before analyzing results of KRAS mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be analyzed. If the test is valid, then may proceed to KRAS mutation analysis.

Briefly, click on Analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable

and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

A) Test is valid only if:

- 1) Water sample is Negative in the Green channel with all Mixes, and
- 2) Water sample is Positive in the Yellow channel with Control Mix with a CT of 28-31, and
- 3) Positive control is Positive in the Green channel with Control Mix with a CT of 26-29, and
- 4) Pos control is Positive in the Yellow channel with Control Mix and CT of 28-31, and
- 5) Positive control is Positive in the Green channel **with all 7 KRAS mutation Mixes** with a CT of 26-40, and
- 6) Sample is Positive in the Green channel with Control Mix with a CT of 22-30 and Positive in Yellow channel with CT of 28-31.

Above steps for Test validation are summarized in Figure 2.

If all above are met, then the Test is valid, and results can be analyzed further for KRAS mutations.


Fig 2. Test validation flowchart

B) Data analysis for KRAS mutations

- 1) Select the patient's samples which are positive in the Green channel with KRAS specific mixes with a CT of 20-40, and document the CTs in the Table 3.
- 2) Calculate ΔCT for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the ΔCT in the Table 3 too.
- 4) If the ΔCT is within the valid range as mentioned in the Table 3, the sample has the KRAS mutation.
- 5) If a sample is positive for two or more KRAS mutations, it is considered positive for the mutation with the lowest ΔCT and negative for the other mutations.

Above steps for detection of KRAS mutations are summarized in Figure 3.

KRAS Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
Ctrl Mix	22-30		NA		
G12A Mix	20-40		≤ 12		
G12C Mix	20-40		≤ 13		
G12D Mix	20-40		≤ 8		
G12R Mix	20-40		≤ 7		
G12S Mix	20-40		≤ 8		
G12V Mix	20-40		≤ 3		
G13D Mix	20-40		≤ 6		

Table 3. Valid CT and ΔCT ranges for KRAS Mixes.

As an example, If CT of a sample is 26.7 with Control Mix, and 32.4 with G12C Mix, 30.2 with G12R Mix, 34.7 with G12V Mix, 40

KRAS RQ (V3.2)

with G13D Mix and negative for other mixes, then ΔCT is 5.7 for G12C (32.4 - 26.7), 3.5 for G12R (30.2-26.7), 8 for G12V and 13.3 for G13D (Table. 4). So, results of G12C and G12R are in valid range. In this case, patient is Positive for G12R since ΔCT is lower for this mutation.

KRAS Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
Ctrl Mix	22-30	26.7	-	-	Neg
G12A Mix	20-40	Neg	≤ 12	-	Neg
G12C Mix	20-40	32.4	≤ 13	5.7	Positive
G12D Mix	20-40	Neg	≤ 8	-	Neg
G12R Mix	20-40	30.2	≤ 7	3.5	Positive
G12S Mix	20-40	Neg	≤ 8	-	Neg
G12V Mix	20-40	34.7	≤ 3	8	Neg
G13D Mix	20-40	40	≤ 6	13.3	Neg

Table 4. KRAS mutation analysis for a specific sample.

Note that a normal sample may cross react with some of the mutation mixes. However, the calculated ΔCT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in green channel for KRAS mutation mixes, will only considered positive only if ΔCT is within the valid range.


Fig 3. Sample analysis flowchart for KRAS mutation detection

Results can also be interpreted as below;

- Sample is **Negative** for KRAS mutations, if it is negative in the Green channel with all mixes or has CT above 40.
- Sample is **Negative** for KRAS mutations if it is positive in the Green channel with CT of 20-40 for one or more KRAS mixes, and a ΔCT above valid range (which mentioned in Table 3).
- Sample is **Positive** for KRAS mutations, if sample has a valid CT and ΔCT (Table 3) in the Green channel for one of the seven KRAS mixes. Briefly,

KRAS RQ (V3.2)

- It is Positive for G12A if G12A Mix $CT \leq 40$ and $\Delta CT \leq 12$.
- It is Positive for G12C if G12C Mix $CT \leq 40$ and $\Delta CT \leq 13$.
- It is Positive for G12D if G12D Mix $CT \leq 40$ and $\Delta CT \leq 8$.
- It is Positive for G12R if G12R Mix $CT \leq 40$ and $\Delta CT \leq 7$.
- It is Positive for G12S if G12S Mix $CT \leq 40$ and $\Delta CT \leq 8$.
- It is Positive for G12V if G12V Mix $CT \leq 40$ and $\Delta CT \leq 3$.
- It is Positive for G13D if G13D Mix $CT \leq 40$ and $\Delta CT \leq 6$.
- If a sample has a valid CT and ΔCT (Table 3) for more than one of the KRAS Mix, sample is Positive for a KRAS mutation which has a lowest ΔCT .

Interpretation of results for KRAS mutation test are summarized in Table 5.

Note that, if a sample is Negative with this kit in terms of KRAS mutations, the following conditions should be considered:

- Sample is Negative only for KRAS mutations which was mentioned.
- Sample is Positive for KRAS mutations in quantities below the sensitivity of the kit.
- Sample could be Positive for other KRAS mutations.

KRAS Mix	Sample CT	Sample Δ CT	Conclusion
G12A Mix	>40	-	Neg for G12A mutation
	20-40	>12	Neg for G12A mutation
	20-40	≤ 12	Pos for G12A mutation
G12C Mix	>40	-	Neg for G12C mutation
	20-40	>13	Neg for G12C mutation
	20-40	≤ 13	Pos for G12C mutation
G12D Mix	>40	-	Neg for G12D mutation
	20-40	>8	Neg for G12D mutation
	20-40	≤ 8	Pos for G12D mutation
G12R Mix	>40	-	Neg for G12R mutation
	20-40	>7	Neg for G12R mutation
	20-40	≤ 7	Pos for G12R mutation
G12S Mix	>40	-	Neg for G12S mutation
	20-40	>8	Neg for G12S mutation
	20-40	≤ 8	Pos for G12S mutation
G12V Mix	>40	-	Neg for G12V mutation
	20-40	>3	Neg for G12V mutation
	20-40	≤ 3	Pos for G12V mutation
G13D Mix	>40	-	Neg for G13D mutation
	20-40	>6	Neg for G13D mutation
	20-40	≤ 6	Pos for G13D mutation

Table 5. Interpretation of results for KRAS mutation.

19. Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated KRAS DNA that can be identified in the background of wild-type DNA and is mentioned in Table 6. The sensitivity depends on the DNA quantity in a sample. Maximum sensitivity of is obtained when a sample with Control Mix has a CT of 22-27 in the Green channel. If CT range is 28-30, the sensitivity will decrease for some mutations.

Reaction	Ctrl Mix CT: 22-27
G12A	1%
G12C	4%
G12D	1%
G12R	2%
G12S	1%
G12V	4%
G13D	4%

Table 6. Sensitivity of the KRAS mutation detection assay.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com

Website: www.novingene.com

23. References

- Luo, J., 2021, February. KRAS mutation in pancreatic cancer. In Seminars in oncology (Vol. 48, No. 1, pp. 10-18). WB Saunders.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Timar, J. and Kashofer, K., 2020. Molecular epidemiology and diagnostics of KRAS mutations in human cancer. Cancer and Metastasis Reviews, 39, pp.1029-1038.
- Tol, J., Dijkstra, J.R., Vink-Börger, M.E., Nagtegaal, I.D., Punt, C.J., Van Krieken, J.H. and Ligtenberg, M.J., 2010. High sensitivity of both sequencing and real-time PCR analysis of KRAS mutations in colorectal cancer tissue. Journal of cellular and molecular medicine, 14(8), pp.2122-2131.

24. Symbols

RUO	Research use only		Manufacturer		Consult instructions for use
LOT	Lot number		Content sufficient for <n> tests		Use-by date
REF	Catalogue number	SN	Serial number		Temperature limit

For more information and resources please visit our website; www.novingene.com