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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Bates, P.A., 2007. Transmission of Leishmania metacyclic promastigotes by phlebotomine sand flies. International journal for parasitology, 37(10), pp.1097-1106.
- Farrell, J.P., 2002 Leishmania. New York, NY: Springer US.

- Feiz Haddad, M.H. et al., 2021 'Epidemiological study of leishmaniasis in Iran and the Middle East in the last two decades', Jundishapur Journal of Medical Sciences, 20(2), pp. 86–100.
- Mackay IM., 2004. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect.; 10 (3): 190 – 212.
- Mann, S., Frasca, K., Scherrer, S., Henao-Martínez, A.F., Newman, S., Ramanan, P. and Suarez, J.A., 2021. A review of leishmaniasis: current knowledge and future directions. Current tropical medicine reports, 8, pp.121-132.

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Leishmania RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of Leishmania DNA
For Research Use Only

 24 (Cat# LeishRQ24)

 48 (Cat# LeishRQ48)

 96 (Cat# LeishRQ96)

 NG-WI-ASL-54-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. Internal Control (IC)	7
14. DNA Isolation	7
15. PCR Protocol	8
16. Devices and software	8
17. Programming Rotor-Gene	8
18. Programming StepOne	9
19. Programming Other Machines	10

20. Data Analysis: Rotor-Gene	10
21. Data Analysis: StepOne	13
22. Sensitivity.....	15
23. Disposal Method	15
24. Technical Support.....	15
25. Contact Information.....	15
26. References	16
27. Symbols.....	16

1. Introduction

Leishmania RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting Leishmania DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

Leishmania RQ kit is intended for detecting Leishmania DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Leishmania is a parasitic protozoan that causes Leishmaniasis in vertebrates. Over 30 species of Leishmania have been known but some of them are more important as they can infect humans, such as; *L. major*, *L. tropica*, *L. aethiopica*, *L. donovani* and *L. infantum*. Phlebotomine sand flies are the intermediate host for Leishmania. Leishmaniasis is most prevalent in tropical and subtropical regions. Based on clinical symptoms, the three main phenotypic categories of disease are Cutaneous Leishmaniasis (CL), Visceral Leishmaniasis (VL), and Mucocutaneous Leishmaniasis (MCL). The infection could be asymptomatic or present with various clinical manifestations. Cutaneous Leishmaniasis is the most prevalent type of infection with Leishmaniasis. Since the symptoms of

Visceral Leishmaniasis are like some common diseases such as tuberculosis, typhoid, and malaria, rapid and specific diagnosis of this parasite is crucial for effective treatment.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
Leish Mix	PCR mix*	360 μ l
Pos Ctrl	Positive Control	150 μ l
Internal Ctrl	Internal control*	250 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- DNA extraction kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**

- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Proper sample to test for leishmania is a sample obtained from tissue samples such as skin lesions (for Cutaneous Leishmaniasis), or biopsy specimens from mucosal areas with abnormalities (for Mucocutaneous Leishmaniasis) or from bone marrow, whole blood (for Visceral Leishmaniasis). Sample can be stored at 2-8°C for a 48 hours or at -20°C or lower for up to a few days.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the Leishmania RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the Leishmania Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to Leishmania Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 28-32 in the Yellow channel on the RotorGene and a CT of 28-34 in the VIC channel on the StepOne.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using either of:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, four for standards and one for the negative control.

If the IC is introduced during the extraction process, pipette 15µl of Leish Mix to each PCR tube.

If the IC is added to the Leish Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.

Then add 10ul of isolated DNA, Control, or water to each tube.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software


Leishmania RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Leishmania template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Leishmania 0.1 is for strip tubes and Leishmania 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Leishmania Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. One Positive control and one negative control, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The Leish Mix contains ROX with the final concentration of 300nM in reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze the data according to the Rotor-Gene manual.

The signal in the **Green** channels is due to **Leishmania** and a signal in the **Yellow** channels is due to **IC**.

Briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green. In the pop-up for Automatic

Leishmania RQ (V1.1)

Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK or simply set the threshold on 0.1 for the Green channel.

Repeat the above for Cycling A. Yellow and manually put threshold on 0.1. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.

The Interpretation of results is summarized in the table on page 15.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for the **Leishmania/FAM** at 0.1 and at 0.05 for the **IC/VIC**. Figures 3 and 4 represent typical graphs for StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical Controls graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/Leishmania channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the FAM/Leishmania channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/Leishmania and VIC/IC channels.

The Interpretation of results is summarized in the following table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

22. Sensitivity

The analytical detection limit of the kit was assessed with a cloned target, and showed a limit of detection equal to 4 copies/ μ l.

23. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

24. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

25. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124




Email: info@novingene.com

Website: www.novingene.com

26. References

- Bates, P.A., 2007. Transmission of Leishmania metacyclic promastigotes by phlebotomine sand flies. International journal for parasitology, 37(10), pp.1097-1106.
- Farrell, J.P., 2002 Leishmania. New York, NY: Springer US.
- Feiz Haddad, M.H. et al., 2021 'Epidemiological study of leishmaniasis in Iran and the Middle East in the last two decades', Jundishapur Journal of Medical Sciences, 20(2), pp. 86–100.
- Mackay IM., 2004. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect.; 10 (3): 190 – 212.
- Mann, S., Frasca, K., Scherrer, S., Henao-Martínez, A.F., Newman, S., Ramanan, P. and Suarez, J.A., 2021. A review of leishmaniasis: current knowledge and future directions. Current tropical medicine reports, 8, pp.121-132.

27. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website: www.novingene.com