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Kit Manual

Autumn 2025, Version 6.0

For Real-Time PCR Quantitative Detection of
BCR-ABL transcripts (p210)
For Research Use Only

 24 (Cat# MBCR210RQ24)

 48 (Cat# MBCR210RQ48)

 96 (Cat# MBCR210RQ96)

 NG-WI-ASL-11-600

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	5
7. Storage and Stability	5
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	7
13. RNA Isolation	7
14. cDNA Synthesis	7
15. PCR Protocol	8
16. Devices and software	8
17. Programming RotorGene	8
18. Programming of StepOne	10
19. Programming Other Machines	10

20. Data Analysis: Rotor-Gene	10
21. Data Analysis: StepOne	13
22. BCR-ABL% calculation	15
23. IS Conversion	16
24. Analytical Sensitivity	17
25. Disposal Method	17
26. Technical Support	18
27. Contact Information	18
28. References	18
29. Symbols	19

1. Introduction

MBCR 210 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying BCR-ABL (p210) Transcripts. All required reagents are included in the PCR Mix provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to calculate BCR-ABL% and prevent false negative results due to failure in extraction. Note that IS-MMR calibrators are also provided with the kit.

This kit is intended for Research Use Only!

Important Note: *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

2. Intended Use

MBCR 210 RQ kit is intended for detecting of BCR-ABL (p210) transcripts (p210, b2a2 or b3a2 break points only) as well as calculation of BCR-ABL% in patients sample undergoing therapy. Additionally, by using IS-MMR, results can be calibrated and reported according to the international standard (IS).

This kit is designed for use with RotorGene, MIC and StepOne machines.

3. Background Information

BCR-ABL also known as Philadelphia chromosome is an abnormality resulted from 9;22 translocation. Consequently, ABL proto-oncogene on chromosome 9 is fused with BCR gene on chromosome 22 (b2a2 or b3a2). This fusion produces mostly 210 or 190 kDa BCR-ABL protein with constitutively active tyrosine kinase activity promoting cell proliferation and inhibition of apoptosis. The fusion gene transcript is detectable in about 95% of Chronic myelogenous leukemia (CML) patients and some cases of

Acute lymphoblastic leukemia (ALL). Also, serial monitoring of patients for identifying and measuring BCR-ABL transcripts provides more precise assessment of response to specific therapies and prediction of those in higher risk of disease progression.

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card with templates and the following reagents:

Label	Content	Quantity
MBCR 210 Mix*	Mix for BCR-ABL	480 μ l
ABL Mix*	Mix for ABL	480 μ l
MA1	Standard 1: 100,000 copies/ μ l	250 μ l
MA2	Standard 2: 10,000 copies/ μ l	250 μ l
MA3	Standard 3: 1,000 copies/ μ l	250 μ l
MA4	Standard 4: 100 copies/ μ l	250 μ l
MA5	Standard 5: 10 copies/ μ l	250 μ l
IS-MMR Calib	IS-MMR Calibrator (RNA Control)	15 μ l
High Pos RNA	RNA Control (High titer)	15 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease-free 1.5ml microtubes and PCR microtubes

- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 72 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used. For optimum results, a sample should include about 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, and Mannheim, Germany).
- TRIzol isolation reagent (Cat. no. 15596026, Invitrogene/Thermo Fisher, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime/Thermo Fisher, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

14. cDNA Synthesis

About 4 μ g of total RNA is required and should be reverse transcribed to cDNA using random hexamers. Different kits are available in the market for this purpose.

Prepare cDNA from patient RNA samples as well as IS-MMR calibrator and High Pos RNA.

Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20 μ l of cDNA add 30 μ l of nuclease free water. Upon successful RNA extraction and cDNA synthesis, ABL titer will reach above 20,000 copy/ μ l.

15. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

Each sample should be examined for both BCR-ABL fusion (p210) gene and for ABL control gene. So, two set of reactions are required. In BCR-ABL set, consider 1 tube for each sample as well as 6 tubes for the 5 standards (MA1 to MA5) and Negative control or NTC. In ABL set, consider 1 tube for each sample and 6 tubes for the 5 standards (MA1 to MA5) and Negative control or NTC. Place required number of tubes on cold block.

Pipette 20µl of MBCR 210 Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA, or standard, or control to each tube.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin the tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.

16. Devices and software

MBCR 210 RQ kit is designed to work with Rotor-Gene, MIC and StepOne.

17. Programming RotorGene


Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the p210 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); M-BCR 0.1 is for strip tubes and M-BCR 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the next page image.

MBCR 210 RQ (V6.0)

Select a tube number containing bcr Mix for the Green channel and tube with ABL Mix for the Yellow channel.

Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition
☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10Fl	15Fl	-10	10
Yellow	20	10Fl	15Fl	-10	10

Annotations in the image:

- An arrow points from the 'Green' row to a box containing 'Tube number containing MBCR 210 Mix'.
- An arrow points from the 'Yellow' row to a box containing 'Tube number containing ABL Mix'.

Buttons on the right: , , ,

Buttons at the bottom: , , ,

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Edit sample names on both p210 and ABL pages. Remember that tubes containing MBCR Mix should only be named in p210 page and tubes containing ABL Mix should only be named in ABL page. Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

18. Programming of StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, 5 standards for MBCR, 4 standards for ABL and a few samples are defined. You may change plate Setup using right-click options (copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished click on the Start Run and save the experiment. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95C x 3 min	1
2	95C x 15 sec	45
	60C x 60 sec	

Fluorescence should be collected at 60C for FAM and VIC dyes. Both of MBCR Mix and ABL Mix contain ROX with final concentration of 300nM in reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC”, respectively.

Analyze the data according to the Rotor-Gene Manual. Perform quantitative analysis for both the **MBCR (Green channel)** and the

ABL (Yellow channel). Briefly, click on Analysis menu and then, under the Quantitation tab double click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold on 0.05. Figures 1 and 2 represent typical graphs for RotorGene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical MBCR Graph in Green Channel for RotorGene


Fig 2. Typical ABL Graph in Yellow Channel for RotorGene

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the Green/MBCR and Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the Green/MBCR channel while it is positive in the Yellow/ABL channel with a sigmoid graph and CT of 20-27 and ABL titer above 20,000 copies/ul.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green/MBCR and the Yellow/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green/MBCR channel while it is positive in the Yellow/ABL channel with sigmoid graph and CT of above 27 and ABL titer less than 20,000 copies/ul. Improper extraction or low RNA input could cause that.

Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted

or less than 1ug RNA has been used. CTs higher than 27 and ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **false negative** reports.

21. Data Analysis: StepOne

Analyze the data according to StepOne manual. Briefly, click on “Analyze” and set the threshold for both the **MBCR/FAM** and **ABL/VIC** on 0.1.

Figures 3 and 4 represent typical graphs for the StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical MBCR Graph in FAM Channel for StepOne


Fig 4: Typical ABL Graph in VIC Channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the FAM/MBCR and the VIC/ABL channels with sigmoid graphs and CT of 20-40 for FAM and 20-30 for the VIC.
- A sample is **Negative** if it is negative in the FAM/MBCR channel while it is positive in the VIC/ABL channel with a sigmoid graph and CT of 20-27 and ABL titer above 20,000 copies/ul.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both the FAM/MBCR and the VIC/ABL channels. Improper extraction or test set up could cause that.
- Results are **inconclusive** and the test should be repeated if a sample is negative in the FAM/MBCR channel while it is positive in the VIC/ABL channel with sigmoid graph and CT of above 27 and ABL titer above 20,000 copies/ul. Improper extraction or low RNA input could cause that.

Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater

than 27 usually happens if not enough cells have been extracted or less than 1ug RNA has been used. CTs higher than 27 and ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **false negative** reports.

22. BCR-ABL% calculation

To assess the response to therapy, BCR-ABL% value for each patient can be calculated. This kit uses NCN method for this purpose (Beillard E. 2003, *Leukemia* 17:2474). In this method BCR-ABL% value is the BCR-ABL expression (titer) normalized by the ABL expression (titer) and then multiplied by 100.

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Usually, ABL gene expression is higher than BCR-ABL fusion gene expression. Therefore, results of above calculation would be a number less than 100%. This value may fall even below 0.001% in optimum cases. However, at diagnosis and in case of relapse, BCR-ABL expression may surpass ABL gene expression. As a result, the ratio of BCR-ABL/ABL would be above 100%. Such results have been reported previously too. For example, please note table 18 of the paper published by J. Gabert *et al* (2003, *Leukemia* 17:2318) for CML patients. He reported a range of 44% to 300% with an average of 86% for peripheral blood and range of 48%-440% with average of 117% for bone marrow samples.

It should be mentioned that a ratio of less than 1% within 6 months or less than 0.1% within 12 months after treatment denote optimum response to treatment while ratio of above 10% 6 months after or above 1% 12 months after treatment, denote failure in treatment (Baccarani M.2003, *Blood* 122: 872).

23. Conversion of Results to the International Scale (IS Conversion)

In 2010, the World Health Organization (WHO) introduced the first International Genetic Reference Panel for BCR-ABL RNA quantitation, establishing it as a global standard. This International Scale (IS) ensures comparable results worldwide, allowing for more precise monitoring of patients' molecular response to treatment. To achieve this standardization, the IS-MMR has been incorporated into NovinGene MBCR 210 RQ kit.

- The IS-MMR Cal and High Pos RNA contains RNA and should be converted to cDNA along with the patient sample RNA.
- cDNAs should be examined with Real-Time PCR assay provided by the kit.
- Calculate the NCN for the patient sample and the IS-MMR calibrator.

Note: The NCN result for the IS-MMR calibrator must fall within the range of 0.3-3x of the IS-MMR value (specified in the Certificate of Analysis (COA)). This NCN value for the calibrator has been determined using the Acrometrix BCR-ABL panel (ThermoFisher Scientific, cat: 956980TS, Lot: 083023).

To calculate IS-NCN for a patient sample, use the following formula:

$$\text{IS-NCN sample} = \frac{\text{NCN sample} \times \text{IS-MMR value}}{\text{NCN (IS-MMR)}}$$

24. Analytical Sensitivity

Analytical sensitivity of the kit has been determined by examining dilution series of positive cDNA in negative cDNA, providing positive samples ranging from 50% to 0.002%. Results were assessed with probit analysis. The analytical detection limit of the kit (LOD) was determined as 4 copies/ul of BCR-ABL or 0.002% in the context of 50,000 copies/ul of ABL mRNA. This means that when a cDNA contains about 50,000 copies/ul of ABL transcripts, there is 95% probability that a sample of 0.002% positive will be detected. Figure 5 shows the probit analysis graph.


Fig 5: Probit Analysis of Test Sensitivity

25. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

26. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

27. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com




Website: www.novingene.com

28. References

- Hochhaus, A., Baccarani, M., Silver, R.T., Schiffer, C., Apperley, J.F., Cervantes, F., Clark, R.E., Cortes, J.E., Deininger, M.W., Guilhot, F., Hjorth-Hansen, H., Hughes, T.P., Lipton, J.H., Mahon, F.X., Mayer, J. and Nicolini, F., 2020. European LeukemiaNet 2020 recommendations for treating chronic myeloid leukemia. *Leukemia*, 34.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3., 2004: 190-212.
- Score, J., Calasanz, M.J., Ottman, O., Pane, F., Yeh, R.F., Sobrinho-Simões, M.A., Kreil, S., Ward, D., Hidalgo-Curtis, C., Melo, J.V. and Wiemels, J., 2010. Analysis of genomic breakpoints in p190 and p210 BCR–ABL indicate distinct mechanisms of formation. *Leukemia*, 24(10), pp.1742-1750.

- Witt, M., Malgorzata Dawidowska and Szczepanski, T., 2012. Molecular Aspects of Hematologic Malignancies. Springer Science & Business Media.

29. Symbols

 Research use only	 Manufacturer	 Consult instructions for use
 Lot number	 Content sufficient for <n> tests	 Use-by date
 Catalogue number	 Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com