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of mutations/polymorphisms. MTHFR polymorphisms and disease, pp.1-20.

- Liu, F., Silva, D., Malone, M.V. and Seetharaman, K., 2018. MTHFR A1298C and C677T polymorphisms are associated with increased risk of venous thromboembolism: a retrospective chart review study. Acta Haematologica, 138(4), pp.208-215.
- Mackay, Ian M., 2004 "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3, 190-212.
- Nihal, U.Ğ.U.Z., Erden, G., Güngör, O., Ceylan, B.A.L., 2012. Determination of the frequency of MTHFR C677T and MTHFR A1298C polymorphisms in persons with polymorphic MTHFR gene. Journal of Clinical and Experimental Investigations, 3(4), pp.472-476.

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## Kit Manual

Autumn 2025, Version 2.1

For Real-Time PCR Detection of MTHFR (A1298C) Variant  
For Research Use Only

 24 (Cat# MTHFR1298RQ24)

 48 (Cat# MTHFR1298RQ48)

 96 (Cat# MTHFR1298RQ96)


 NG-WI-ASL-44-201

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction.....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle.....	3
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	4
8. Product Use Limitations.....	4
9. Additionally Required Materials .....	5
10. General Precautions .....	5
11. Specimen, Storage and Transport.....	6
12. Interfering Substances.....	6
13. Internal Control.....	6
14. DNA Isolation .....	6
15. PCR Protocol.....	7
16. Devices and software .....	7
17. Programming Rotor-Gene .....	7
18. Programming StepOne.....	8
19. Programming Other Machines.....	9

20. Data Analysis: Rotor-Gene .....	9
21. Data Analysis: StepOne .....	12
22. Disposal Method.....	14
23. Technical Support.....	15
24. Contact Information .....	15
25. References .....	15
26. Symbols.....	16

## **1. Introduction**

MTHFR-1298 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting of MTHFR A1298C mutation. All required reagents are included in the PCR Mix provided in the kit.

This kit is for Research Use Only!

## **2. Intended Use**

MTHFR-1298 RQ kit is intended for detecting of MTHFR (A1298C) variant in human DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## **3. Background Information**

MTHFR stands for Methylene tetrahydrofolate reductase which is a key regulatory enzyme in folate and homocysteine metabolism. Since, the C677T and A1298C variants show reduced activity (30-65% and 15-30% respectively), these polymorphisms are among the most important genetic determinants of plasma total homocysteine. Due to their influence on homocysteine metabolism, they have been linked to several complex disorders including arterial thrombosis, atherosclerosis, coronary heart disease, neural tube defects, several types of cancer, adverse drug interactions and osteoporosis though, findings are still controversial. Depending on the ethnic background up to 50% of population may carry the variant allele of MTHFR.

## **4. Test Principle**

The target sequence is detected using PCR, where primers specific to the target sequence amplify it. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled

probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

## 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
MTHFR Mix	PCR Mix*	480 µl
Thrmb MM	Homozygous Positive Control	100 µl
Thrmb WM	Heterozygous Positive Control	100 µl
Thrmb WW	Negative Control	100 µl
Water	PCR Grade Water	200 µl

\*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

## 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.

- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## **9. Additionally Required Materials**

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipetters and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## **10. General Precautions**

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for:  
a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice after.**
- Keep PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.

- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## **11. Specimen, Storage and Transport**

Whole blood (0.5ml) is the preferred sample. Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for few days). For longer terms, sample should be aliquoted and stored at -20°C which is stable for a few months.

## **12. Interfering Substances**

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## **13. Internal Control**

Each sample carries wild type, mutant or both alleles of MTHFR. These alleles serve as both the target and the Internal Control of the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, the reaction is inhibited, and the test should be repeated.

## **14. DNA Isolation**

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).

- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

## 15. PCR Protocol

Thaw the reagents on crushed ice completely, followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for each control and one for the NTC.

**Pipette 20ul of MTHFR Mix directly to each tube followed by adding 5ul of controls, or sample DNA.**

Cap the tubes and visually inspect to ensure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring too.*

***Please note that genotype of the samples can be called only if all the four controls provided with kit are used in each test!***

## 16. Devices and software


MTHFR-1298 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

## 17. Programming Rotor-Gene

***- Before you start the machine, make sure you have attached the locking ring on the rotor!***

Open the MTHFR A1298C template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); MTHFR 1298 0.1 is for strip tubes and MTHFR 1298 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file. Edit sample names.

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on "Plate Setup". Controls and few samples are defined. You may change plate set up using right click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on "Start

Run” and save the experiment on the desired location. Instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	50
	<b>57°C x 40 sec</b>	

Fluorescence should be collected at 57°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

## 20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as “Positive control”. Patient samples are defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze data according to Rotor-Gene manual. Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both the Green and Yellow channels and then click on “Show”.


**Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test! Also, make sure that in the Setup section, the alleles MM, WM, and WW are defined.**

The “Scatter Graph Window” is now shown in the lower left window. Each dot represents one sample. The vertical axis shows the Green fluorescence and horizontal axis is for the Yellow

fluorescence. Note that M allele is detected in the Green channel and W allele in the Yellow channel.

Therefore, samples are located in four regions of the graph. MM samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.


In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 1).

**Fig 1.** Defining genotypes of different controls of the test on Rotor-Gene

Now turn on all the samples and the results are displayed in the tab “Scatter Analysis Results” (Figure 2).

**Fig 2.** Typical scatter graph for controls and samples on Rotor-Gene

**Note!** If location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also, if a sample is located in between regions, it should be re-examined.

## 21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and select “Allelic Discrimination” tab. Genotype of the samples are automatically determined by the software based on comparison of their fluorescence with the controls.


Please note that the device software can only detect the genotype

of the samples if all three kit controls and water or a DNA-free control have been used in the test!

please ensure that the alleles MM, WM, and WW are defined in the Setup section.

Figure 3 represents the typical graph. Each dot represents one sample. All samples are color coded based on their genotype as blue for MM, Red for WW, Green for WM, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "X".

Note that the vertical axis shows the FAM fluorescence and horizontal axis is for the VIC fluorescence. Accordingly, M allele is detected in the FAM channel and W allele in the VIC channel. Therefore, samples are located in four regions of the graph. MM samples with high FAM and low VIC fluorescence gather in upper left; WW samples with low FAM fluorescence and high VIC gather in lower right; WM samples with high FAM and high VIC fluorescence gather in upper right; and finally, NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence. Genotypes of samples are automatically determined by the software based on comparison with the controls (Figure 3).

**Fig 3.** Typical graph for controls and samples on StepOne

**Note!** If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid, and test should be repeated!

## 22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## **23. Technical Support**

For technical support, contact us via

Phone: +98-9936223241

Email: [info@novingene.com](mailto:info@novingene.com)

## **24. Contact Information**

### **NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: [info@novingene.com](mailto:info@novingene.com)






Website: [www.novingene.com](http://www.novingene.com)

## **25. References**

- Leclerc, D., Sibani, S. and Rozen, R., 2005. Molecular biology of methylenetetrahydrofolate reductase (MTHFR) and overview of mutations/polymorphisms. MTHFR polymorphisms and disease, pp.1-20.
- Liu, F., Silva, D., Malone, M.V. and Seetharaman, K., 2018. MTHFR A1298C and C677T polymorphisms are associated with increased risk of venous thromboembolism: a retrospective chart review study. *Acta Haematologica*, 138(4), pp.208-215.
- Mackay, Ian M., 2004 "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3, 190-212.
- Nihal, U.Ğ.U.Z., Erden, G., Güngör, O., Ceylan, B.A.L., 2012. Determination of the frequency of MTHFR C677T and MTHFR A1298C polymorphisms in persons with polymorphic MTHFR

gene. Journal of Clinical and Experimental Investigations, 3(4), pp.472-476.

## 26. Symbols

 <b>RUO</b> Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b> Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b> Catalogue number	 <b>SN</b> Serial number	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**