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- Mackay, Ian M., 2004. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3: 190-212.
- Teruya, J., 2021. Management of bleeding patients. Cham: Springer.
- Vaughan, D.E., 2005. PAI-1 and atherothrombosis. Journal of Thrombosis and Haemostasis, 3(8), pp.1879-1883.
- Zhang, Q., Jin, Y., Li, X., Peng, X., Peng, N., Song, J. and Xu, M., 2020. Plasminogen activator inhibitor-1 (PAI-1) 4G/5G promoter polymorphisms and risk of venous thromboembolism—a meta-analysis and systematic review. Vasa.

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PAI-1 RQ Kit Manual

Autumn 2025, Version 2.0

For Real-Time PCR Detection of PAI-1 4G/5G polymorphism
For Research Use Only

 24 (Cat# PAIRQ24)

 48 (Cat# PAIRQ48)

 96 (Cat# PAIRQ96)

 NG-WI-ASL-22-200

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	4
9. Additionally Required Substances Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering	6
13. Internal Control (IC)	6
14. DNA Isolation	6
15. PCR Protocol	7
16. Devices and software	7
17. Programming Rotor-Gene	7
18. Programming StepOne	9
19. Programming Other Machines	9

20. Data Analysis: Rotor-Gene	9
21. Data Analysis: StepOne	12
22. Disposal Method	14
23. Technical Support	15
24. Contact Information	15
25. References	15
26. Symbols	16

1. Introduction

PAI-1 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting of PAI-1. All required reagents are included in the PCR Mix provided in the kit.

This kit is for Research Use Only!

2. Intended Use

PAI-1 RQ kit is intended for detecting 4G/5G polymorphism in promoter region of PAI-1 gene in human genomic DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene or StepOne and MIC machines.

3. Background Information

PAI-1 (Plasminogen Activator Inhibitor 1) is a 52 kDa glycoprotein. It is the main inhibitor of activators of plasminogen, and hence fibrinolysis and blood clot degradation. PAI-1 gene is located on chromosome 7 and consists of 9 exons. It includes many variations. The most significant variation, however, lies in the promoter region at -675, which is a single guanine base insertion/deletion known as 4G/5G polymorphism. While 5G allele binds both transcription activator and repressor, 4G allele only binds activator and leads to higher expression of PAI-1. 4G allele is considered a predisposing factor and has been related to various clinically important conditions including deep vein thrombosis, atherosclerosis and myocardial infarction, cancer, sepsis, miscarriage, and asthma.

4. Test Principle

The target sequence is detected using PCR, where primers specific to the target sequence amplify its unique sequence. Real-Time

PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
PAI-1 RQ Mix	Ready to use PCR Mix*	480 µl
4G/4G Ctrl	4G/4G Homozygous Control	100 µl
4G/5G Ctrl	4G/5G Heterozygous Control	100 µl
5G/5G Ctrl	5G/5G Homozygous Control	100 µl
Water	PCR Grade Water	200 µl

* 1, 2 or 4 tubes for 24, 48 or 96 reaction kits, respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.

- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by

flickering followed by a quick spin and **store on crushed ice after.**

- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

The required sample for this test is Whole blood. Sample should be collected in sterile condition in proper tubes containing EDTA or citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 3 days). For long term storage, sample should be aliquoted and stored at -20°C which remains stable for few months. DNA can be extracted from whole blood or buffy coat.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

Either wild type or mutant alleles may function as internal control for this assay as any patients carry 4G or 5G or both alleles. If a patient is negative for both, then the test is not valid and should be repeated.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)

- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus four for controls and NTC.

Pipette 20µl of PAI-1 RQ Mix directly to each tube followed by adding 5µl of sample or control or water.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Please note that genotype of the samples can be called only if all the four controls provided with kit are used in each test!

16. Devices and software


PAI-1 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the PAI-1 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); PAI-1 0.1 is for strip tubes and PAI-1 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Click on the Start button (Green button on the top menu). On the pop up window click start again and save the run file. Edit sample names.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click on Template (provided in the flash card, or accessible by kit QR code). One negative control, 3 positive controls and some samples are defined. You may change the plate set up by using right-click options (copy, paste, clear). You may also add /remove samples or

change sample name on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. Instrument will start shortly.

Alternatively, in Setup menu, Genotyping option could be selected. Assign controls as directed by StepOne manual and set the run parameters according to the following table:

Step	Temperature and time	Cycles
1	65°C x 30 sec	Pre-PCR Read
2	95°C x 3 min	1
3	95°C x 15 sec	45
	65°C x 40 sec	
4	65°C x 30 sec	Post-PCR Read

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	50
	65°C x 40 sec	

Fluorescence should be collected at 65°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as “Positive control”. Patient samples are defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC” respectively.

Analyze data according to Rotor-Gene manual. Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both the Green and Yellow channels and then click on “Show”.

Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test! Also, make sure that in the Setup section, the alleles MM, WM, and WW are defined.

In the “Scatter Graph” window, each dot represents one sample. The vertical axis shows the Green fluorescence and horizontal axis is for the Yellow fluorescence. Note that 4G allele is detected in the Green channel and 5G allele in Yellow channel.

Therefore, samples are located in four regions of the graph. 4G/4G samples with high Green and low Yellow fluorescence gather in upper left; 5G/5G samples with low Green fluorescence and high Yellow gather in lower right; 4G/5G samples with average Green and Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then, holding the right click, draw a rectangle around each control. Each box represents one of the above mentioned regions. Label 4G/4G genotype region as “Mutant”, 4G/5G as “Heterozygous”, 5G/5G as “Wild Type” and Negative Control as “None”(Figure 1).


Fig 1. Defining genotypes of different controls of the test on Rotor-Gene

Now turn on all the samples and the results are displayed in the tab "Scatter Analysis Results" (Figure 2).


Fig 2. Typical scatter graph for controls and samples on Rotor-Gene

Note! If the location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also, if a sample is located in between regions, it should be re-examined.

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and select "Allelic Discrimination" tab. Genotype of the samples are automatically determined by the software based on comparison of their fluorescence with the controls.

Please note that the device software can only detect the genotype of the samples if all three kit controls and water or a DNA-free control have been used in the test!

please ensure that the alleles MM, WM, and WW are defined in the Setup section.

Figure 3 represents the typical graph. Each dot represents one sample. In order to analyze samples, In Plate Setup, select Assign, and define Task for each Control based on their alleles. All samples are color coded based on their genotype as Red for **4G/4G**, blue for **5G/5G**, Green for **5G/4G**, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "x".

Note that the vertical axis shows the VIC fluorescence and horizontal axis is for the FAM fluorescence. Accordingly, 4G allele is detected in the FAM channel and 5G allele in the VIC channel. Therefore, samples are located in four regions of the graph. 5G/5G samples with high VIC and low FAM fluorescence gather in upper left; 4G/4G samples with high FAM fluorescence and low VIC gather in lower right; 5G/4G samples with high FAM and high VIC fluorescence gather in upper right; and finally, NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence (Figure 3).


Fig 3. Typical graph for controls and samples on StepOne

Note! If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid and test should be repeated!

22. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

23. Technical Support

For technical support, contact us via

Phone +98-9936223241

Email: info@novingene.com

24. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124


Email: info@novingene.com

Website: www.novingene.com

25. References

- Teruya, J., 2021. Management of bleeding patients. Cham: Springer.
- Heiman, M., Gupta, S., Khan, S.S., Vaughan, D.E. and Shapiro, A.D., 2017. Complete plasminogen activator inhibitor 1 deficiency.
- Vaughan, D.E., 2005. PAI-1 and atherothrombosis. Journal of Thrombosis and Haemostasis, 3(8), pp.1879-1883.
- Zhang, Q., Jin, Y., Li, X., Peng, X., Peng, N., Song, J. and Xu, M., 2020. Plasminogen activator inhibitor-1 (PAI-1) 4G/5G promoter polymorphisms and risk of venous thromboembolism—a meta-analysis and systematic review. Vasa.
- Mackay, Ian M., 2004. "Real-time PCR in the microbiology laboratory." Clinical microbiology and infection 10, no. 3: 190-212.

26. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com