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Kit Manual

Autumn 2025, Version 1.0

For Real-Time PCR Detection of PML-RARA (bcr1, bcr2, bcr3)
Transcripts
For Research Use Only

 24 (Cat# PML-RARA3XRQ24)

 48 (Cat# PML-RARA3XRQ48)

 96 (Cat# PML-RARA3XRQ96)

 NG-WI-ASL-65-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.

Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability.....	5
8. Product Use Limitations	5
9. Additionally Required Materials.....	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	7
13. RNA Isolation	7
14. cDNA Synthesis	7
15. PCR Protocol	7
16. Devices and software.....	8
17. Programming RotorGene	8
18. Programming of StepOne	9

19. Programming Other Machines	10
20. Data Analysis: RotorGene	10
21. Data Analysis: StepOne	13
22. Analytical Sensitivity	16
23. Disposal Method	17
24. Technical Support.....	17
25. Contact Information.....	17
26. References	17
27. Symbols	18

1. Introduction

PML-RARA 3X RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting 3 isoforms of PML-RARA transcripts named bcr1, bcr2 and bcr3. All required reagents are included in the PCR Mixes provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

Important Note: *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

2. Intended Use

PML-RARA 3X RQ kit is intended for detecting 3 isoforms of PML-RARA, bcr1, bcr2 and bcr3 transcripts. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC.

3. Background Information

PML-RARA is an abnormality resulted from t (15;17) (q22; q21) translocation. This translocation results in fusion of PML (promyelocytic) gene with RARA (retinoic acid receptor alpha) gene, and production of chimeric PML-RARA protein which is a transcription repressor and impairs the myeloid differentiation. While RARA breakpoints always occur in intron 2, PML breakpoints involve three different regions of intron 6 (55%), exon 6 (5%) and intron 3 (40%). These isoforms of PML-RARA are respectively called bcr1/Long/L, bcr2/Variant/V and bcr3/Short/S.

PML-RARA accounts for more than 90% of APL (acute progranulocytic leukemia) cases, 10-15% of AML (acute myeloid leukemia) cases.

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
bcr1 Mix*	Mix for bcr1	480 µl
bcr2 Mix*	Mix for bcr2	480 µl
bcr3 Mix*	Mix for bcr3	480 µl
ABL Mix*	Mix for ABL	480 µl
PML/ABL Pos Ctrl	Equal to 1,000 copies/µl bcr, and 10,000 copies/µl ABL	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2 ml PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant.

RNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used. For optimum results, a sample should include about 10 million WBC per 150 μ l.

Whole blood or buffy coat should be shipped and stored at +4°C for 48 hrs. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, and Mannheim, Germany).
- TRIzol isolation reagent (Cat. no. 15596026, Invitrogene/ Thermo Fisher, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime/ Thermo Fisher, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

14. cDNA Synthesis

about 4ug of total RNA is required and should be reverse transcribed to cDNA using random hexamers. Different kits are available on the market for this purpose.

Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20ul of cDNA add 30ul of nuclease free water.

Upon successful RNA extraction and cDNA synthesis, CT of ABL for sample should be less than CT of ABL Positive Control.

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin.

Each sample should be examined for bcr1,2,3 fusion gene and for ABL control gene. So, four sets of reactions are required for each sample. In each bcr set, consider one tube for each sample as well as one tube for Pos Control and one for Negative control or NTC. In ABL set, consider one tube for each sample and one tube for ABL Pos Control and one tube for Negative control or NTC. Place the required number of the tubes on cold block.

Pipette 20µl of bcr 1 Mix to the first series of tubes,

Pipette 20µl of bcr 2 Mix to the second series of tubes,

Pipette 20µl of bcr 3 Mix to the third series of tubes,

Pipette 20µl of ABL Mix to the fourth series of tubes.

Continue by adding 5µl of cDNA, and Controls to each tube.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.

16. Devices and software

PML-RARA 3X RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming RotorGene


Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the PML-RARA template file for RotorGene (provided in the flash card, or accessible by kit QR code); PML-RARA 3X 0.1 is for strip tubes and PML-RARA 3X 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation.

Adjust the setting according to the image.

Select a tube number containing bcr Mix for the Green channel and tube with ABL Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file.

Edit sample names.

Make sure in the "Type" column, Positive control have been defined as "Positive Control". Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

18. Programming of StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. Positive and Negative controls for bcr and ABL and a few samples are defined. You may change the plate setup using right-click options (copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished click on Start Run and save the experiment. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Both bcr Mix and ABL Mix contain ROX with a final concentration of 300nM in reaction.

20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu Patient samples should be defined as "unknown" and Negative control, or no template control as "Negative Control" or "NTC", respectively.

Analyze the data according to Rotor-Gene manual. Perform quantitative analysis for **bcr (the Green channel)** and **ABL (the Yellow channel)**. Briefly, click on analysis menu and then, under Quantitation tab double click on “Cycling A. Green” and Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold at 0.05. Figures 1 and 2 represent typical graphs for the RotorGene machine.


Fig 1. Typical Controls graph with bcr Mixes in Green channel for Rotor-Gene


Fig 2. Typical Controls graph with ABL Mix in Yellow channel for Rotor-Gene

Consider the following points when analyzing:

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT, if present, is not reliable.

- A sample is **Positive for bcr1** if it is positive in the Green channel with bcr1 Mix and the Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for Yellow.
- A sample is **Positive for bcr2** if it is positive in the Green channel with bcr2 Mix and the Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for Yellow.
- A sample is **Positive for bcr3** if it is positive in the Green channel with bcr3 Mix and the Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for Yellow.
- A sample is **Negative** if it is negative in the Green channel with bcr1, bcr2 and bcr3 Mixes while it is positive in the Yellow/ABL channel with CT of equal or less than Positive Control CT with ABL Mix.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green channel with the all Mixes bcr1, bcr2 and bcr3 and in the Yellow/ABL channel with ABL Mix. Improper extraction or test set up could cause that.
 - Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green/bcr channel while it is

positive in the Yellow/ABL channel with sigmoid graph and CT of above Positive Control CT with ABL Mix. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of less than Positive Control CT with ABL Mix. CTs above Positive Control CT in Yellow/ABL channel usually happens if not enough cells have been extracted or less than 1ug RNA has been used, high ABL CTs leads to reduce the sensitivity of the test and may result in **False negative** reports.*

- Note: It is possible for a sample to be positive with a sigmoid graph in the Green channel for all three mixes; bcr1, bcr2, and bcr3 and also positive with a sigmoid graph and a CT of 20-30 in the ABL/Yellow channel. In such cases, the sample is considered positive for the mix that shows the lowest CT in the Green channel, which is most often bcr1.
- Note: It is possible for a sample to be positive with a sigmoid graph in the Green channel for two mixes; bcr2 and bcr3 and also show a positive result with a sigmoid graph and a CT of 20-30 in the ABL/Yellow channel. In such cases, the sample is considered positive for the mix that shows the lowest CT in the Green channel, which is most often bcr2.

21. Data Analysis: StepOne

Analyze the data according to StepOne Manual. Briefly, click on “Analyze” and set the threshold on 0.05 for **bcr/FAM** and 0.1 for **ABL/VIC**. Figures 3 and 4 represent typical graphs for StepOne.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT if present, is not reliable.


Fig 3. Typical Controls graph with bcr Mixes in FAM channel for StepOne


Fig 4. Typical Controls graph with ABL Mix in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive for bcr1** if it is positive in the FAM channel with bcr1 Mix and in VIC/ABL channel with sigmoid graphs and CT of 20-40 for the FAM and 20-30 for the VIC.
- A sample is **Positive for bcr2** if it is positive in the FAM channel with bcr2 Mix and in VIC/ABL channel with sigmoid graphs and CT of 20-40 for the FAM and 20-30 for the VIC.
- A sample is **Positive for bcr3** if it is positive in the FAM channel with bcr3 Mix and in VIC/ABL channel with sigmoid graphs and CT of 20-40 for the FAM and 20-30 for the VIC.
- A sample is **Negative** if it is negative in the FAM channel with bcr1, bcr2 and bcr3 Mixes while it is positive in the VIC/ABL channel with CT of equal or less than Positive Control CT with ABL Mix.

- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM channel with all Mixes bcr1, bcr2 and bcr3 and in the VIC/ABL channel with ABL Mix. Improper extraction or error in test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM/bcr channel while it is positive in the VIC/ABL channel with sigmoid graph and CT of above Positive Control CT with ABL Mix. Improper extraction or cDNA synthesis and low RNA input could cause that.

*Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of less than ABL Positive Control. CTs above Positive Control CT with ABL Mix, usually happens if not enough cells have been extracted or less than 1ug RNA has been used, which leads to reduce the sensitivity of the test and may result in **False negative** reports.*

- Note: It is possible for a sample to be positive with a sigmoid graph in the FAM channel for all three mixes bcr1, bcr2, and bcr3 and also positive with a sigmoid graph and a CT of 20-30 in the ABL/VIC channel. In such cases, the sample is considered positive for the mix that shows the lowest CT in the FAM channel, which is most often bcr1.
- Note: It is possible for a sample to show positive results with a sigmoid graph in the FAM channel for two mixes bcr2 and bcr3 and also show a positive result with a sigmoid graph and a CT of 20-30 in the ABL/VIC channel. In such cases, the sample is considered positive for the mix that shows the lowest CT in the FAM channel, which is most often bcr2.

22. Analytical Sensitivity

Analytical sensitivity of the kit has been determined by examining dilution series of cloned cDNA and was estimated about 2

copies/ul or 0.02% to reach this sensitivity each sample should contain higher than 10,000 copies/ul ABL transcripts.

23. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

24. Technical Support

For technical support, contact us via
Phone: +98 993-6223241
Email: info@novingene.com

25. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124

Email: info@novingene.com

Website: www.novingene.com


26. References

- Bain, B.J., 2017. Leukaemia diagnosis. Hoboken, Nj John Wiley & Sons, Inc.
- De Braekeleer, E., Douet-Guilbert, N. and De Braekeleer, M., 2014. RARA fusion genes in acute promyelocytic leukemia: a review. Expert review of hematology, 7(3), pp.347-357.
- Liquori, A., Ibañez, M., Sargas, C., Sanz, M.Á., Barragán, E. and Cervera, J., 2020. Acute Promyelocytic Leukemia: A

Constellation of Molecular Events around a Single PML-RARA Fusion Gene. Cancers, 12(3).

- Mackay IM., 2004. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect.; 10 (3): 190 – 212.

27. Symbols

RUO	Research use only		Manufacturer		Consult instructions for use
LOT	Lot number		Content sufficient for <n> tests		Use-by date
REF	Catalogue number	SN	Serial number		Temperature limit

For more information and resources please visit our website; www.novingene.com