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
Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition
☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10Fl	15Fl	1	10
Yellow	1	5Fl	10Fl	1	10

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Parvovirus B19 RQ Kit Manual

Autumn 2025, Version 1.0

For Real-Time PCR Detection and Quantitation of Parvovirus B19
For Research Use Only

 24 (Cat# ParvoRQ24)

 48 (Cat# ParvoRQ48)

 96 (Cat# ParvoRQ96)

 NG-WI-ASL-55-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	5
8. Product Use Limitations	5
9. Additionally Required Materials.....	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	7
13. Internal Control (IC)	7
14. DNA Isolation	8
15. PCR Protocol	8
16. Devices and software.....	9
17. Programming Rotor-Gene.....	9
18. Programming StepOne	10
19. Programming Other Machines	11

20. Data Analysis: Rotor-Gene	11
21. Data Analysis: StepOne	14
22. Quantitation	16
23. Linear Range.....	16
24. Sensitivity.....	17
25. Disposal Method	17
26. Technical Support.....	17
27. Contact Information.....	17
28. References	17
29. Symbols.....	18

1. Introduction

Parvovirus B19 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying Parvovirus B19 DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

Parvovirus B19 RQ kit is intended for detecting and quantifying Parvovirus B19 virus DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Parvovirus B19 is a member of Parvoviridae family with a non-enveloped icosahedral capsid and single-stranded DNA genome. This virus causes Erythema infectious or Fifth Disease, which is more common in children and usually has mild symptoms like mild rash on the cheeks. The main transmission route is through respiratory droplets while, it is also transmitted to the fetus during pregnancy. Since fetal infection can lead to serious consequences including fetal loss, early diagnosis during pregnancy is important.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
PV-B19 Mix	PCR mix*	360 µl
PV1	Standard 1: 100,000 copies/µl	150 µl
PV2	Standard 2: 10,000 copies/µl	150 µl
PV3	Standard 3: 1,000 copies/µl	150 µl
PV4	Standard 4: 100 copies/µl	150 µl
PV5	Standard 5: 10 copies/µl	150 µl
Internal Ctrl	Internal control	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components **on ice** completely, mix by flickering followed by a quick spin and **store on crushed ice** after.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Peripheral blood plasma or serum should be collected in sterile condition in proper and sterile containers. Samples can be stored at +4°C for 48 hours or stored at -20°C for up to a few weeks. The minimum sample volume recommended for testing is 200µl of plasma/serum, which requires 0.5ml of whole blood.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the Parvovirus B19 RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the PV-B19 Mix.

To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If IC is added to PV-B19 Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC reagent should be added to each reaction. For example, for 10 PCR reactions, 10ul of the Internal control should be added to 150ul of PV-B19 Mix before it is added to the tubes.

In a successful DNA extraction and PCR test, the Internal Control should generate a CT of 28-32 in the Yellow Channel of Rotor-Gene and 28-34 in the VIC Channel of StepOne.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample plus one for each standard and one for the negative control.

If the IC is introduced during the extraction process, pipette 15µl of PV-B19 Mix to each PCR tube.

If the IC is added to the PV-B19 Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.

Then add 10ul of extracted DNA, standard, or water to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software


Parvovirus B19 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Parvovirus template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); PV-B19 0.1 is for strip tubes and PV-B19 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Parvovirus Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. Five standards, one negative control and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add/remove samples or change sample names on "Define Targets and

Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The PV-B19 Mix contains ROX. Final concentration of ROX in reaction is 300nM.

20. Data Analysis: Rotor-Gene

To analyze data briefly, click on the Analysis menu, and then under the Quantitation tab double click on Cycling A. Green. Manually put threshold at 0.1. Repeat the above for the Yellow Channel. Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.


Fig 1. Typical PV-B19 graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for PV-B19/FAM at 0.1 and IC/VIC at 0.05. Figures 3 and 4 represent typical graphs for the StepOne machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical PV-B19 graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/Parvovirus B19 channel with a sigmoid graph and a CT of less than 40.
- A sample is **Negative** if it is negative in the FAM/Parvovirus B19 channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/Parvovirus B19 and VIC/IC channels.

The interpretation of the results is summarized in the following table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

22. Quantitation

The kit provides five quantitation standards with defined titers to generate a standard curve for the quantification of samples, viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently, using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ μ l. To convert the result to copy/ml following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result (copy/\mu l)} \times \text{elution volume (\mu l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

23. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target, and showed to be linear in the range of 10,000,000 copies/ μ l to 10 copies/ μ l.

24. Sensitivity

The analytical detection limit of the kit was assessed with a dilution series of the cloned target and showed a limit of detection equal to 4 copies/ μ l.

25. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

26. Technical Support

For technical support, contact us via
Phone +98 993-6223241
Email: info@novingene.com

27. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393
+98-990 1813124

Email: info@novingene.com

Website: www.novingene.com

28. References

- Bale Jr, J.F., 2014. Measles, mumps, rubella, and human parvovirus B19 infections and neurologic disease. Handbook of clinical neurology, 121, pp.1345-1353.

- Kerr, J.R., 2006. Parvoviruses. London: Hodder Arnold ; New York.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Rogo, L.D., Mokhtari-Azad, T., Kabir, M.H. and Rezaei, F., 2014. Human parvovirus B19: a review. Acta Virol, 58(3), pp.199-213.
- Ugai, S., Aizawa, Y., Kanayama, T. and Saitoh, A., 2018. Parvovirus B19: a cause of sepsis like syndrome in an infant. Pediatrics, 141(6).

29. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com