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ایمیل: info@novingene.com

وبسایت: www.novingene.ir

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- Griffiths, C., Drews, S.J. and Marchant, D.J., 2017, 'Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment', *Clinical Microbiology Reviews*, 30(1), pp. 277–319.
- Ji, T., Liu, Z., Wang, G., Guo, X., Lai, C., Chen, H., Huang, S., Xia, S., Chen, B., Jia, H. and Chen, Y., 2020. Detection of COVID-19: A review of the current literature and future perspectives. *Biosensors and Bioelectronics*, 166, p.112455.
- Mackay IM. Real-time PCR in the microbiology laboratory. *Clin. Microbiol. Infect.* 2004; 10 (3): 190 – 212.
- Palese, P. and Kingsbury, D.W., 2011. *Genetics of Influenza Viruses*. Springer.
- Paules, C. and Subbarao, K., 2017, 'Influenza', *The Lancet*, 390(10095), pp. 697–708.
- Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022, *Influenza*. *Lancet* (London, England), 400(10353), 693–706.

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Resp I RQ Kit Manual

Autumn 2025, Version 1.0

For Real-Time RT-PCR Detection of SARS-CoV-2,
Human Influenza A/B and RSV (RNA)
For Research Use Only

 24 (Cat#RespIRQ24)

 48 (Cat#RespIRQ48)

 96 (Cat#RespIRQ96)

 NG-WI-ASL-61-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Internal Control (IC)	6
13. RNA Isolation	6
14. RT-PCR Protocol	7
15. Devices and software	7
16. Programming of the Rotor-Gene	7
17. Programming Other Machines	8
18. Data Analysis: Rotor-Gene	9
19. Sensitivity	13
20. Disposal Method	13

21. Technical Support.....	13
22. Contact Information.....	13
23. References	13
24. Symbols.....	14

1. Introduction

Resp I RQ kit provides a ready-to-use Real-Time RT-PCR test designed for detecting RNA of the following viruses, SARS-CoV-2, Influenza A virus/Influenza B virus, and RSV. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting RNase P as an Internal Control (IC).

This kit is intended for Research Use Only!

2. Intended Use

Resp I RQ kit is intended for detecting RNA of following viruses, SARS-CoV-2, Human Influenza A or B viruses (without differentiation between A and B), and RSV RNA. Detection is achieved using One-Step Multiplex Real-Time RT-PCR system with Rotor-Gene and MIC machine.

3. Background Information

Coronaviruses are enveloped viruses with single stranded RNA genome which lead to respiratory illnesses.

With the emergence of SARS (Sever Acute Respiratory Syndrome by SARS CoV), then MERS (Middle East Respiratory Syndrome by MERS CoV), and recently COVID-19 (by SARS CoV-2), they are considered major health problem issues.

Human influenza A and B are enveloped RNA viruses of the Orthomyxoviridae family. While both viruses cause seasonal flu epidemics, Influenza A virus is responsible for the majority of them and may also cause pandemics.

Respiratory Syncytial Virus (RSV) is also an enveloped RNA virus and belongs to Paramyxoviridae family. RSV mostly causes lower respiratory tract infections with cold like symptoms in healthy

individuals, but it may also cause severe respiratory infections in some.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
Resp A Mix	RT-PCR Master mix* for SARS-CoV-2, Flu A/Flu B, RSV	480 µl
Resp A Pos Ctrl	Positive Control	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit and required equipments/items
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

This kit does not require cDNA synthesis reagents!

10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**

- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted RNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep RT-PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

Proper sample to test is upper respiratory tract sample mainly nasopharyngeal swab. Samples can be stored at 2-8°C for a few hours or at -20°C or lower for up to a few months.

12. Internal Control (IC)

To assess the possibility of RNA extraction failure and PCR inhibition and to prevent false negative results, The Resp A Mix contains primers and probe of a housekeeping gene (RNase P) as an internal control (IC). In a successful RNA extraction and RT-PCR test, IC should generate a CT of 20-30 in Orange/Rox Channel.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)

- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

14. RT-PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for positive sample and one for water.

Aliquot 20µl of Resp A Mix directly to each PCR tube. Then add 5ul of extracted RNA, Positive Control or water.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

15. Devices and software


Resp I RQ kit is designed to work with Rotor-Gene and MIC.

16. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Resp I template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Resp I 0.1 is for strip tubes and Resp I 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Resp A Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

In this experiment, ROX should not be selected as the reference dye.

18. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

To analyze data briefly, click on the Analysis menu and then under Quantitation tab double click on Cycling A. Green. Manually put threshold, at 0.1. Repeat the above for Yellow, Red and Orange Channels.

Figures 1, 2, 3 and 4 represent typical graphs for Rotor-Gene.

To interpret the results, please note that:

A Signal in the **Green** channel indicates **SARS-CoV-2**, a signal in the **Yellow** channel is due to **Influenza A or B**, **Red** channel is due to **RSV** and **Orange** channel indicates **IC or RNase P**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** for SARS-CoV-2 if it is positive in the **Green** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for Influenza A or B if it is positive in the **Yellow** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for RSV if it is positive in the **Red** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Negative** for Influenza A, B, RSV and SARS-CoV-2 if it is negative in the **Green** and the **Yellow** and the **Red** channels while it is positive in the **Orange** channel with a sigmoid graph and CT of 20-30.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in **Orange** channel.

The interpretation of results is summarized in bellow Table.

Green	Yellow	Red	Orange	Result
+	-	-	+	Pos for SARS-CoV-2
-	+	-	+	Pos for Flu A/B
-	-	+	+	Pos for RSV
-	-	-	+	Negative
- / +	- / +	- / +	-	Inconclusive


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene


Fig 3. Typical Controls graph in Red channel for Rotor-Gene


Fig 4. Typical Controls graph in Orange channel for Rotor-Gene

19. Sensitivity

The analytical detection limit of the kit was assessed with the cloned targets genomes and showed a limit of detection equal to 13 copies/ μ l for SARS-CoV-2 and 10 copies/ μ l for Influenza A and B, and 9 copies/ μ l for RSV.

20. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: info@novingene.com





Website: www.novingene.com

23. References

- Griffiths, C., Drews, S.J. and Marchant, D.J., 2017 'Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment', Clinical Microbiology Reviews, 30(1), pp. 277–319.

- Ji, T., Liu, Z., Wang, G., Guo, X., Lai, C., Chen, H., Huang, S., Xia, S., Chen, B., Jia, H. and Chen, Y., 2020. Detection of COVID-19: A review of the current literature and future perspectives. Biosensors and Bioelectronics, 166, p.112455.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Palese, P. and Kingsbury, D.W. (2011). Genetics of Influenza Viruses. Springer.
- Paules, C. and Subbarao, K., 2017. 'Influenza', The Lancet, 390(10095), pp. 697–708.
- Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022. Influenza. Lancet (London, England), 400(10353), 693–706.

24. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com