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- COVID-19: A review of the current literature and future perspectives. Biosensors and Bioelectronics, 166, p.112455.
- Liu, W.-K., Liu, Q., Chen, D.-H., Liang, H.-X., Chen, X.-K., Huang, W.-B., Qin, S., Yang, Z.-F. and Zhou, R., 2013. Epidemiology and clinical presentation of the four human parainfluenza virus types. BMC Infectious Diseases, 13(1).
 - Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
 - Palese, P. and Kingsbury, D.W., 2011. Genetics of Influenza Viruses. Springer.
 - Paules, C. and Subbarao, K., 2017 'Influenza', The Lancet, 390 (10095), pp. 697–708.
 - Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022. Influenza. Lancet (London, England), 400 (10353), 693–706.

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Resp II RQ Kit Manual

Fall 2024, Version 1.0

For Real-Time RT-PCR Detection of SARS-CoV-2,
Human Influenza Virus A/B, Respiratory Syncytial Virus (RSV)
and Human Parainfluenza Viruses 1-4 (HPIV1-4) RNA
For use with Rotor-Gene
For Research Use Only

 24 (Cat# RespIIRQ24)

 48 (Cat# RespIIRQ48)

 96 (Cat# RespIIRQ96)

 NG-WI-ASL-63-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

| | |
|---|----|
| 1. Introduction | 3 |
| 2. Intended Use | 3 |
| 3. Background Information | 3 |
| 4. Test Principle | 4 |
| 5. Kit Contents | 4 |
| 6. Packaging models | 5 |
| 7. Storage and Stability | 5 |
| 8. Product Use Limitations | 5 |
| 9. Additionally Required Materials | 5 |
| 10. General Precautions | 6 |
| 11. Specimen, Storage and Transport | 6 |
| 12. Internal Control (IC) | 7 |
| 13. RNA Isolation | 7 |
| 14. RT-PCR Protocol | 7 |
| 15. Devices and software | 8 |
| 16. Programming of the Rotor-Gene | 8 |
| 17. Programming Other Machines | 9 |
| 18. Data Analysis: Rotor-Gene | 10 |
| 19. Sensitivity | 17 |
| 20. Disposal Method | 18 |

| | |
|------------------------------|----|
| 21. Technical Support..... | 18 |
| 22. Contact Information..... | 18 |
| 23. References | 18 |
| 24. Symbols..... | 19 |

1. Introduction

Resp II RQ kit provides a ready-to-use Real-Time RT-PCR test designed for detecting RNA of the following viruses, SARS-CoV-2, Influenza A virus/Influenza B virus, Respiratory Syncytial Virus (RSV) and Human Parainfluenza Viruses types 1 to 4 (HPIVs 1-4). All required reagents are included in the PCR Mixes provided in the kit. Resp A Mix also contains a different series of primers and probe for detecting RNase P as an Internal Control (IC).

This kit is intended for Research Use Only!

2. Intended Use

The Resp II RQ kit includes two distinct Mixes, Resp A Mix and Resp B Mix. Resp A Mix is designed for detecting RNA from SARS-CoV-2, Human Influenza A or B viruses (without distinguishing between A and B), RSV and RNase P gene as IC. Resp B Mix is intended for the detection of HPIVs 1 to 4. Mixes are used with a One-Step Multiplex Real-Time RT-PCR system with the Rotor-Gene machine.

3. Background Information

Coronaviruses are enveloped viruses characterized with single stranded RNA genome, which cause respiratory illnesses.

The emergence of SARS (Severe Acute Respiratory Syndrome by SARS CoV), MERS (Middle East Respiratory Syndrome by MERS CoV), and recently COVID-19 (by SARS CoV-2), has highlighted their significant impact on global health.

Human influenza A and B are enveloped RNA viruses of the Orthomyxoviridae family. While both viruses cause seasonal flu epidemics, Influenza A virus is responsible for the majority of them and may also cause pandemics.

Respiratory Syncytial Virus (RSV) is also an enveloped RNA virus and belongs to Paramyxoviridae family. RSV mostly causes lower respiratory tract infections with cold like symptoms in healthy individuals, but it may also cause severe respiratory infections in some.

Human parainfluenza viruses are single-stranded RNA viruses classified under the paramyxoviridae family. These viruses are divided into four distinct types HPIV-1, HPIV-2, HPIV-3 and HPIV-4, each associated with different clinical symptoms and varying degrees of severity. Signs and symptoms are similar to common cold, but it can also cause more serious infection especially in children.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual a flash card with Rotor-Gene templates and the following reagents:

| Label | Content | Quantity |
|-----------------|---|----------|
| Resp A Mix | RT-PCR mix* for SARS-CoV-2, Flu A/Flu B, RSV | 480 µl |
| Resp B Mix | RT-PCR mix* for Human Parainfluenza Viruses (HPIV1-4) | 480 µl |
| Resp A Pos Ctrl | Positive Control for Resp A | 150 µl |

| | | |
|-----------------|-----------------------------|--------|
| Resp B Pos Ctrl | Positive Control for Resp B | 150 µl |
| Water | PCR Grade Water | 200 µl |

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips

- RNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

This kit does not require cDNA synthesis reagents!

10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted RNA to the PCR tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Place kit components on ice till completely thawed, mix by flickering followed by a quick spin and store on crushed ice after.**
- **Keep RT-PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.**
- Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

Proper samples for testing SARS CoV-2 and RSV can be collected from upper and lower respiratory tract and for Influenza from upper respiratory tract, mainly nasopharyngeal swabs. These samples should be stored at 2-8°C for a few hours and at -20°C or lower for up to a few weeks.

12. Internal Control (IC)

To assess the possibility of RNA extraction failure and PCR inhibition and to prevent false negative results, the Resp A Mix contains primers and a probe targeting the housekeeping gene RNase P as an internal control (IC). In a successful RNA extraction and RT-PCR test, the reaction of sample with Resp A Mix should generate a CT of 20-30 in Orange/Rox Channel. So, all patient samples must be examined with Resp A Mix to control RNA extraction.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

14. RT-PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. This kit includes two mixes, [Resp A Mix](#) which detects SARS-CoV-2, Influenza A/B virus and RSV and [Resp B Mix](#) which detects HPIVs 1-4. Each sample should be examined with both mixes in two separate series of reactions. In each series, consider one tube for each sample plus one for positive control and one for water.

Aliquot 20µl of Resp A Mix to the first series of PCR tube followed by adding 5ul of extracted RNA, Positive Control A or water.

Aliquot 20µl of Resp B Mix to each tube of the second series followed by adding 5ul of extracted RNA, Positive Control B or water.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

15. Devices and software

Resp II RQ kit is designed for use with Rotor-Gene.

16. Programming of the Rotor-Gene


Before starting the machine, ensure you have attached the locking ring on the rotor!

Open the memory card provided in the kit and double-click on Respiratory template "Resp II 0.2" or "Resp II 0.1" depending on the tubes used.

Note: To adjust the gain in the software, go to the View menu. Then, select the Gain Optimisation. Under the Auto-Gain Optimisation Setup first, select the Optimise Acquiring option.

In the Channel Setting, set the Tube Position to number 1 for all channels (consider that Tube number 1 should contain Mix Resp A). and set the Target Sample range within the range of 10 to 15. Enable the option Perform Optimisation Before 1st Acquisition option and close the window.

Adjust the setting according to the bellow image.

Gain Optimisation Setting

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

| Step | Temperature and time | Cycles |
|------|----------------------|--------|
| 1 | 50°C x 10 min | 1 |
| 2 | 95°C x 3 min | 1 |
| 3 | 95°C x 15 sec | 45 |
| | 60°C x 60 sec | |

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

18. Data Analysis: Rotor-Gene

The results are interpreted based on the reaction of the samples with two mixes (Resp A Mix and Resp B Mix) in four channels. It is specified in the next table that the four channels of each mix are related to which virus.

| Mix/Channel | Green/FAM | Yellow/HEX | Orange/ROX | Red/Cy5 |
|-------------|------------|------------------|---------------|---------|
| Resp A Mix | SARS-CoV-2 | Influenza A or B | Internal Ctrl | RSV |
| Resp B Mix | HPIV-1 | HPIV-2 | HPIV-3 | HPIV-4 |

A) Analysis of the sample results with [Resp A Mix](#):

To analyze data of [Resp A Mix](#) briefly, click on the analysis menu and then under Quantitation tab double click on cycling A. Green. Manually put threshold, at 0.1. Repeat the above for Yellow, Red and Orange Channels.

Figures 1, 2, 3 and 4 represent typical graphs with Resp A Mix for Rotor-Gene.

To interpret the results, please note that:

In the [Resp A Mix](#), a signal detected in the **Green** channel indicates **SARS-CoV-2**, a signal in the **Yellow** channel is due to **Influenza A or B**, a signal in the **Red** channel represents **RSV**, and a signal in the **Orange** channel indicates **IC or RNase P**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical Controls graph with Resp A Mix in Green channel for Rotor-Gene (SARS-CoV-2)


Fig 2. Typical Controls graph with Resp A Mix in Yellow channel for Rotor-Gene (Influenza A or B)


Fig 3. Typical Controls graph with Resp A Mix in Red channel for Rotor-Gene (RSV)


Fig 4. Typical Controls graph with Resp A Mix in Orange channel for Rotor-Gene (Internal Control)

Consider the following points when analyzing results with [Resp A Mix](#):

- A sample is **Positive** for SARS-CoV-2 if it is positive in the **Green** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for Influenza A or B if it is positive in the **Yellow** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for RSV if it is positive in the **Red** channel with sigmoid graphs and CT of less than 40 and also positive in the **Orange** channel with a CT of 20-30.
- A sample is **Negative** for Influenza A, B, RSV and SARS-CoV-2 if it is negative in the **Green** and the **Yellow** and the **Red** channels while it is positive in the **Orange** channel with a sigmoid graph and CT of 20-30.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in **Orange** channel with Resp A Mix.

The interpretation of results is summarized in the following Table.

| Resp A Mix | | | | |
|------------|--------|------|--------|--------------------|
| Green | Yellow | Red | Orange | Result |
| + | - | - | + | Pos for SARS-CoV-2 |
| - | + | - | + | Pos for Flu A/B |
| - | - | + | + | Pos for RSV |
| - | - | - | + | Negative |
| - /+ | - /+ | - /+ | - | Inconclusive |

B) Analysis of the sample results with Resp B Mix:

To analyze data of **Resp B Mix** briefly, click on the analysis menu and then under Quantitation tab double click on cycling A. Green. Manually put threshold, at 0.1 for Green, Yellow and Red. Repeat the above for Orange Channel and put threshold at 0.05.

Figures 5, 6, 7 and 8 represent typical graphs with Resp B Mix for Rotor-Gene.

Signal with **Resp B Mix** in the **Green** channel indicates **HPIV-1**, a signal in the **Yellow** channel is due to **HPIV-2**, **Orange** channel indicates **HPIV-3** and **Red** channel is due to **HPIV-4**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 5. Typical Controls graph with Resp B Mix in Green channel for Rotor-Gene (HPIV-1)


Fig 6. Typical Controls graph with Resp B in Yellow channel for Rotor-Gene (HPIV-2)


Fig 7. Typical Controls graph with Resp B Mix in Orange channel for Rotor-Gene (HPIV-3)


Fig 8. Typical Controls graph with Resp B in Red channel for Rotor-Gene (HPIV-4)

Consider the following points when analyzing with **Resp B Mix:**

- A sample is **Positive** for HPIV-1 if it is positive in the **Green** channel with sigmoid graphs and CT of less than 40 and also positive for IC with Resp A Mix in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for HPIV-2 if it is positive in the **Yellow** channel with sigmoid graphs and CT of less than 40 and also positive for IC with Resp A Mix in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for HPIV-3 if it is positive with Resp B Mix in the **Orange** channel with sigmoid graphs and CT of less than 40 and also positive for IC with Resp A Mix in the **Orange** channel with a CT of 20-30.
- A sample is **Positive** for HPIV-4 if it is positive in the **Red** channel with sigmoid graphs and CT of less than 40 and

also positive for IC with Resp A Mix in the **Orange** channel with a CT of 20-30.

- A sample is **Negative** for Human Para Influenza viruses types 1 to 4 if it is negative in the **Green** and the **Yellow** and the **Red** channels while IC is positive with Resp A Mix in the **Orange** channel with a sigmoid graph and CT of 20-30.
- Results are **Inconclusive** and the test should be repeated if a sample IC is negative with Resp A Mix in the **Orange** channel.

The interpretation of results is summarized in the following Table.

| Resp B Mix | | | | Result |
|--|--------|-----|--------|---|
| Green | Yellow | Red | Orange | |
| + | - | - | - | Pos for HPIV-1 |
| - | + | - | - | Pos for HPIV-2 |
| - | - | - | + | Pos for HPIV-3 |
| - | - | + | - | Pos for HPIV-4 |
| - | - | - | - | Negative (if IC is Pos with Resp A in Orange channel) |
| If sample is Neg in Orange channel with Resp A Mix | | | | Inconclusive |

19. Sensitivity

The analytical detection limit of the kit was assessed with the cloned targets genomes (partial) and showed a limit of detection equal to 13 copies/μl for SARS-CoV-2 and 10 copies/μl for Influenza A and B, 9 copies/μl for RSV, 8 copies/μl for HPIV-1, 8.4 copies/μl for HPIV-2, 13.6 copies/μl for HPIV-3 and 9.2 copies/μl for HPIV-4.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: info@novingene.com

Website: www.novingene.com


23. References

- Griffiths, C., Drews, S.J. and Marchant, D.J., 2017 'Respiratory syncytial virus: Infection, detection, and new options for prevention and treatment', *Clinical Microbiology Reviews*, 30(1), pp. 277–319.
- Ji, T., Liu, Z., Wang, G., Guo, X., Lai, C., Chen, H., Huang, S., Xia, S., Chen, B., Jia, H. and Chen, Y., 2020. Detection of COVID-19: A review of the current literature and future perspectives. *Biosensors and Bioelectronics*, 166, p.112455.
- Liu, W.-K., Liu, Q., Chen, D.-H., Liang, H.-X., Chen, X.-K., Huang, W.-B., Qin, S., Yang, Z.-F. and Zhou, R., 2013.

Epidemiology and clinical presentation of the four human parainfluenza virus types. BMC Infectious Diseases, 13(1).

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Palese, P. and Kingsbury, D.W., 2011. Genetics of Influenza Viruses. Springer.
- Paules, C. and Subbarao, K., 2017 'Influenza', The Lancet, 390 (10095), pp. 697–708.
- Uyeki, T. M., Hui, D. S., Zambon, M., Wentworth, D. E., & Monto, A. S., 2022. Influenza. Lancet (London, England), 400 (10353), 693–706.

24. Symbols

| | | |
|------------------------------|--|--|
| RUO Research use only |  Manufacturer |  Consult instructions for use |
| LOT Lot number |  Content sufficient for <n> tests |  Use-by date |
| REF Catalogue number | SN Serial number |  Temperature limit |

For more information and resources please visit our website; www.novingene.com