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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Agyeman, A.A. and Ofori-Asenso, R., 2017. Tuberculosis—an overview. *Journal of Public Health and Emergency*, 1(7), pp.1-11 .
- Gagneux, S., 2018. Ecology and evolution of *Mycobacterium tuberculosis*. *Nature Reviews Microbiology*, 16(4), pp.202-213.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212 .
- Pieters, J. and McKinney, J.D., 2013. Pathogenesis of *Mycobacterium tuberculosis* and its Interaction with the Host Organism. Springer Science & Business Media.

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# TB RQ Kit Manual

Autumn 2025, Version 3.5

For Real-Time PCR Detection and Quantitation of Mycobacterium  
Tuberculosis Complex  
For Research Use Only

 24 (Cat# TBRQ24)

 48 (Cat# TBRQ48)

 96 (Cat# TBRQ96)


 NG-WI-ASL-10-305

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	4
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	4
8. Product Use Limitations .....	5
9. Additionally Required Materials .....	5
10. General Precautions .....	5
11. Proper Specimen .....	6
12. Internal control (IC) .....	6
13. DNA Isolation .....	7
14. PCR Protocol .....	7
15. Devices and software .....	8
16. Programming Rotor-Gene .....	8
17. Programming StepOne .....	9
18. Programming Other Machines .....	10
19. Data Analysis: Rotor-Gene .....	10

20. Data Analysis: StepOne .....	12
21. Quantitation .....	14
22. Linear Range .....	15
23. Sensitivity.....	15
24. Disposal Method .....	15
25. Technical Support.....	15
26. Contact Information.....	15
27. References .....	16
28. Symbols .....	16

## 1. Introduction

TB RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying Mycobacterium Tuberculosis Complex DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

## 2. Intended Use

TB RQ kit is intended for detecting and quantifying Mycobacterium tuberculosis complex DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

## 3. Background Information

Tuberculosis (TB) is one of the major public health problems worldwide and results in about 2 million deaths annually. This infection is caused by a group of closely related species and subspecies of intracellular pathogens called Mycobacterium tuberculosis complex (MTC) including *M. tuberculosis*, *M. bovis* and *M. bovis* BCG.

While bacterial culturing is still the gold standard for TB diagnosis, Real-Time PCR is increasingly used as an alternative due to its rapid, sensitive, and specific results. It also provides the highest sensitivity and widest dynamic range among other methods.

#### 4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

#### 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
TB Mix	PCR mix*	360 µl
TB S1	Standard 1: 100,000 copies/µl	150 µl
TB S2	Standard 2: 10,000 copies/µl	150 µl
TB S3	Standard 3: 1,000 copies/µl	150 µl
TB S4	Standard 4: 100 copies/µl	150 µl
TB S5	Standard 5: 10 copies/µl	150 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

#### 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

#### 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaws more than 3 times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) application.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- Mycobacterium/TB DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- **Decontaminate samples prior to DNA extraction.**

- Within the pre-PCR work area, assigns three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the TB Mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

## 11. Proper Specimen

Sputum, BAL (bronchoalveolar lavage), bronchial secretion, CSF and peritoneal puncture are proper patient samples. Samples can be stored at 2-8°C for a few hours or at -20°C or lower for up to a few days.

## 12. Internal control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the TB RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the TB Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and the patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

**Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its**

efficiency.

If the IC is added to TB Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 28-34 in the Yellow/VIC Channel.

### 13. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. Make sure the kit is especially manufactured for Mycobacterium/TB DNA extraction. If not available, you may use:

- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany) with some modification.

For more details, please refer to "artus® M. tuberculosis RG PCR Kit Handbook (Mat# 1046960EN, Qiagen)".

To monitor DNA extraction process, internal control should be applied to the extraction process. For more details, please, refer to section 12 of this handbook.

### 14. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, five for standards and one for the negative control.

**If the IC is introduced during the extraction process, pipette 15ul of TB Mix directly to each PCR tube.**

**If the IC is added to the TB Mix, add 15ul of the prepared mix (as described in section 12) to each PCR tube.**

**Then add 10ul of extracted DNA, standard, or water to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## **15. Devices and software**

TB RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

## **16. Programming Rotor-Gene**

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the TB template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); TB 0.1 is for strip tubes and TB 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain TB Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

**Auto-Gain Optimisation Setup**

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to  degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10FI	15FI	1	10
Yellow	1	5FI	10FI	1	10

## 17. Programming StepOne

Open the StepOne software (V 2.\*). in the Set-Up menu, Open the TB template file for StepOne (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, five standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment to the desired location. The instrument will start shortly.

## 18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	40
	60°C x 60 sec	


Fluorescence should be collected at 60°C for FAM and VIC dyes. The TB Mix contains ROX. The final reaction concentration of ROX is 300nM.

## 19. Data Analysis: Rotor-Gene


Before analyzing results, make sure that in the sample menu, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples are defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC”, respectively. the data according to the Rotor-Gene manual. Perform the quantitative analysis for the **TB (Green channel)** and qualitative analysis for the **Internal Control (Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on Cycling A. Green.

In the pop-up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then, click on OK. Repeat the above for “Cycling A. Yellow” but, cancel the Automatic Threshold and manually put threshold on 0.1. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine. **Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 1.** Typical TB graph in Green channel for Rotor-Gene

**Fig 2.** Typical IC graph in Yellow channel for Rotor-Gene

Consider the following points when analyzing:


- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 35. The titer and quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.

The interpretation of results is summarized in the table on page 14.


## 20. Data Analysis: StepOne

Analyze obtained data according to StepOne manual. Briefly, click on “Analyze” and set the threshold for the **TB/FAM** at 0.1 and at 0.1 for the **IC/VIC**. Figures 3 and 4 represent typical graphs for the StepOne machine.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 3.** Typical TB graph in FAM channel for StepOne

**Fig 4.** Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/TB channel with a sigmoid graph and a CT of less than 35. The titer and quantitation results are valid.
- A sample is **Negative** if it is negative in the FAM/TB channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/TB and VIC/IC channels.

Interpretation of results is summarized in the following table.

Green	Yellow	Result
+	+/-	Positive
-	+ (CT: 28-34)	Negative
-	-	Inconclusive

## 21. Quantitation

The kit provides five quantitation standards with defined titers to generate a standard curve for the quantification of samples load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently, using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ $\mu$ l. To convert the result to copy/ml the following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result (copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

## **22. Linear Range**

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copies/μl to 10 copies/μl.

## **23. Sensitivity**

The analytical detection limit of the kit was assessed with a cloned target and showed a limit of detection equal to 1 copy/μl.

## **24. Disposal Method**

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## **25. Technical Support**

For technical support, contact us via  
Phone: +98-9936223241  
Email: [info@novingene.com](mailto:info@novingene.com)

## **26. Contact Information**

### **NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990-1813124





Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)

## 27. References

- Agyeman, A.A. and Ofori-Asenso, R., 2017. Tuberculosis—an overview. Journal of Public Health and Emergency, 1(7), pp.1-11 .
- Gagneux, S., 2018. Ecology and evolution of Mycobacterium tuberculosis. Nature Reviews Microbiology, 16(4), pp.202-213.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212 .
- Pieters, J. and McKinney, J.D., 2013. Pathogenesis of Mycobacterium tuberculosis and its Interaction with the Host Organism. Springer Science & Business Media.

## 28. Symbols

 <b>RUO</b> Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b> Lot number	 <b>Content sufficient for &lt;n&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b> Catalogue number	 <b>SN</b> Serial number	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**