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
Auto-Gain Optimisation Setup

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 58 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	10Fl	15Fl	1	10
Yellow	1	5Fl	10Fl	1	10
Orange	1	2Fl	5Fl	1	10
Red	1	5Fl	8Fl	1	10

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Thrombo4X Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of Factor II, Factor V Leiden,
MTHFR C677T and MTHFR A1298C mutations
For use with Rotor-Gene
For Research Use Only

 24 (Cat#T4X24)

 48 (Cat#T4X48)

 96 (Cat#T4X96)

 NG-WI-ASL-60-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use.....	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	5
6. Packaging models.....	5
7. Storage and Stability	5
8. Product Use Limitations	5
9. Additionally Required Materials	6
10. General Precautions	6
11. Specimen, Storage and Transport	7
12. InterferingSubstances.....	7
13. Internal Control.....	7
14. DNA Isolation	7
15. PCR Protocol	8
16. Devices and software	8
17. Programming Rotor-Gene	8
18. Programming Other Machines	10
19. Data Analysis: Rotor-Gene.....	10
20. Disposal Method	28

21. Technical Support.....	28
22. Contact Information	28
23. References.....	29
24. Symbols	29

1. Introduction

Thrombo4X kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting Factor II (G2021A), Factor V Leiden (G1691A), MTHFR (C677T) and MTHFR (A1298C). All required reagents are included in the PCR Mix provided in the kit.

This kit is intended for Research Use Only!

2. Intended Use

Thrombo4X kit is intended for the detection of Factor II (G2021A), Factor V Leiden (G1691A), MTHFR (C677T) and MTHFR (A1298C) mutations in human DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene and MIC machine.

3. Background Information

Thrombophilia (hypercoagulability) is an acquired or inherited disorder in which abnormality in blood coagulation increases the risk of thrombosis. Several point mutations are known to play role in thrombophilia including Factor II (Prothrombin), Factor V, Factor VIII, MTHFR, Fibrinogen, Plasminogen, Antithrombin III, Protein S and Protein C.

Factor V Leiden is an autosomal dominant inherited disorder of blood clotting and refers to specific G1691A substitution in Factor V gene (R506Q). About 3 to 8 percent of Caucasian population and 20 percent of patients with Venous thromboembolism (VTE) carry one copy of this mutation (heterozygous) and about 1 in 5000 and 0.1 percent of patients with VTE carry two copies of it (homozygous). Heterozygosity for factor V Leiden is associated with 7-fold increased risk of venous thromboembolism (VTE)

including deep venous thrombosis (DVT) and pulmonary embolism. The risk of VTE is increased 80 fold in Factor V Leiden homozygous. Factor V Leiden is also associated with a 2 to 3 fold increase in relative risk of pregnancy loss.

Prothrombin or factor II is the precursor for thrombin in the coagulation cascade. G20210A mutation in prothrombin gene is the second most prevalent genetic risk factor for venous thromboembolism (VTE) after Factor V Leiden. This mutation increases the risk of venous thrombosis about 2 to 7 times. It has also been linked to increased risk of pregnancy loss. Prevalence of this mutation ranges between 0.5 to 5% depending on the race and ethnic background.

MTHFR stands for Methylenetetrahydrofolate reductase which is a key regulatory enzyme in folate and homocysteine metabolism. Since, the C677T and A1298C variants show reduced activity (30-65% and 15-30% reduction respectively), these polymorphisms are among the most important genetic determinants of plasma total homocysteine.

Depending on the ethnic background up to 50% of population may carry the variant allele of MTHFR.

4. Test Principle

The target sequence is detected using PCR, where primers specific to the target sequence amplify it. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
F2-F5 Mix	PCR Mix* for Factor 2 and Factor 5	480 µl
MTHFR Mix	PCR Mix* for MTHFR A1298C and MTHFR C677T	480 µl
Thrmb MM Ctrl	Homozygous Positive Control	100 µl
Thrmb WM Ctrl	Heterozygous Positive Control	100 µl
Thrmb WW Ctrl	Negative Control	100 µl
Water	PCR Grade Water	200 µl

*1, 2 or 4 tubes for 24, 48 or 96 reaction kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.

- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Whole blood and peripheral blood (0.5ml) are the preferred samples. Whole blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for few days). For longer terms, sample should be aliquoted and stored at -20°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control

This kit detects four mutations in Factor II, Factor V and MTHFR gene. Since each person carries wild type, mutant or both alleles, it serves as both the target and the Internal Control for the assay. Any sample should always be positive for at least one of them. If a sample is negative for both alleles, then the test should be repeated.

14. DNA Isolation

DNA isolation can be performed using different kits from various

manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. This kit includes two mixes, **F2-F5 Mix** which detects Factor V and Factor II alleles and **MTHFR Mix** which detects MTHFR A1298C and MTHFR C677T. Each sample should be examined with both F2-F5 Mix and MTHFR Mix in two separate set of reactions. In each set, consider 1 tube for each sample plus four tubes for controls and NTC.

Pipette 20ul of F2-F5 Mix to the first series of tubes and 20ul of MTHFR Mix to each tube of second series followed by adding 5ul of controls or sample DNA and NTC.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software

Thrombo4X RQ kit is designed to work with Rotor-Gene and MIC.


17. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Thrombo4X template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Thrombo4X 0.1 is

for strip tubes and Thrombo4X 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the run file. Edit sample names. in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

18. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	57°C x 30 sec	

Fluorescence should be collected at 57°C for FAM, VIC, Rox and Cy5 dyes.

Please note that, Mixes does not contain ROX dye as normalizer!

19. Data Analysis: Rotor-Gene

Results are analyzed based on presence or absence of reaction in each channel with each mix.

Briefly with F2-F5 mix, Factor V Leiden wild type allele will have a reaction in VIC channel and the mutant allele in FAM channel. With the same mix, Factor II wild type allele will have a reaction in ROX channel and the mutant allele in CY5 channel.

Similar to the above, with MTHFR mix, MTHFR1298 wild type allele will show a reaction in VIC channel and the mutant allele in FAM channel. With the same mix, MTHFR677 wild type allele will show cross reaction in ROX channel and the mutant allele in CY5 channel. Results are summarized in the following table.

Mix	Target	FAM	VIC	CY5	ROX	Results
F2-F5 Mix	F5	+	-			F5 MM
		+	+			F5 WM
		-	+			F5 WW

	F2			+	-	F2 MM
				+	+	F2 WM
				-	+	F2 WW
MTHFR Mix	1298	+	-			1298 MM
		+	+			1298 WM
		-	+			1298 WW
	677			+	-	677 MM
				+	+	677 WM
				-	+	677 WW
		-	-	-	-	NTC

Based on above point, two methods can be used for data analysis which will be described afterwards.

A) Analysis with scatter graph

B) Analysis with endpoint fluorescence

A) Analysis with scatter graph

To analyze the results with scatter graph for Factor V:

Briefly, click on Analysis menu, select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Green and Yellow channels and then click on “Show” and select F2-F5.

Please note that the genotyping of the samples can only be determined if all three controls of the kits were used in the experiment.

The “Scatter Graph Window” is now shown. Each dot represents one sample. The vertical axis shows the Green fluorescence and belongs to M allele and horizontal axis is for the Yellow fluorescence represent W allele.

Therefore, samples are located in four regions of the graph. MM samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Neg Control as “None” (Fig 1). Then turn on patient’s sample and determine its genotype.


Fig 1. Factor V Leiden-Typical scatter graph for Controls with F2-F5 Mix for Rotor-Gene

A) To analyze the results with scatter graph for Factor II:

Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Red and Orange channels and then click on “Show” and select F2-F5.

Please note that the genotyping of the samples can only be determined if all three controls of the kits were used in the experiment.

The “Scatter Graph Window” is now shown. Each dot represents one sample. The vertical axis shows the Orange fluorescence

and belongs to W allele and horizontal axis is for the Red fluorescence and M allele.

Therefore, samples are located in four regions of the graph. MM samples with high Red and low Orange fluorescence gather in lower right; WW samples with low Red fluorescence and high Orange gather in upper left; WM samples with high Red and high Orange fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Red and low Orange fluorescence.

In order to determine genotypes of the samples, the regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 2). Then, turn on patient’s sample and determine its genotype.


Fig 2. Factor II-Typical scatter graph for Controls with F2-F5 Mix For Rotor-Gene

To analyze the results with scatter graph for MTHFR A1298C:

Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Green and Yellow channels and then click on “Show” and select MTHFR.

Please note that the genotyping of the samples can only be determined if all three controls of the kits were used in the experiment.

The “Scatter Graph Window” is now shown. Each dot represents one sample. The vertical axis shows the Green fluorescence and

M allele. horizontal axis is for the Yellow fluorescence and W allele.

Therefore, samples are located in four regions of the graph. MM samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 3). Then, turn on patient’s sample and determine its genotype.


Fig 3. MTHFR A1298C-Typical scatter graph for Controls with MTHFR Mix in Rotor-Gene

A) To analyze the results with scatter graph for MTHFR C677T:

Briefly, click on Analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Red and Orange channels and then click on “Show” and select MTHFR.

Please note that the genotyping of the samples can only be determined if all three controls of the kits were used in the experiment.

The “Scatter Graph Window” is now shown. Each dot represents one sample. The vertical axis shows the Orange fluorescence

and W allele. horizontal axis is for the Red fluorescence and belong to M allele. Note that M allele is detected in Red channel and W allele in Orange channel for MTHFR C677T.

Therefore, samples are located in four regions of the graph. MM samples with high Red and low Orange fluorescence gather in lower right; WW samples with low Red fluorescence and high Orange gather in upper left; WM samples with high Red and high Orange fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Red and low Orange fluorescence.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above-mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None” (Figure 4). Then, turn on patient’s sample and determine its genotype.


Fig 4. MTHFR C677T-Typical scatter graph for Controls with MTHFR Mix in Rotor-Gene

Note! If location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also, if a sample is located in between regions, it should be re-examined!

B) To analyze by the endpoint fluorescence:

To analysis the results of each genetic factor, each sample should be checked with the corresponding mix in two channels belonging to the Wild type and Mutant alleles. It is specified in the

bellow table that the four channels of each mix are related to the alleles of which genetic factors.

Mix/Channel	Green/FAM	Yellow/HEX	Red/Cy5	Orange/ROX
F2-F5 Mix	F5 Mutant	F5 Wild type	F2 Mutant	F2 Wild type
MTHFR	1298 Mutant	1298 Wild type	677 Mutant	677 Wild type

To assess Factor V Leiden with endpoint fluorescence:

To analyze data briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green and F2-F5, manually put threshold at 0.2. the Green signal is due to Mutant or M allele. Repeat the above for the Yellow Channel and F2-F5 and put the threshold at 0.2 which due to wild or W allele. Figures 5 and 6 represent typical graphs for the Rotor-Gene.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase.


Fig 5. Factor V Leiden-Typical Control graph with F2-F5 Mix in Green channel for Rotor-Gene


Fig 6. Factor V Leiden-Typical Control graph with F2-F5 Mix in Yellow channel for Rotor-Gene

Consider following points when analysis:

- If a sample is Pos in **Green** channel with F2-F5 Mix, and Negative in **Yellow** channel it is **Homozygote / MM for Factor V Leiden**.
- If a sample is Pos in both **Green** and **Yellow** channel with F2-F5 Mix, sample is **Heterozygote / WM for Factor V Leiden**.
- If a sample is Pos in **Yellow** and Negative in **Green** channel with F2-F5 Mix, sample is **Homozygote Wild type or WW**.
- If a sample is negative for both alleles, then results are **inconclusive**, and the test should be repeated. Improper extraction could cause that.

The interpretation of the results is summarized in the table of page 9.

To assess Factor II with Quantitation graph:

To analyze data briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Red and F2-F5, manually put threshold at 0.2. the Red signal is due to Mutant or M allele. Repeat the above for the Orange Channel and F2-F5 and put the threshold at 0.2. The Orange signal is due to wild or W allele. Figures 7 and 8 represent typical graphs for the Rotor-Gene.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase.


Fig 7. Factor II -Typical Control graph with F2-F5 Mix in Red channel for Rotor-Gene


Fig 8. Factor II -Typical Control graph with F2-F5 Mix in Orange channel for Rotor-Gene

- If a sample is Pos in **Red** channel with F2-F5 Mix, and Negative in **Orange** channel it is **Homozygote / MM for Factor II**.
- If a sample is Pos in both **Red** and **Orange** channel with F2-F5 Mix, sample is **Heterozygote / WM for Factor II**.
- If a sample is Pos in **Orange** and Negative in **Red** channel with F2-F5 Mix, sample is **Homozygote Wild type or WW Factor II**.
- If a sample is negative for both alleles, then results are **inconclusive**, and the test should be repeated. Improper extraction could cause that.

The interpretation of the results is summarized in the table of page 9.

To assess MTHFR A1298C with Quantitation graph:

To analyze data briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Green and MTHFR, manually put threshold at 0.1. the Green signal is due to Mutant or M allele. Repeat the above for the Yellow Channel and MTHFR and put the threshold at 0.1 which due to wild or W allele. Figures 9 and 10 represent typical graphs for the Rotor-Gene.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase.


Fig 9. MTHFR A1298C -Typical Control graph with MTHFR Mix in Green channel for Rotor-Gene


Fig 10. MTHFR A1298C -Typical Control graph with MTHFR Mix in Yellow channel for Rotor-Gene

- If a sample is Pos in **Green** channel with MTHFR Mix, and Negative in **Yellow** channel it is **Homozygote / MM for MTHFR A1298C**.
- If a sample is Pos in both **Green** and **Yellow** channel with MTHFR Mix, sample is **Heterozygote / WM for MTHFR A1298C**.
- If a sample is Pos in **Yellow** and Negative in **Green** channel with MTHFR Mix, sample is **Homozygote Wild type or WW**.
- If a sample is negative for both alleles, then results are **inconclusive**, and the test should be repeated. Improper extraction could cause that.

The interpretation of the results is summarized in the table of page 9.

To assess MTHFR C677T with Quantitation graph:

To analyze data briefly, click on the Analysis menu and then under the Quantitation tab double-click on Cycling A. Red and MTHFR, manually put threshold at 0.1. the Red signal is due to Mutant or M allele. Repeat the above for the Orange Channel and MTHFR and put the threshold at 0.1. The Orange signal is due to wild or W allele. Figures 11 and 12 represent typical graphs for the Rotor-Gene.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase.


Fig 11. MTHFR C677T -Typical Control graph with MTHFR Mix in Red channel for Rotor-Gene


Fig 12. MTHFR C677T - Typical Control graph with MTHFR Mix in Orange channel for Rotor-Gene

- If a sample is Pos in **Red** channel with MTHFR Mix, and Negative in **Orange** channel it is **Homozygote MM for MTHFR C677T**.
- If a sample is Pos in both **Red** and **Orange** channel MTHFR Mix, sample is **Heterozygote WM for MTHFR C677T**.
- If a sample is Pos in **Orange** and Negative in **Red** channel with MTHFR Mix, sample is **Homozygote Wild type or WW**.
- If a sample is negative for both alleles, then results are **inconclusive** and the test should be repeated. Improper extraction could cause that.

The interpretation of the results is summarized in the Table of page 9.

20. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98-9936223241

email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: info@novingene.com

Website: www.novingene.com

23. References

- Andrea Luigi Tranquilli, 2011. Thrombophilia. Rijeka: Intech.
- Christiansen, S.C., Cannegieter, S.C., Koster, T., Vandenbroucke, J.P. and Rosendaal, F.R., 2005. Thrombophilia, clinical factors, and recurrent venous thrombotic events. *Jama*, 293(19), pp.2352-2361.
- Khan, S. and Dickerman, J.D., 2006. Hereditary thrombophilia. *Thrombosis journal*, 4(1), pp.1-17.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3, 2004: 190-212.

24. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com