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Toxo RQ

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Auto-Gain Optimisation Setup [X]

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 Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

Channel Settings :

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	5FI	10FI	1	10
Yellow	1	5FI	10FI	1	10

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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Adem, D.M. and Ame, M.M., 2023. Toxoplasmosis and its Significance in Public Health: A Review. Journal of Biomedical and Biological Sciences, 2(1), pp.1-20.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. Clinical microbiology and infection, 10(3), pp.190-212.
- Smith, N.C., Goulart, C., Hayward, J.A., Kupz, A., Miller, C.M. and van Dooren, G.G., 2021. Control of human toxoplasmosis. International Journal for Parasitology, 51(2-3), pp.95-121.

- Weiss, L.M. and Dubey, J.P., 2009. Toxoplasmosis: A history of clinical observations. International journal for parasitology, 39(8), pp.895-901.

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Toxo RQ Kit Manual

Autumn 2025, Version 1.1

For Real-Time PCR Detection of Toxoplasma Gondii
For Research Use Only

 24 (Cat# ToxoRQ24)

 48 (Cat# ToxoRQ48)

 96 (Cat# ToxoRQ96)

 NG-WI-ASL-53-101

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	4
9. Additionally Required Materials	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. Internal Control (IC)	6
14. DNA Isolation	7
15. PCR Protocol	7
16. Devices and software	8
17. Programming of the Rotor-Gene	8
18. Programming StepOne	9
19. Programming Other Machines	10

20. Data Analysis: Rotor-Gene	10
21. Data Analysis: StepOne	12
22. Sensitivity	14
23. Disposal Method	14
24. Technical Support	15
25. Contact Information	15
26. References	15
27. Symbols	16

1. Introduction

Toxo RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting *Toxoplasma gondii* DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps. This kit is intended for Research Use Only!

2. Intended Use

Toxo RQ kit is intended for detecting *Toxoplasma gondii* in DNA samples extracted from human. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Toxoplasma gondii is a parasitic protozoan that causes toxoplasmosis in warm-blooded animals. Cats are the definitive host for *Toxoplasma gondii*, but several intermediate hosts also exist. It is usually transmitted from infected meat or cat feces. *Toxoplasma gondii* is asymptomatic in most humans but it causes various symptoms in immunocompromised patients, pregnant women, and congenitally infected children with toxoplasmosis.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence

provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
Toxo Mix	PCR mix	360 μ l*
Pos Ctrl	Positive Control	150 μ l
Internal Ctrl	Internal Control*	250 μ l*
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.

- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.**

- Keep PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

We recommend amniotic fluid, blood and CSF. Sample can be stored at 2-8°C for 48 hours or at -20°C or lower for up to a few days.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the Toxo RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the Toxo Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to Toxo Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC reagent should be added to each reaction. For example, for 10 PCR reactions, 15ul of the internal control should be added to 150ul of Toxo Mix before it is added to the tubes.

In a successful DNA extraction and PCR test, the IC should generate a CT of 28-32 in the Yellow Channel on RotorGene and 28-34 in the VIC channel on StepOne.

14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample plus one for positive control and one for negative control (water).

If the IC is introduced during the extraction process, pipette 15µl of the [Toxo Mix](#) to each PCR tube.

If the IC is added to the Toxo Mix, add 15ul of the [prepared mix](#) (as described in section 13) to each PCR tube.

Then add 10ul of isolated DNA, [Positive Control](#), or water to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software


Toxo RQ kit is designed to work with Rotor-Gene, StepOne and MIC.

17. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the Toxo template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); Toxo 0.1 is for strip tubes and Toxo 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Toxo Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on desired location.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click on Template (provided in the flash card, or accessible by kit QR code). Click on "Plate Setup". Controls and a few samples are defined. You may change plate setup using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment to the desired location. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

20. Data Analysis: Rotor-Gene

To analyze data briefly, click on Analysis menu and then under Quantitation tab double-click on Cycling A. Green. Manually put threshold at 0.1. Repeat the above for the Yellow and Orange Channels. Figures 1 and 2 represent typical graphs for the Rotor-Gene. To interpret the results, please note that:

Signal in the **Green** channels is due to **Toxoplasma** and a signal in the **Yellow channel** is due to **IC**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene

Consider the following points when analyzing:

- A sample is **Positive** for *Toxoplasma gondii* if it is positive in the Green channel with sigmoid graphs and CT of less than 40.
- A sample is **Negative** for *Toxoplasma gondii* if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green and Yellow channels.

21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for the **FAM** channel at 0.1 and 0.05 for the **VIC** channel. To interpret the results, please note that signal in the **FAM** channels is due to **Toxoplasma**, and a signal in the **VIC** channel is due to **IC**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present is not reliable.


Fig 3. Typical Control graph in FAM channel for StepOne


Fig 4. Typical IC graph in VIC channel for StepOne

Consider the following points when analyzing:

- A sample is **Positive** for *Toxoplasma gondii* if it is positive in the FAM channel with sigmoid graphs and CT of less than 40.
- A sample is **Negative** for *Toxoplasma gondii* if it is negative in the FAM channel while it is positive in the VIC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM and VIC channels.

The interpretation of the results is summarized in the following Table.

Green	Yellow	Result
+	+/-	Positive
-	+	Negative
-	-	Inconclusive

22. Sensitivity

The analytical detection limit of the kit was assessed with cloned target, of *Toxoplasma gondii* genome and showed a limit of detection equal to 4 copy/ μ l.

23. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

24. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

25. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124




Email: info@novingene.com

Website: www.novingene.com

26. References

- Adem, D.M. and Ame, M.M., 2023. Toxoplasmosis and its Significance in Public Health: A Review. *Journal of Biomedical and Biological Sciences*, 2(1), pp.1-20.
- Mackay, I.M., 2004. Real-time PCR in the microbiology laboratory. *Clinical microbiology and infection*, 10(3), pp.190-212.
- Smith, N.C., Goulart, C., Hayward, J.A., Kupz, A., Miller, C.M. and van Dooren, G.G., 2021. Control of human toxoplasmosis. *International Journal for Parasitology*, 51(2-3), pp.95-121.
- Weiss, L.M. and Dubey, J.P., 2009. Toxoplasmosis: A history of clinical observations. *International journal for parasitology*, 39(8), pp.895-901.

27. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com