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## Ureaplasma RG (V1.0)

**Auto-Gain Optimisation Setup**

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 Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to  degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

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## Ureaplasma RG (V1.0)

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

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## Ureaplasma RG (V1.0)

- Marovt, M., Keše, D., Kotar, T., Kmet, N., Miljković, J., Šoba, B. and Maticič, M., 2015. Ureaplasma parvum and Ureaplasma urealyticum detected with the same frequency among women with and without symptoms of urogenital tract infection. European Journal of Clinical Microbiology & Infectious Diseases, 34(6), pp.1237-1245.

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# Ureaplasma RG Kit Manual

Summer 2025, Version 1.0

For Multiplex Real-Time PCR Detection of *Ureaplasma parvum*  
and *Ureaplasma urealyticum* DNA

For Use with Rotor-Gene or MIC

For Research Use Only

 24 (Cat# UreaplasmaRQ24)

 48 (Cat# UreaplasmaRQ48)

 96 (Cat# UreaplasmaRQ96)


 NG-WI-ASL-67-100

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	4
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	4
8. Product Use Limitations .....	5
9. Additionally Required Materials .....	5
10. General Precautions .....	5
11. Specimen, Storage and Transport .....	6
12. Internal Control (IC) .....	6
13. DNA Isolation .....	7
14. PCR Protocol .....	7
15. Devices and software .....	8
16. Programming of the Rotor-Gene .....	8
17. Programming Other Machines .....	9
18. Data Analysis: Rotor-Gene .....	10
19. Sensitivity .....	13

20. Disposal Method .....	13
21. Technical Support.....	13
22. Contact Information.....	13
23. References .....	13
24. Symbols.....	14

## **1. Introduction**

Ureaplasma RG kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting DNA of *Ureaplasma parvum* and *Ureaplasma urealyticum*.

All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

## **2. Intended Use**

Ureaplasma RG kit is intended for detecting DNA of *Ureaplasma parvum* and *Ureaplasma urealyticum*. Detection is achieved using Real-Time PCR system with Rotor-Gene and MIC machine.

## **3. Background Information**

*Ureaplasma parvum* and *Ureaplasma urealyticum* are classified under Mycoplasmataceae. These cell wall-deficient bacteria are commonly found in the urogenital tract of many individuals. Although they are often asymptomatic, in some cases they are associated with complications such as infertility, genital tract inflammation, and pregnancy-related conditions including preterm birth and chorioamnionitis.

Since, culturing *Ureaplasma* species is both difficult and time-consuming, molecular diagnostic methods particularly PCR and Real-Time PCR are recognized as the most reliable approaches for their detection. These techniques provide rapid, sensitive, and accurate identification, ensuring precise differentiation.

#### 4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

#### 5. Kit Contents

The kit contains a manual, flash card, with Rotor-Gene templates and the following reagents:

Label	Content	Quantity
Ureaplasma Mix	PCR mix* for Ureaplasma parvum and Ureaplasma urealyticum	360 µl
STI Pos Ctrl	Positive Control for	150 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

#### 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

#### 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for in vitro diagnostics.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and

- c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
  - Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
  - Keep PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
  - Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

## **11. Specimen, Storage and Transport**

Proper samples to test are urogenital swabs (endocervical/urethral) and urine. In special clinical contexts amniotic fluid (in pregnant women) and some other samples if requested by the physician are appropriate samples. Samples can be stored at 2-8°C for 48 hours or at -20°C or lower for a few months.

## **12. Internal Control (IC)**

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the Ureaplasma RG kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the Ureaplasma Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to Ureaplasma Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the Ureaplasma Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 27-32 in the Yellow/VIC Channel.

### 13. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- QIAamp DNA Blood Mini Kit (Cat# 51104, Qiagen, Germany)
- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)

### 14. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus one for positive sample and one for water.

**If the IC is introduced during the extraction process, pipette 15ul of the [Ureaplasma Mix](#) to the PCR tube.**

**If the IC is added to the [Ureaplasma Mix](#), add 15ul of the [prepared mix](#) (as described in section 12) to the PCR tube.**

**Then add 10ul of extracted DNA, [Positive Control](#) or water.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: If using Rotor-Gene attach the locking ring.*

## 15. Devices and software

Ureaplasma RG kit is designed to work with Rotor-Gene and MIC.

## 16. Programming of the Rotor-Gene


*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the memory card provided in the kit and double-click on Ureaplasma template " Ureaplasma 0.1" or " Ureaplasma 0.2" depending on the tubes used.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the image.

Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Ureaplasma Mix).

Then close the Auto Gain Optimisation window.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

## 17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table.

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and dyes.


*Note! ROX should not be selected as reference dye in this test.*

## 18. Data Analysis: Rotor-Gene


To analyze data Ureaplasma briefly, click on the analysis menu and then under Quantitation tab double click on cycling A. Green. Manually put threshold, at 0.1. Repeat the above for Yellow and Orange Channels.

Figures 1, 2 and 3 represent typical graphs for Rotor-Gene. To interpret the results, please note that; a signal in the **Green** channel indicates **Ureaplasma Urealyticum**, a signal in **Orange** channel indicates **Ureaplasma parvum** and a signal in the **Yellow** channel indicates **Internal Control**.


**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 1.** Typical Controls graph in Green channel for Rotor-Gene (Ureaplasma urealyticum)

**Fig 2.** Typical Controls graph in Orange channel for Rotor-Gene (Ureaplasma parvum)

**Fig 3.** Typical Controls graph in Yellow channel for Rotor-Gene (Internal Control)

Consider the following points for analyzing results:

- A sample is **Positive** for *Ureaplasma urealyticum* if it is positive in the **Green** channel with sigmoid graphs and CT of less than 40 and in the **Yellow** channel with a CT of 20-32.
- A sample is **Positive** for *Ureaplasma parvum* if it is positive in the **Orange** channel with sigmoid graphs and CT of less than 40 and in the **Yellow** channel with a CT of 20-32.
- A sample is **Negative** for *Ureaplasma parvum*, and *Ureaplasma urealyticum* if it is negative in the **Green** and the **Orange** channels while it is positive in the **Yellow** channel with a sigmoid graph and CT of 20-32.
- Results are **Inconclusive** and the test should be repeated if a sample IC is negative in **Yellow** channel.

The interpretation of results is summarized in the following table.

Green	Orange	Yellow	Result
+	-	+	Pos for <i>Ureaplasma urealyticum</i>
-	+	+	Pos for <i>Ureaplasma parvum</i>
-	-	+	Negative
-	-	-	Inconclusive

## 19. Sensitivity

The analytical detection limit of the kit was assessed with a cloned targets genomes and showed a limit of detection equal to 7.5 copies/μl for *Ureaplasma urealyticum* and 7.7 copies/μl for *Ureaplasma parvum*.

## **20. Disposal Method**

The contents of the kit do not require any special method of disposal and can be directly discarded. Infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## **21. Technical Support**

For technical support, contact us via

Phone: +98 993-6223241

email: [info@novingene.com](mailto:info@novingene.com)

## **22. Contact Information**

### **NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: [info@novingene.com](mailto:info@novingene.com)

Website: [www.novingene.com](http://www.novingene.com)


## **23. References**

- Buder, S., Schöfer, H., Meyer, T., Bremer, V., Kohl, P.K., Skaletz-Rorowski, A. and Brockmeyer, N., 2019. Bacterial sexually transmitted infections. JDDG: Journal der Deutschen Dermatologischen Gesellschaft, 17(3), pp.287–315.
- Hunjak, B., Sabol, I., Vojnović, G., Fistončić, I., Erceg, A.B., Peršić, Z. and Grce, M., 2014. Ureaplasma urealyticum and Ureaplasma parvum in women of reproductive age. Archives of gynecology and obstetrics, 289(2), pp.407-412.
- Mackay IM., 2004. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect.; 10 (3): 190 – 212.

- Mahon, C.R. and Lehman, D.C., 2023. Textbook of Diagnostic Microbiology - E-Book. Elsevier Health Sciences.
- Marovt, M., Keše, D., Kotar, T., Kmet, N., Miljković, J., Šoba, B. and Matičić, M., 2015. Ureaplasma parvum and Ureaplasma urealyticum detected with the same frequency among women with and without symptoms of urogenital tract infection. European Journal of Clinical Microbiology & Infectious Diseases, 34(6), pp.1237-1245.

## 24. Symbols

<b>RUO</b> Research use only	 Manufacturer	 Consult instructions for use
<b>LOT</b> Lot number	 Content sufficient for <n> tests	 Use-by date
<b>REF</b> Catalogue number	<b>SN</b> Serial number	 Temperature limit

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**