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ایمیل: info@novingene.com

وبسایت: www.novingene.com

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- Cilloni, D., Renneville, A., Hermitte, F., Hills, R.K., Daly, S., Jovanovic, J.V., Gottardi, E., Fava, M., Schnittger, S., Weiss, T. and Izzo, B., 2009. Real-time quantitative polymerase chain reaction detection of minimal residual disease by standardized WT1 assay to enhance risk stratification in acute myeloid leukemia: a European LeukemiaNet study. *Journal of Clinical Oncology*, 27(31), pp.5195-5201.
- Hastie, N., 2016. The Wilms' tumor (WT1) gene: methods and protocols. New York, NY: Humana Press.
- Hochhaus, A., Baccarani, M., Silver, R.T., Schiffer, C., Apperley, J.F., Cervantes, F., Clark, R.E., Cortes, J.E., Deininger, M.W., Guilhot, F. and Hjorth-Hansen, H., 2020. European LeukemiaNet 2020 recommendations for treating chronic myeloid leukemia. *Leukemia*, 34(4), pp.966-984.

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Yang, L., Han, Y., Saurez Saiz, F. and Minden, M.D., 2007. A tumor suppressor and oncogene: the WT1 story. Leukemia, 21(5), pp.868-876.

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WT1 RQ Kit Manual

Autumn 2025, Version 1.2

For Real-Time PCR Quantitation of WT1 Transcripts
For Research Use Only

 24 (Cat# WT1RQ24)

 48 (Cat# WT1RQ48)

 96 (Cat# WT1RQ96)

 NG-WI-ASL-46-102

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	5
8. Product Use Limitations	5
9. Additionally Required Materials	5
10. General Precautions	6
11. Specimen, Storage and Transport	6
12. Interfering Substances	7
13. RNA Isolation	7
14. cDNA Synthesis	7
15. PCR Protocol	7
16. Devices and software	8
17. Programming RotorGene	8
18. Programming of StepOne	9

19. Programming Other Machines	10
20. Data Analysis: RotorGene.....	10
21. Data Analysis: StepOne	13
22. WT1 Normalization	15
23. Analytical Sensitivity	15
24. Disposal Method	15
25. Technical Support.....	16
26. Contact Information.....	16
27. References	16
28. Symbols	17

1. Introduction

WT1 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying WT1 Transcripts. All required reagents are included in the PCR Mix provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

Important Note: *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

2. Intended Use

WT1 RQ kit is intended for detecting and quantifying of WT1 expression in sample. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene, StepOne and MIC machines.

3. Background Information

Measurable Residual Disease (MRD) is an important factor for risk stratification and treatment planning for AML patients. MRD could be evaluated according to different parameters including AML related gene expression. According to current LeukemiaNet guidelines, when no other parameter is available, WT1 expression may be used as MRD marker (A Hochhaus, 2020, Leukemia, 34:966).

Clinical studies revealed, WT1 is expressed in normal hematopoietic cells. The average expression level of this gene has been reported up to 50 copies per 10^4 ABL copies in peripheral blood and up to 250 copies/ 10^4 ABL in bone marrow. (D Cilloni, 2009, Journal of Clinical Oncology 27:5195).

This kit provides a ready-to-use system for quantitation of WT1 and ABL transcripts.

4. Test Principle

The Translocation is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a flash card, and the following reagents:

Label	Content	Quantity
WT1 Mix	Mix for WT1*	480 µl
ABL Mix	Mix for ABL*	480 µl
WT1 S1	WT1 Standard 1: 100,000 copies/µl	150 µl
WT1 S2	WT1 Standard 2: 10,000 copies/µl	150 µl
WT1 S3	WT1 Standard 3: 1,000 copies/µl	150 µl
WT1 S4	WT1 Standard 4: 100 copies/µl	150 µl
WT1 S5	WT1 Standard 5: 10 copies/µl	150 µl
A1	ABL Standard 1: 100,000 copies/µl	150 µl
A2	ABL Standard 2: 10,000 copies/µl	150 µl
A3	ABL Standard 3: 1,000 copies/µl	150 µl
A4	ABL Standard 4: 100 copies/µl	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reactions kits.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease-free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- **Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 48 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used. For optimum results, a sample should include about 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for two days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, and Mannheim, Germany).
- TRIzol isolation reagent (Cat. no. 15596026, Invitrogene/Thermo Fisher, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime/Thermo Fisher, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

14. cDNA Synthesis

about 4ug of total RNA is required and should be reverse transcribed to cDNA using random hexamers. Different kits are available on the market for this purpose.

Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20ul of cDNA add 30ul of nuclease free water.

Upon successful RNA extraction and cDNA synthesis, ABL titer will reach above 20,000 copy/ul.

15. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

WT1 RQ (V1.2)

Each sample should be examined for both in WT1 fusion gene and for ABL control gene. So, two sets of reactions are required. In WT1 set, consider 1 tube for each sample as well as 6 tubes for the 5 standards (WT1 S1 to WT1 S5) and Negative control or NTC. In ABL set, consider 1 tube for each sample and 5 tubes for the 4 standards (A1 to A4) and Negative control or NTC. Place required number of tubes on cold block.

Pipette 20µl of WT1 Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA, or standard, or control to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using RotorGene attach the locking ring.

16. Devices and software

WT1 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.


17. Programming RotorGene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the WT1 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); WT1 0.1 is for strip tubes and WT1 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Select a tube number containing WT1 Mix for the Green channel and tube with ABL Mix for the Yellow channel.

WT1 RQ (v1.2)

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Edit sample names on both WT1 and ABL pages. Remember that tubes containing WT1 Mix should only be named in WT1 page and tubes containing ABL Mix should only be named in ABL page.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

18. Programming of StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click Template (provided in the flash card, or accessible by kit QR code). Click on Plate Setup. One negative control, five standards for WT1, four standards for ABL and a few samples are defined. You may change plate set up using right-click options (copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished click on Start Run and save the experiment. Instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Both of WT1 Mix and ABL Mix contain ROX with final concentration of 300nM in reaction.

20. Data Analysis: RotorGene

Before analyzing results, make sure that in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as “unknown” and Negative control or no template control as “Negative Control” or “NTC”, respectively.

Analyze the data according to RotorGene manual. Perform quantitative analysis for both **WT1 (the Green channel)** and **ABL**

WT1 RQ (V1.2)

(the Yellow channel). Briefly, click on analysis menu and then, under Quantitation tab double-click on “Cycling A. Green”. Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1. Figures 1 and 2 represent typical graphs for the RotorGene machine.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 1. Typical WT1 graph in Green channel for Rotor-Gene


Fig 2. Typical ABL graph in Yellow channel for Rotor-Gene

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the Green/WT1 and Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the Green/WT1 channel while it is positive in the Yellow/ABL channel with a sigmoid graph and CT of 20-27 and ABL titer above 20,000 copies/ul.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the Green/WT1 and Yellow/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the Green/WT1 channel while it is positive in the Yellow/ABL channel with sigmoid graph and CT of above 27 and ABL titer less than 20,000 copies/ul. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 1ug RNA has been used. CTs higher than 27 and ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **alse negative** reports.*

21. Data Analysis: StepOne

Analyze the data according to StepOne manual. Briefly, click on “Analyze” and set the threshold for the **WT1/FAM** and the **ABL/VIC** on 0.1. Figures 3 and 4 represent typical graphs for the StepOne machine.


Fig 3. Typical WT1 graph in FAM channel for StepOne


Fig 4. Typical ABL graph in VIC channel for StepOne

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the FAM/WT1 and VIC/ABL channels with sigmoid graphs and CT of 20-40 for the FAM and 20-30 for the VIC.
- A sample is **Negative** if it is negative in the FAM/WT1 channel while it is positive in the VIC/ABL channel with a sigmoid graph and CT of 20-27 with ABL titer above 20,000 copies/ul.

- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/WT1 and VIC/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM/WT1 channel while it is positive in the VIC/ABL channel with sigmoid graph and CT of above 27 and ABL titer less than 20,000 copies/ul. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 1ug RNA has been used. CTs higher than 27 or ABL titer less than 20,000 copies/ul reduce the sensitivity of the test and may result in **false negative** reports.*

22. WT1 Normalization

To assess the response to therapy, WT1 value for each patient can be calculated with NCN method (D Cilloni, 2009, Journal of Clinical Oncology 27:5195).

In this method, the WT1 transcripts (titer) normalized by the ABL transcripts (titer) and then multiplied by 10,000 copies of ABL.

23. Analytical Sensitivity

Analytical sensitivity of the kit has been determined by examining dilution series of cloned cDNA and was estimated about 2 copies/ul.

24. Disposal Method

The contents of the kit do not require any special treatment before

disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

25. Technical Support

For technical support, contact us via:

Phone: +98 993-6223241

Email: info@novingene.com

26. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Website: www.novingene.com


27. References

- Cilloni, D., Renneville, A., Hermitte, F., Hills, R.K., Daly, S., Jovanovic, J.V., Gottardi, E., Fava, M., Schnittger, S., Weiss, T. and Izzo, B., 2009. Real-time quantitative polymerase chain reaction detection of minimal residual disease by standardized WT1 assay to enhance risk stratification in acute myeloid leukemia: a European LeukemiaNet study. *Journal of Clinical Oncology*, 27(31), pp.5195-5201.
- Hastie, N., 2016. The Wilms' tumor (WT1) gene: methods and protocols. New York, NY: Humana Press.
- Hochhaus, A., Baccarani, M., Silver, R.T., Schiffer, C., Apperley, J.F., Cervantes, F., Clark, R.E., Cortes, J.E., Deininger, M.W., Guilhot, F. and Hjorth-Hansen, H., 2020.

European LeukemiaNet 2020 recommendations for treating chronic myeloid leukemia. Leukemia, 34(4), pp.966-984.

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.
- Yang, L., Han, Y., Saurez Saiz, F. and Minden, M.D., 2007. A tumor suppressor and oncogene: the WT1 story. Leukemia, 21(5), pp.868-876.

28. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com