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**Auto-Gain Optimisation Setup**

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Auto-Gain Optimisation will read the fluorescence on the inserted sample at different gain levels until it finds one at which the fluorescence levels are acceptable. The range of fluorescence you are looking for depends on the chemistry you are performing.

Set temperature to  degrees.

☒ Perform Optimisation Before 1st Acquisition

☐ Perform Optimisation At 60 Degrees At Beginning Of Run

**Channel Settings :**

Name	Tube Position	Min Reading	Max Reading	Min Gain	Max Gain
Green	1	5FI	10FI	-10	10
Yellow	20	10FI	15FI	-10	10

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Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

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## **Kit Manual**

Autumn 2025, Version 4.0

For Real-Time PCR Detection and Quantitation of  
BCR-ABL Transcripts (p190)  
For Research Use Only

 24 (Cat# MBCR190RQ24)

 48 (Cat# MBCR190RQ48)

 96 (Cat# MBCR190RQ96)


 NG-WI-ASL-12-400

**RUO**

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	4
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability.....	5
8. Product Use Limitations .....	5
9. Additionally Required Materials.....	5
10. General Precautions .....	6
11. Specimen, storage and transport.....	6
12. Interfering Substances .....	7
13. RNA isolation .....	7
14. cDNA synthesis .....	7
15. PCR Protocol .....	7
16. Devices and software.....	8
17. Programming Rotor-Gene.....	8
18. Programming StepOne .....	10
19. Programming Other Machines .....	10

20. Data Analysis: Rotor-Gene .....	10
21. Data Analysis: StepOne .....	13
22. BCR-ABL% calculation .....	15
23. Analytical Sensitivity .....	16
24. Disposal Method .....	16
25. Technical Support.....	16
26. Contact Information.....	16
27. References .....	17
28. Symbols .....	17

## 1. Introduction

mbcr 190 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying BCR-ABL (p190) Transcripts. All required reagents are included in the PCR Mix provided in the kit. The kit also contains ABL Mix for the detection of *ABL* gene transcripts to prevent false negative results due to failure in extraction.

This kit is intended for Research Use Only!

**Important Note:** *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

## 2. Intended Use

mbcr 190 RQ kit is intended for the detection and quantitation of BCR-ABL (p190) transcripts (p190, e1a2 break point only) in peripheral blood and BCR-ABL% calculation in patients. This kit is designed for use with Rotor-Gene, StepOne and MIC machines.

## 3. Background Information

Philadelphia chromosome is an abnormality resulted from 9;22 translocation. Consequently, *ABL* proto-oncogene on chromosome 9 is fused with *BCR* gene on chromosome 22. This fusion produces BCR-ABL protein, mostly 210kDa (b2a2 or b3a2) or 190kDa (e1a2), with constitutively active tyrosine kinase activity promoting cell proliferation and inhibition of apoptosis. The fusion gene transcript is detectable in about 95% of CML patients and some cases of ALL. Also, serial monitoring of patients for identifying and measuring BCR-ABL transcripts provides more precise assessment of response to specific therapies and prediction of those in higher risk of disease progression.

#### 4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

#### 5. Kit Contents

The kit contains a manual, a flash card and the following reagents:

Label	Content	Quantity
mbcr 190 Mix*	Mix for BCR-ABL	480 $\mu$ l
ABL Mix*	Mix for ABL	480 $\mu$ l
m1	mbcr Standard 1: 100,000 copies/ $\mu$ l	150 $\mu$ l
m2	mbcr Standard 2: 10,000 copies / $\mu$ l	150 $\mu$ l
m3	mbcr Standard 3: 1,000 copies/ $\mu$ l	150 $\mu$ l
m4	mbcr Standard 4: 100 copies/ $\mu$ l	150 $\mu$ l
m5	mbcr Standard 5: 10 copies/ $\mu$ l	150 $\mu$ l
A1	ABL Standard 1: 100,000 copies/ $\mu$ l	150 $\mu$ l
A2	ABL Standard 2: 10,000 copies/ $\mu$ l	150 $\mu$ l
A3	ABL Standard 3: 1,000 copies/ $\mu$ l	150 $\mu$ l
A4	ABL Standard 4: 100 copies/ $\mu$ l	150 $\mu$ l
Water	PCR Grade Water	200 $\mu$ l

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

#### 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit
- cDNA synthesis kit
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Take utmost care to avoid RNase contamination during RNA extraction and cDNA synthesis.
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes, and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Always keep PCR Mix tube at -20°C. Take it out just before use and return it to freezer immediately after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

## 11. Specimen, storage and transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 48 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity buffy coat can be used. For optimum results, a sample should include 5 to 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise, should be aliquoted and stored at -70°C which is stable for a few months.

## 12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 13. RNA isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using:

- TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, and Mannheim, Germany).
- TRIzol isolation reagent (Cat. no. 15596026, Invitrogene/Thermo Fisher, USA)
- Isol-RNA isolation reagent (Cat. no. 2302700, 5 prime/Thermo Fisher, Germany)
- Accuzol isolation reagent (Cat. no. K-3090, Bioneer, Korea)

## 14. cDNA synthesis

1 $\mu$ g of total RNA is required and should be reverse transcribed to cDNA. Different kits are available on the market for this purpose. Dilute prepared cDNA 2.5x with nuclease free water. For example, to 20 $\mu$ l of cDNA add 30 $\mu$ l of nuclease free water.

## 15. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

Each sample should be examined for both BCR-ABL fusion (p190) gene and for ABL control gene in two separate sets of reactions. In BCR-ABL set consider 1 tube for each sample as

well as 6 tubes for the 5 standards (m1 to m5) and Negative control or NTC. In ABL set consider 1 tube for each sample and 5 tubes for the 4 standards (A1 to A4) and Negative control or NTC. Place required number of tubes on cold block.

**Pipette 20µl of mbcr 190 Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA or standard or control to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## 16. Devices and software

mbcr 190 RQ kit is designed to work with Rotor-Gene, StepOne and MIC.


## 17. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the p190 template file for Rotor-Gene (provided in the flash card, or accessible by kit QR code); mbcr 0.1 is for strip tubes and mbcr 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the following image.

Select a tube number containing p190 Mix for the Green channel and tube with ABL Mix for the Yellow channel.

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

Edit sample names on both mbcr and ABL pages. Remember that tubes containing mbcr mix should only be named in mbcr page and tubes containing ABL Mix should only be named in ABL page.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been

entered. Patient samples should be defined as "unknown" and no template control as "NTC" respectively.

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template (provided in the flash card, or accessible by kit QR code).

Click on Plate Setup. One negative control, five standards for mbr, four standards for ABL and a few samples are defined. You may change plate Setup using right-click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on the Start Run and save the experiment to the desired location. Instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Both mbr Mix and ABL Mix contain ROX with a final concentration of 300nM in reaction.

## 20. Data Analysis: Rotor-Gene


Before analyzing results, make sure in the sample menu, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be

defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.


Analyze data according to Rotor-Gene manual. Perform quantitative analysis for both the **mbcr (the Green channel)** and the **ABL (the Yellow channel)**. Briefly, click on Analysis menu and then under the Quantitation tab double click on cycling A. Green.

Close the pop-up window and manually set threshold at 0.05. Repeat above for the Yellow channel and set the threshold on 0.05. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

**Fig 1.** Typical mbcr graph in Green channel for Rotor-Gene

**Fig 2.** Typical ABL graph in Green channel for Rotor-Gene

Consider the following points when analyzing:


- A sample is **Positive** if it is positive in both the Green/mbcr and the Yellow/ABL channels with sigmoid graphs and CT of 20-40 for Green and 20-30 for the Yellow.
- A sample is **Negative** if it is negative in the Green/mbcr channel while it is positive in the Yellow/ABL channel with a sigmoid graph and CT of 20-27.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both the Green/mbcr and the Yellow/ABL channels. Improper extraction or test Setup could cause that.
- Results are **inconclusive** and the test should be repeated if a sample is negative in the Green/mbcr channel while it is positive in Yellow/ABL channel with sigmoid graph and CT of above 27. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the Yellow/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used. CTs higher than 27 reduce the sensitivity of the test and may result in **false negative** reports.*


## 21. Data Analysis: StepOne

Analyze data according to StepOne manual. Briefly, click on Analyze and set the threshold for both the **mbcr/FAM** and the **ABL/VIC** on 0.1.

Figures 3 and 4 represent typical graphs for the StepOne machine.

**Fig 3.** Typical mbcr graph in FAM channel for StepOne

**Fig 4.** Typical ABL graph in VIC channel for StepOne

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

Consider the following points when analyzing:

- A sample is **Positive** if it is positive in both the FAM/mbcr and the VIC/ABL channels with sigmoid graphs and CT of 20-40 for FAM and 20-30 for the VIC.
- A sample is **Negative** if it is negative in the FAM/mbcr channel while it is positive in the VIC/ABL channel with a sigmoid graph and CT of 20-27.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both the FAM/mbcr and VIC/ABL channels. Improper extraction or test set up could cause that.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in the FAM/mbcr channel while it is positive in VIC/ABL channel with sigmoid graph and CT of above 27. Improper extraction or low RNA input could cause that.

*Note: All patient samples should be positive in the VIC/ABL channel with a sigmoid graph and CT of 27 or less. ABL CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used. CTs higher than 27 reduce the sensitivity of the test and may result in **false negative** reports.*

## 22. BCR-ABL% calculation

To assess the response to therapy, BCR-ABL% value for each patient can be calculated. This kit uses NCN method for this

purpose (Beillard E. 2003, *Leukemia* 17:2474). In this method BCR-ABL% value is the BCR-ABL expression (titer) normalized by the ABL expression (titer) and then multiplied by 100.

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## 23. Analytical Sensitivity

The analytical detection limit of the kit (LOD) was determined as 7 copies/ul of BCR-ABL or 0.14% in the context of 5,000 copies/ul of ABL mRNA. This means that when a cDNA contains about 5,000 copies/ul of ABL transcripts, there is 95% probability that a sample of 0.14% positive will be detected.

## 24. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 25. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: [info@novingene.com](mailto:info@novingene.com)

## 26. Contact Information

**NovinGene ParsVira**

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Website: [www.novingene.com](http://www.novingene.com)

## 27. References

- Hochhaus, A., Baccarani, M., Silver, R.T., Schiffer, C., Apperley, J.F., Cervantes, F., Clark, R.E., Cortes, J.E., Deininger, M.W., Guilhot, F., Hjorth-Hansen, H., Hughes, T.P., Lipton, J.H., Mahon, F.X., Mayer, J. and Nicolini, F., 2020. European LeukemiaNet 2020 recommendations for treating chronic myeloid leukemia. *Leukemia*, 34.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3 2004: 190-212.
- Score, J., Calasanz, M.J., Ottman, O., Pane, F., Yeh, R.F., Sobrinho-Simões, M.A., Kreil, S., 2010. Analysis of genomic breakpoints in p190 and p210 BCR–ABL indicate distinct mechanisms of formation. *Leukemia*, 24(10), pp.1742-1750.
- Witt, M., Malgorzata Dawidowska and Szczepanski, T., 2012. *Molecular Aspects of Hematologic Malignancies*. Springer Science & Business Media.

## 28. Symbols

<b>RUO</b> Research use only	 Manufacturer	 Consult instructions for use
<b>LOT</b> Lot number	 Content sufficient for <n> tests	 Use-by date
<b>REF</b> Catalogue number	<b>SN</b> Serial number	 Temperature limit

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**