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NRAS RQ Kit Manual

Spring 2026, Version 1.0

For Real-Time PCR Detection of NRAS Mutations
For use with Rotor-Gene and StepOne
For Research Use Only

 24 (Cat# NRASRQ24)

 48 (Cat# NRASRQ48)

 NG-WI-ASL-73-100

RUO

NovinGene ParsVira

No. 20, 4th St, Gisha, Tehran, Iran 1446843434.


Table of Contents

- 1. Introduction 3
- 2. Intended Use 3
- 3. Background Information 3
- 4. Test Principle 5
- 5. Kit Contents 5
- 6. Packaging models 6
- 7. Storage and Stability 6
- 8. Product Use Limitations 6
- 9. Additionally Required Items 6
- 10. General Precaution 7
- 11. Specimen 7
- 12. DNA Extraction 7
- 13. Protocol: DNA Sample Assessment 8
- 14. Devices and software 8
- 15. Programming Real-Time PCR 9
- 16. Analysis: Sample Assessment 10
- 17. Protocol: Detection of NRAS mutations 12
- 18. NRAS Mutation Detection Analysis 13

19. Analytical Sensitivity	20
20. Disposal Method	21
21. Technical Support.....	21
22. Contact Information.....	21
23. References	22
24. Symbols.....	22

1. Introduction

NRAS RQ Kit is a ready-to-use real-time PCR assay for the detection of NRAS gene mutations in codons 12, 13, 59, and 61 in human genomic DNA. All required reagents are included in the provided PCR Mixes. All Mutation Mixes contains specific primers and probes for detecting mutations as well as an internal control based on a synthetic DNA sequence.

This kit is intended for Research Use Only (RUO).

2. Intended Use

NRAS RQ kit is intended for detecting of NRAS gene mutations in codons 12, 13, 59 and 61 in human genomic DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene and StepOne machines.

3. Background Information

The RAS gene family, consisting of HRAS, KRAS, and NRAS, functions as critical molecular switches in intracellular signaling pathways. These proteins regulate a wide range of cellular processes, including proliferation, differentiation, migration, survival, and apoptosis. Activating mutations in RAS genes impair their intrinsic GTPase activity, resulting in persistent activation of downstream signaling pathways, particularly MAPK/ERK pathway and other kinase-dependent pathways. Consequently, aberrant RAS signaling promotes uncontrolled cell growth and contributes to tumor initiation and progression in numerous human malignancies. In clinical studies, the presence of RAS mutations has been linked to distinct clinicopathological features, prognosis, and responses to therapeutic interventions.

Among the members of RAS family, NRAS plays a prominent oncogenic role in several cancers, most notably melanoma, acute myeloid leukemia (AML), and thyroid carcinoma, colorectal cancer, and endometrial cancer. NRAS gene is located on chromosome 1, and its mutations are primarily found in exon 2 (codons 12 and 13), exon 3 (codons 59 and 61), and exon 4 (codons 117 and 146). Among these alterations, mutations affecting codon Q61 are of particular clinical significance, and can be associated with reduced treatment efficacy and, in some cases, the development of drug resistance. The most important NRAS mutations are listed in Table 1.

Exon	Codon	Protein Change	Base change
2	12	G12A	35G>C
2	12	G12C	34G>T
2	12	G12D	35G>A
2	12	G12R	34G>C
2	12	G12S	34G>A
2	12	G12V	35G>T
2	13	G13A	38G>C
2	13	G13C	37G>T
2	13	G13D	38G>A
2	13	G13R	37G>C
2	13	G13S	37G>A
2	13	G13V	38G>T
3	59	A59T	175G>A
3	59	A59D	176C>A
3	61	Q61E	181C>G
3	61	Q61H	183A>C
3	61	Q61H	183A>T
3	61	Q61K	181C>A
3	61	Q61L	182A>T

3	61	Q61R	182A>G
4	117	K117N	351G>C
4	117	K117N	351G>T
4	117	K117R	350A>G
4	146	A146P	436G>C
4	146	A146T	436G>A
4	146	A146V	437C>T

Table1. NRAS Mutations

4. Test Principle

The target is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual and following reagents:

Label	Content	Quantity	Volume
Ctrl Mix	PCR Mix for quality control	2	480 µl
G12 Mix	PCR Mix to check G12 mutations	1	480 µl
G13 Mix	PCR Mix to check G13 mutations	1	480 µl
A59D Mix	PCR Mix to check A59D mutations	1	480 µl
Q61H Mix	PCR Mix to check A61H mutations	1	480 µl
Q61K Mix	PCR Mix to check A61K mutations	1	480 µl
Q61L Mix	PCR Mix to check A61L mutations	1	480 µl
Q61R Mix	PCR Mix to check A61R mutations	1	480 µl
NRAS Pos	Positive control	1	250 µl
NRAS Neg	Negative control	1	250 µl
Water	PCR Grade Water	1	200 µl

Table 2. Kit Contents

6. Packaging models

The kit is available in 24 and 48 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

9. Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and the accessory computer
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments
- Nuclease-free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves

- Cold block

10. General Precaution

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

11. Specimen

Sections of paraffin block are the proper sample for patients with solid tumor such as colorectal cancer and thyroid. If tissue is not feasible, plasma sample can be used. In hematological malignancies such as AML, peripheral blood or bone marrow aspirate are standard specimens.

12. DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)

- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 50 microliters. We recommend extracting 100 microliters or higher volume of DNA, enough to repeat the test if necessary.

Note: Extracted sample should contain 10-50 ng/ul DNA.

13. Protocol: DNA Sample Assessment

Before examining a sample for NRAS mutations, quality of DNA should be assessed. If the results are within the desired range, then the second test for detecting NRAS mutations will be performed. To qualify DNA extraction, follow steps below.

Note! Check the volume of extracted DNA sample and make sure it is more than 50 microliters and preferably about 100 microliters. Volumes less than 50 microliters are not enough to proceed.

First, thaw the **Ctrl Mix** on ice completely, mix by flickering followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for negative controls and water.

Pipette 20µl of Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and water to each tube.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.


14. Devices and software

NRAS RQ kit is designed to work with Rotor-Gene and StepOne.

15. Programming Real-Time PCR

Open the NRAS template file for Rotor-Gene or StepOne machine (provided in the flash card, or accessible by kit QR code). In RotorGene template NRAS 0.1 is for strip tubes and NRAS 0.2 is for 0.2ml tubes. Program starts.

Note: For Gain Optimisation, in the View menu select the Gain Optimisation. Adjust the setting according to the below image. Make sure to set the Tube Position to number 1 for all channels (note that Tube number 1 should contain Ctrl Mix).

Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the file to start the machine.

You can also use table 3:

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Table 3. Temperature Profile

Fluorescence should be collected at 60°C for FAM and VIC dyes. Mix contain ROX with final concentration of 300nM in reaction. Make sure that in the sample menu, the Positive controls have been defined as "Positive control". Patient samples are defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

16. Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. In the first step, results of positive and negative controls will be evaluated.

Perform quantitative analysis for both **DNA sample (the Green channel)** and **Internal control (the Yellow channel)**. Briefly, click on Analysis menu and then, under Quantitation tab double click on "Cycling A. Green". Close the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

A) Test is valid only if:

- 1) Water sample is **Negative** in the Green channel, and it is **Positive** in Yellow channel with CT of 28-31.

2) Positive control/Negative control is **Positive** in the Green channel with CT of 24-28. And it is **Positive** in the Yellow channel with CT of 28-31.

Expected results for a valid test are summarized in the Table 4.

Sample	Yellow	Green
NTC	Pos CT: 28-31	Neg
Pos/Neg Ctrl	Pos CT: 28-31	Pos CT: 24-28

Table 4. Expected results for Sample Assessment Test validity.

If above are met, then the Test results are valid and may proceed to DNA sample assessment.

B) DNA sample Quality Analysis:

- **DNA sample is qualified only if** sample is positive in the Green channel with CT of 22-28 and also positive in the Yellow channel with CT of 28-31. This sample can be further examined for NRAS mutations.

- **Sample should be diluted with water**, if a sample in is positive in the Green channel with CT of less than 22 and also positive in the Yellow channel with CT of 28-31. Dilute the sample to reach CT of 22-28. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the NRAS mutation test.

- **Result is invalid**, if a sample in the Green channel has CT above 28. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 5.

	FAM/Green	VIC/Yellow	Result
1	+ CT: 22-28	+ CT: 28-31	Valid
2	+ CT<22	+ CT: 28-31	Sample dilution
3	+ CT>28	+ CT: 27-33	Invalid
4	+ CT: 22-28	+ CT>33	Invalid

Table 5. DNA Quality Analysis

17. Protocol: Detection of NRAS mutations

To detect NRAS mutations, each sample should be examined with 8 mixes, one for DNA quantitation and one for each of NRAS mutations, each in a separate microtube. Therefore, to examine only one sample, 32 microtubes are required; 8 for the sample, 8 for the Positive control, 8 for the Negative control and 8 for water (NTC). Respectively for each extra samples 8 microtubes will be added. Therefore for 2 samples, 40 microtubes and for 3 samples 48 microtubes are required. Figure1 shows tube setup for 3 samples. Place required number of tubes on cold block organized in series of 3 for each sample.

Pipette 20µl of Ctrl Mix to each tube in the first series.

Pipette 20µl of G12 Mix to each tube in the 2nd series.

Pipette 20µl of G13 Mix to each tube in the 3th series.

Pipette 20µl of A59D Mix to each tube in the 4th series.

Pipette 20µl of Q61H Mix to each tube in the 5th series.

Pipette 20µl of Q61K Mix to each tube in the 6th series.

Pipette 20µl of Q61L Mix to each tube in the 7th series.

Pipette 20µl of Q61R Mix to each tube in the 8th series.

Continue by adding 5µl of extracted DNA, Positive control, Negative Control and water to each tube. Consider the first, second and third tube in each series for positive and Negative control and water. Then next would be for samples.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine and program it according to the section 15.

Note: If using Rotor-Gene attach the locking ring.


Fig 1. Tubes setup for Mixes and samples.

18. NRAS Mutation Detection Analysis

Before analyzing results of NRAS mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be analyzed. If the test is valid, then may proceed to NRAS mutation analysis.

Briefly, for RotorGene software click on Analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close

the pop-up window and manually set threshold at 0.1. Repeat above for the Yellow channel and set the threshold at 0.1. Also put threshold at 0.1 for FAM and VIC in StepOne.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

A) Test is valid only if:

- 1) Water sample is Negative in the Green channel with all Mixes, and
- 2) Water sample is Positive in the Yellow channel with Control Mix with a CT of 28-31, and
- 3) Positive control is Positive in the Green channel with Control Mix with a CT of 26-29, and
- 4) Pos control is Positive in the Yellow channel with Control Mix and CT of 28-31, and
- 5) Positive control is Positive in the Green channel **with all severe NRAS mutation Mixes** with a CT of 26-40, and
- 6) Sample is Positive in the Green channel with Control Mix with a CT of 22-30 and Positive in Yellow channel with CT of 28-31.

Above steps for Test validation are summarized in Figure 2.

If all above are met, then the Test is valid, and results can be analyzed further for NRAS mutations.


Fig 2. Test validation flowchart

B) Data analysis for NRAS mutations

- 1) Select the patient’s samples which are positive in the Green channel with NRAS specific mixes with a CT of 20-40, and document the CTs in the Table 6.
- 2) Calculate ΔCT for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the ΔCT in the Table 6 too.
- 4) If the ΔCT is within the valid range as mentioned in the Table 3, the sample has the NRAS mutation.
- 5) If a sample is positive for two or more NRAS mutations, it is considered positive for the mutation with the lowest ΔCT and negative for the other mutations.

Above steps for detection of NRAS mutations are summarized in Figure 3.

NRAS Mix	Valid CT	Sample CT	Valid ΔCT in RotorGene	Valid ΔCT in StepOne	Sample ΔCT	Result
Ctrl Mix	22-28		NA	NA		
G12 Mix	20-40		≤ 7.4	≤ 7.4		
G13 Mix	20-40		≤ 11.2	≤ 11.8		
A59D Mix	20-40		≤ 8.4	≤ 9.7		
Q61H Mix	20-40		≤ 8	≤ 10.2		
Q61K Mix	20-40		≤ 9.6	≤ 10.9		
Q61L Mix	20-40		≤ 12.4	≤ 12.8		
Q61R Mix	20-40		≤ 12.4	≤ 12.8		

Table 6. Valid CT and ΔCT ranges for NRAS Mixes.

As an example, If CT of a sample is 25.4 with Control Mix, and 35.1 with G12 Mix, 32.5 with G13 Mix, 36.8 with A59D Mix, 37.8 with Q61H Mix, 38.1 with Q61K Mix, 40 with Q61L Mix and 39 with Q61R Mix, then ΔCT is 9.7 for G12 (35.1 – 25.4), 7.1 for G13 (32.5-25.4), 11.4 for A59D (36.8-25.4), 12.4 for Q61H (37.8-25.4), 12.7 for Q61K (38.1-25.4), 14.6 for Q61L (40-25.4), 13.6 for Q61R (39-25.4), (Table. 7). So, results of G12 and G13 are in valid range. In this case, patient is Positive for G13.

NRAS Mix	Valid CT	Sample CT	Valid ΔCT in RotorGene	Valid ΔCT in StepOne	Sample ΔCT	Result
Ctrl Mix	22-28	25.4	NA	NA		
G12 Mix	20-40	35.1	≤ 7.4	≤ 7.4	9.7	Negative
G13 Mix	20-40	32.5	≤ 11.2	≤ 11.8	7.1	Positive
A59D Mix	20-40	36.8	≤ 8.4	≤ 9.7	11.4	Negative
Q61H Mix	20-40	37.8	≤ 8	≤ 10.2	12.4	Negative
Q61K Mix	20-40	38.1	≤ 9.6	≤ 10.9	12.7	Negative
Q61L Mix	20-40	40	≤ 12.4	≤ 12.8	14.6	Negative
Q61R Mix	20-40	39	≤ 12.4	≤ 12.8	13.6	Negative

Table 7. NRAS mutation analysis for a specific sample.

Note that a normal sample may cross react with some of the mutation mixes. However, the calculated ΔCT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in green channel for NRAS mutation mixes, will only be considered positive only if ΔCT is within the valid range.


Fig 3. Sample analysis flowchart for NRAS mutation detection

Results can also be interpreted as below;

- Sample is **Negative** for NRAS mutations, if it is negative in the Green channel with all mixes or has CT above 40.
- Sample is **Negative** for NRAS mutations if it is positive in the Green channel with CT of 20-40 for one or more NRAS mixes, and a Δ CT above valid range (which mentioned in Table 6).
- Sample is **Positive** for NRAS mutations, if sample has a valid CT and Δ CT (Table 6) in the Green channel for one of the seven NRAS mixes. Briefly,
- If a sample has a valid CT and Δ CT for more than one of the NRAS Mix, sample is Positive for a NRAS mutation which has a lowest Δ CT.

Interpretation of results for NRAS mutation test are summarized in Table 8.

Note that, if a sample is Negative with this kit in terms of NRAS mutations, the following conditions should be considered:

- Sample is Negative only for NRAS mutations which was mentioned.
- Sample is Positive for NRAS mutations in quantities below the sensitivity of the kit.
- Sample could be Positive for other NRAS mutations.

NRAS Mix	Sample CT	Δ CT In RotorGene	Sample Δ CT In StepOne	Conclusion
G12 Mix	>40	-	-	Neg for G12 mutations
	20-40	>7.4	>7.4	Neg for G12 mutations
	20-40	\leq 7.4	\leq 7.4	Pos for G12 mutations
G13 Mix	>40	-	-	Neg for G13 mutations
	20-40	>11.2	>11.8	Neg for G13 mutations
	20-40	\leq 11.2	\leq 11.8	Pos for G13 mutations
A59D Mix	>40	-	-	Neg for A59D mutations
	20-40	>8.4	>9.7	Neg for A59D mutations
	20-40	\leq 8.4	\leq 9.7	Pos for A59D mutations
Q61H Mix	>40	-	-	Neg for Q61H mutation
	20-40	>8	>10.2	Neg for Q61H mutation
	20-40	\leq 8	\leq 10.2	Pos for Q61H mutation
Q61K Mix	>40	-	-	Neg for Q61K mutation
	20-40	>9.6	>10.9	Neg for Q61K mutation
	20-40	\leq 9.6	\leq 10.9	Pos for Q61K mutation
Q61L Mix	>40	-	-	Neg for Q61L mutation
	20-40	>12.4	>12.8	Neg for Q61L mutation
	20-40	\leq 12.4	\leq 12.8	Pos for Q61L mutation
Q61R Mix	>40	-	-	Neg for Q61R mutation
	20-40	>12.4	>12.8	Neg for Q61R mutation
	20-40	\leq 12.4	\leq 12.8	Pos for Q61R mutation

Table 8. Interpretation of results for NRAS mutation.

19. Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated NRAS DNA that can be identified in the background of wild-type DNA and is mentioned in Table 9. The sensitivity depends on the DNA quantity in a sample. Maximum sensitivity is obtained when an extracted sample contains 10-50 ng/ul DNA. Also, Sample

with Control Mix has a CT of 22-28 in the Green channel. If CT range is 28-30, the sensitivity will decrease for some mutations.

Reaction	Sensitivity
G12	~0.8
G13	~0.3
A59D	~1.4
Q61H	~0.7
Q61K	~0.5
Q61L	~0.5
Q61R	~1.2

Table 9. Sensitivity of the NRAS mutation detection assay.

20. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via

Phone: +98 993-6223241

Email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 20, 4th St, Gisha St, Tehran, Iran 1446843434.

Tel: +98 21-88837393




+98 990-1813124

Email: info@novingene.com
 Website: www.novingene.com

23. References

- Kelleher, F.C. and McArthur, G.A., 2012. Targeting NRAS in melanoma. *The Cancer Journal*, 18(2), pp.132-136.
- Mackay IM. Real-time PCR in the microbiology laboratory. *Clin. Microbiol. Infect.* 2004; 10 (3): 190 – 212.
- Mulligan, G., Lichter, D.I., Di Bacco, A., Blakemore, S.J., Berger, A., Koenig, E., 2014. Mutation of NRAS but not KRAS significantly reduces myeloma sensitivity to single-agent bortezomib therapy. *Blood, The Journal of the American Society of Hematology*, 123(5), pp.632-639.
- Muñoz-Couselo, E., Adelantado, E.Z., Ortiz, C., García, J.S. and Perez-Garcia, J., 2017. NRAS-mutant melanoma: current challenges and future prospect. *OncoTargets and therapy*, pp.3941-3947.
- Yang, X. and Wu, H., 2024. ‘RAS signaling in carcinogenesis, cancer therapy and resistance mechanisms’, *Journal of Hematology & Oncology*, 17(1).

24. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com